




BTL Biotechno Labs Pvt. Ltd.

Authorized Channel Partner of VIDIA

VIDIA is a private Czech biotechnology company founded in 1991 by experienced research scientists. Specializing in the development, manufacture, and sale of in vitro diagnostic (IVD) kits, VIDIA serves both human medicine and environmental monitoring sectors. All VIDIA IVD products are CE certified, designed for high stability, long shelf life, and ease of use. Their product portfolio includes ELISA kits, immunofluorescence assays (IFAs), line immunoassays (LIAs), apid tests, peptide and antigen conjugate synthesis services.


BTL Biotechno Labs Pvt. Ltd. is an authorised channel partner of VIDIA India, dedicated to being the most reliable supplier of life science research products across the Indian market. We are committed to supporting the scientific community by offering innovative, high-quality, and affordable research solutions. As a distributor of VIDIA in India, we provide their complete product range to researchers at the best price, ensuring access to cutting-edge tools for scientific advancements.

VARICELLA ZOSTER VIRUS



VIDIA kits

REF	Product	Evaluation	Incubation	Sample	Number of tests
ODZ-087	anti-VZV IgG (CSF)	semiquant., 5ST quant.	30' / 30' / 15'	serum, plasma, cerebrospinal fluid	96
ODZ-168	anti-VZV IgG	semiquant.	30' / 30' / 15'	serum, plasma	96
ODZ-197	anti-VZV IgM	semiquant.	30' / 30' / 15'	serum, plasma	96
ODZ-233	anti-VZV IgG (CSF) a avidity IgG	semiquant., quant.	30' / 10' / 30' / 15'	serum, plasma, cerebrospinal fluid	96
ODZ-284	anti-VZV IgA	semiquant.	30' / 30' / 15'	serum, plasma	96



TICK-BORNE ENCEPHALITIS VIRUS – TBEV

VIDIA kits

REF	Products
ODZ-192	ELISA-VIDITEST anti-TBEV IgG (CSF)
ODZ-170	ELISA-VIDITEST anti-TBEV IgG and IgG avidity
ODZ-194	ELISA-VIDITEST anti-TBEV IgM


BORRELIA BURGDORFERI SENSU LATO



VIDIA kits


REF	Products
ODZ-398	ELISA -VIDITEST anti-Borrelia recombinant IgG + VIse (CSF)
ODZ-398/5ST	ELISA -VIDITEST anti-Borrelia recombinant IgG + VIse (CSF)
ODZ-399	ELISA -VIDITEST anti-Borrelia recombinant IgM (CSF)
ODZ -399/5ST	ELISA -VIDITEST anti-Borrelia recombinant IgM (CSF)

CYTOMEGALOVIRUS



VIDIA kits


REF	Product	Evaluation	Incubation	Sample	Number of tests
ODZ-102	anti-CMV IgG	semiquant.	30' / 30' / 15'	serum, plasma	96
ODZ-102/5ST	anti-CMV IgG (CSF)	semiquant., quant.	30' / 30' / 15'	serum, plasma, cerebrospinal fluid	96
ODZ-177	anti-CMV IgG a avidity IgG	semiquant.	30' / 10' / 30' / 15'	serum, plasma	48
ODZ-402	anti-CMV IgM	semiquant.	30' / 30' / 15'	serum, plasma	96
ODZ-164	anti-CMV IgA	semiquant.	30' / 30' / 15'	serum, plasma	96



HUMAN HERPESVIRUS 6

VIDIA kits


REF	Product	Evaluation	Incubation	Sample	Number of tests
ODZ-235	anti-HHV-6 IgG	semiquant.	60' / 60' / 20'	serum	96
ODZ-344	anti-HHV-6 IgG (CSF)	semiquant., 5 ST quant.	60' / 60' / 20'	serum, cerebrospinal fluid	96
ODZ-345	anti-HHV-6 IgM	semiquant.	60' / 60' / 20'	serum	96



POLYOMAVIRUS JC

VIDIA kits

REF	Product	Evaluation	Incubation	Sample	Number of tests
ODZ-450	ELISA-VIDITEST anti-JCV IgG	semiquant., quant.	60' / 60' / 10'	serum, plasma	96



POLYOMAVIRUS BK


VIDIA kits

REF	Product	Evaluation	Incubation	Sample	Number of tests
ODZ-405	anti-BKV IgG	semiquant.	30' / 30' / 10'	serum, plasma	96
ODZ-450	anti-JCV IgG	semiquant., quant.	60' / 60' / 10'	serum, plasma	96

COMPLEMENT FACTOR H

VIDIA kits

REF	Product	Evaluation	Incubation	Sample	Number of tests
ODZ-166	anti-complement factor H	quant.	60' / 60' / 20'	serum, plasma	48




ASPERGILLUS FUMIGATUS



VIDIA kits

REF	Product	Evaluation	Incubation	Sample	Number of tests
ODZ-512	anti-Aspergillus fumigatus IgG	semiquant., quant.	30' / 30' / 15'	serum, plasma	96
ODZ-513	anti-Aspergillus fumigatus IgM	semiquant., quant.	30' / 30' / 15'	serum, plasma	96
ODZ-514	anti-Aspergillus fumigatus IgA	semiquant., quant.	30' / 30' / 15'	serum, plasma	96




CHLAMYDIA PNEUMONIAE

VIDIA kits

REF	Product	Evaluation	Incubation	Sample	Number of tests	VIDIMAT
ODZ-457	anti-Chlamydia pneumoniae IgG	semiquant.	30' / 30' / 15'	serum, plasma	96	✓
ODZ-458	anti-Chlamydia pneumoniae IgM	semiquant.	30' / 30' / 15'	serum, plasma	96	✓
ODZ-459	anti-Chlamydia pneumoniae IgA	semiquant.	30' / 30' / 15'	serum, plasma	96	✓


PSEUDOMONAS AERUGINOSA



VIDIA kits

REF	Products	Evaluation	Incubation	Sample	Numbers of tests	THUNDERBOLT
ELISA-VIDITEST						
ODZ-526	anti-Pseudomonas aeruginosa IgG	quant.	60' / 60' / 10'	serum, plasma	96	✓ NEW

STEP-BY-STEP PRINCIPLE



BTL Biotechno Labs Pvt. Ltd.

613-A, Ansal Chamber-II, Bhikaji Cama Place, New Dehi - 110066
 Phone : +91 8860924629, 7291852429
 E-mail : info@biotechnolabs.com
 web : www.biotechnolabs.com

VIDIA spol. s r.o.

Nad Safinou II 365, 252 50 Vestec
 Czech Republic
 phone: +420 261 090 565
 fax: +420 261 090 566
 e-mail: info@vidia.cz
 www.vidia.cz

