

Neuroscience Research

- Recombinant Proteins
- Antibodies
- Elisa Kits
- Cytokines for Neural Cell Culture



Authorized Distributor:

BTL Biotechno Labs Pvt. Ltd.

613-A, Ansal Chamber-II, Bhikaji Cama Place, New Delhi-110066

Mob : +91- 8860924349; 07291852429

E-mail : info@biotechnolabs.com , Web : www.biotechnolabs.com



YOUR GOOD PARTNER IN BIOLOGY RESEARCH!

Postal Address : 7505 Fannin St. Ste 610-322, Houston, TX 77054, USA

Tel : 301-363-4651 (Available 9 a.m. to 5 p.m. CST from Monday to Friday)

Email: support@cusabio.com Web: www.cusabio.com

Neuroscience research is important for understanding human brain function, treating neurological disorders, and advancing related drug development. Neurological disorders are diseases that affect the brain, spinal cord, peripheral nerves, and other components of the nervous system. These diseases include, but are not limited to, cerebrovascular disease, Parkinson's disease, Alzheimer's disease, multiple sclerosis, epilepsy, and others.

CUSABIO provides neuroscience research related proteins, antibodies, ELISA kits products, as well as cytokines commonly used in neural cell culture, to help related mechanisms of functional research and drug discovery and development.

● Neurological Diseases and Popular Research Targets

Alzheimer's Disease (AD)

Alzheimer's Disease is a chronic, progressive neurodegenerative disease characterized by memory loss, cognitive dysfunction, and behavioral and psychological symptoms. It is the most common type of dementia and is thought to be caused by the abnormal accumulation of amyloid plaques in the brain, resulting in massive neuronal death, and abnormal phosphorylation of Tau proteins.

Currently, drug research targets for AD are mainly focused on A β , Tau proteins, BACE1, etc. Other potential targets are also being explored, such as α -synuclein and RAGE inhibitors.

ACHE	AGER	AKT1	APOE	APP	BACE1
BCHE	BDNF	CAR1	GFAP	IER3	IL6
INS	MAPT	PSEN1	RENBP	SNCA	TNF

Parkinson's Disease (PD)

Parkinson's Disease is a neurodegenerative disorder that primarily affects motor function, with characteristic symptoms including resting tremor, muscle rigidity, bradykinesia, and balance disorders. The pathogenesis of Parkinson's disease is complex and involves multiple factors, including loss of dopamine neurons, protein misfolding and aggregation, mitochondrial dysfunction, neuroinflammation, and oxidative stress.

Currently, drug research targets for Parkinson's disease focus on α -synuclein aggregation, LRRK2 gene mutations, PAAN/MIF nuclease activity, mitochondrial dysfunction, autophagy mechanisms, and specific dopaminergic neuron subpopulations discovered based on single-cell sequencing technology.

AKT1	APOE	APP	BDNF	GBA1	IER3
IL1A	IL6	INS	LRRK2	MAOB	MAPT
NFKB1	PARK7	PINK1	PRKN	SLC6A3	SNCA
SOD1	TH	TNF			

Amyotrophic Lateral Sclerosis (ALS)

Amyotrophic Lateral Sclerosis, also known as acromegaly, is a fatal neurodegenerative disease characterized by the progressive loss of motor neurons in the brain and spinal cord, resulting in gradual muscle weakness and atrophy that ultimately affects motor, communication, swallowing, and respiratory functions, and ultimately death.

Currently, its research targets focus on genes and proteins such as SOD1, C9ORF72, TDP-43, FUS, OPTN, TBK1, and RIPK1, as well as exploring new therapies that utilize a combination of antisense oligonucleotides (ASOs), stem cells and genetic engineering.

APP	ATXN2	BDNF	C9orf72	FUS	GFAP
-----	-------	------	---------	-----	------

GRN	HTT	IER3	IL6	MAPT	OPTN
SNCA	SNCA	SOD1	SQSTM1	TARDBP	TBK1
TNF	VCP				

Multiple Sclerosis/Multiple Sclerosis (MS)

Multiple Sclerosis is a chronic inflammatory disease of the central nervous system (CNS) resulting in demyelination and neurodegeneration. The clinical manifestations of MS are varied and include motor deficits, sensory abnormalities, visual disturbances, psychiatric symptoms, and cognitive impairment.

Drug research targets in MS have focused on B and T cells, sphingosine-1-phosphate receptor (S1PR), Nrf2 pathway, potassium channels, α 2-adrenergic receptor, Syk kinase, S1P1 receptor, CD20-positive B cells, immunosuppressive drugs, mitochondrial protective agents, ion channels, and glutamate-mediated excitotoxicity, which involve immunomodulation, neuroprotection, and cellular damage.

CD4	CD8A	CSF2	FOXP3	IER3	IFNB1
IFNG	IL10	IL17A	IL1A	IL1B	IL2
IL6	KRT20	MBP	MOG	TNF	

Huntington's Disease (HD)

Huntington's Disease is a rare autosomal dominant neurodegenerative disorder with typical manifestations including progressively progressive choreiform movements, cognitive deficits, and psychiatric behavioral abnormalities. The pathogenesis of Huntington's disease is primarily related to an abnormal repeat of the CAG trinucleotide on the Huntington gene (HTT), which leads to the production of polyglutamine chains, which in turn encode an abnormal Huntington's protein (mHTT), which accumulates in the central nervous system and leads to neuronal damage and death.

Pharmacological research in Huntington's disease has focused on reducing the expression and levels of mutant Huntington's protein (mHTT), activating cellular autophagy pathways to promote mHTT degradation, and symptomatic treatment to improve motor and psychiatric symptoms.

AKT1	APP	BDNF	GFAP	HTT	IER3
INS	MAPT	MTOR	NFKB1	PPP1R1B	PSEN1
SNCA	SOD1	TARDBP	TNF	TP53	


CUSABIO Featured Products Related to Neuroscience

➔ Recombinant Proteins	03
➔ Antibodies	17
➔ Elisa Kits	60
➔ Cusabio Cytokines for Neural Cell Culture	72

CUSABIO Featured Products Related to Neuroscience


Recombinant Proteins

Recombinant Human Interleukin-17A (IL17A) CSB-BP624104HU(M)




Measured by its binding ability in a functional ELISA. Immobilized Human IL17A at 2 µg/ml can bind Anti-IL17A recombinant antibody (CSB-RA624104MA1HU), the EC₅₀ is 1.818-2.170 ng/mL.

Recombinant Human Interleukin-6(IL6) CSB-YP011664HU




Measured by its binding ability in a functional ELISA. Immobilized Human IL6 at 2 µg/ml can bind Anti-IL6 recombinant antibody (CSB-RA011664MA1HU), the EC₅₀ is 35.80-41.82 ng/mL.

Recombinant Human Microtubule-associated protein tau (MAPT) CSB-MP013481HU(F8)




Measured by its binding ability in a functional ELISA. Immobilized Human MAPT at 2 µg/ml can bind Anti-MAPT recombinant antibody (CSB-RA013481A1HU), the EC₅₀ is 4.547-6.284 ng/mL.

Recombinant Human Angiopoietin-2 (ANGPT2) CSB-MP001707HU(A4)




Measured by its binding ability in a functional ELISA. Immobilized Human ANGPT2 at 2 µg/ml can bind anti-ANGPT2 recombinant antibody (CSB-RA001707MA01HU), the EC₅₀ is 0.6666-0.8876 ng/mL.

Recombinant Human Programmed cell death 1 ligand 1 (CD274) CSB-MP878942HU1




Measured by its binding ability in a functional ELISA. Immobilized PD-L1 at 2 µg/ml can bind Anti-PD-L1 mouse monoclonal antibody (CSB-MA878942A1m, antigen from E.coli), the EC₅₀ of human PD-L1 protein is 1.252-1.653 ng/mL.

Recombinant Human CD276 antigen (CD276) CSB-MP733578HU




Measured by its binding ability in a functional ELISA. Immobilized PD-1 at 2 µg/ml can bind Anti-PD-1 recombinant antibody, the EC₅₀ of human PD-1 protein is 6.087-7.854 ng/mL.

Recombinant Human B-lymphocyte antigen CD20 (MS4A1) CSB-MP015007HU



Measured by its binding ability in a functional ELISA. Immobilized human CD20 at 2 µg/ml can bind Anti-CD20 recombinant antibody (CSB-RA015007A1HU), the EC₅₀ is 3.243-7.085 ng/mL.

Recombinant Human Glucagon-like peptide 1 receptor (GLP1R) CSB-MP009514HUb1



Measured by its binding ability in a functional ELISA. Immobilized Human GLP1R at 2 µg/ml can bind Anti-GLP1R recombinant antibody (CSB-RA009514MA1HU), the EC₅₀ is 54.54-94.23 ng/mL.

Target	Organism	Source	Tag Info	Code
ABCD1	Human	in vitro E.coli expression system	C-terminal 10xHis-tagged Active	CSB-CF001068HU
Ache	Rat	E.coli	N-terminal 6xHis-tagged	CSB-EP001154RA
ACVR2B	Human	E.coli	N-terminal 6xHis-B2M-tagged Active	CSB-EP623829HU
ACVR2B	Human	Mammalian cell	C-terminal 6xHis-Fc-tagged	CSB-AP005701HU
AGER	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP001441HU
AKT1	Human	E.coli	N-terminal GST-tagged	CSB-EP001553HU
Akt1	Rat	E.coli	N-terminal 10xHis-tagged	CSB-EP001553RA
Alb	Mouse	Baculovirus	C-terminal 9xHis-tagged	CSB-BP001561MO
ALB	Cynomolgus monkey	Mammalian cell	N-terminal 6xHis-Myc-tagged	CSB-MP001561MOV
ALB	Dog	Yeast	N-terminal 6xHis-tagged	CSB-YP001561DO
Alb	Mouse	Yeast	N-terminal 6xHis-tagged	CSB-YP001561MO
ALB	Cynomolgus monkey	Yeast	N-terminal 6xHis-tagged	CSB-YP001561MOV
ALB	Bovine	Yeast	C-terminal 6xHis-Myc-tagged	CSB-YP001561BO
ALB	Chicken	Yeast	C-terminal 6xHis-tagged	CSB-YP001561CH
ALB	Chicken	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP001561CH
Alb	Rat	E.coli	N-terminal 6xHis-tagged and C-terminal Avi-tagged	CSB-EP001561RA
ANGPT2	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP001707HU
Angpt2	Rat	Baculovirus	C-terminal 6xHis-tagged	CSB-BP001707RA
ANGPT2	Human	Mammalian cell	C-terminal 10xHis-tagged Active	CSB-MP001707HU(A4)
ANGPT2	Dog	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-MP001707DO
Angpt2	Rat	E.coli	N-terminal 6xHis-tagged	CSB-EP001707RAa0
Apoe	Mouse	Baculovirus	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-BP001936MO
Apoe	Mouse	E.coli	N-terminal 6xHis-tagged	CSB-EP001936MO
APOE	Rabbit	E.coli	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged	CSB-EP001936RBb3
APOE	Rabbit	E.coli	N-terminal 10xHis-B2M-tagged and C-terminal Myc-tagged	CSB-EP001936RBb7
Apoe	Mouse	Mammalian cell	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-MP001936MO
Apoe	Mouse	Yeast	N-terminal 6xHis-tagged	CSB-YP001936MO
Apoe	Mouse	E.coli	C-terminal 6xHis-tagged	CSB-EP001936MOc7

Target	Organism	Source	Tag Info	Code
APOE	Cynomolgus monkey	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP001936MOV
ApoE	Rat	E.coli	N-terminal 6xHis-tagged	CSB-EP001936RA
ATXN2	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP859520HU
BDNF	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP002655HU
BDNF	Human	E.coli	N-terminal GST-tagged	CSB-EP002655HU1
Bdnf	Rat	E.coli	N-terminal 6xHis-tagged	CSB-RP166394r
BDNF	Human	E.coli	Tag-Free Active	CSB-AP003781HU
C1S	Human	E.coli	N-terminal 10xHis-tagged	CSB-EP003657HU
C3	Mouse	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP365507MO
C3	Bovine	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP647263BO
C3	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP355590HU
C5	Rat	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP003995RA
C5	Pig	Yeast	N-terminal 6xHis-tagged	CSB-YP003995PI
CD19	Human	Mammalian cell	C-terminal Fc-tagged Active	CSB-AP005061HU
CD274	Human	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP878942HU2
CD274	Human	Mammalian cell	C-terminal hFc-tagged Active	CSB-MP878942HU1
CD274	Dog	Mammalian cell	C-terminal 10xHis-tagged	CSB-MP5015DO1
CD276	Human	Mammalian cell	C-terminal hFc-Myc-tagged Active	CSB-MP733578HU
CD276	Human	Mammalian cell	C-terminal hFc-Avi-tagged	CSB-MP733578HU-B
CD33	Human	E.coli	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged	CSB-EP004925HU
CD33	Human	Mammalian cell	C-terminal hFc-Myc-tagged Active	CSB-MP004925HU
CD38	Human	Mammalian cell	C-terminal hFc-tagged	CSB-MP004929HU1
CD38	Human	Yeast	C-terminal 6xHis-tagged	CSB-YP004929HU1
CD38	Human	E.coli	N-terminal GST-tagged	CSB-EP004929HU2
CD3E	Human	E.coli	N-terminal GST-tagged	CSB-EP004931HU(A4)
CD3E	Human	Yeast	N-terminal 6xHis-tagged	CSB-YP004931HU
Cd3e	Mouse	Yeast	N-terminal 6xHis-tagged	CSB-YP004931MO
Cd3e	Mouse	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP004931MO

Target	Organism	Source	Tag Info	Code
CD3E	Pig	E.coli	N-terminal 6xHis-tagged	CSB-EP801781PI
Cd3e	Golden hamster	E.coli	N-terminal 6xHis-tagged and C-terminal Strep II-tagged	CSB-EP4828MRG
CD40	Human	Mammalian cell	C-terminal Fc-tagged Active	CSB-AP005181HU
CD40	Human	Mammalian cell	C-terminal hFc-tagged Active	CSB-MP004936HU1
CD40	Pig	E.coli	C-terminal 6xHis-tagged	CSB-EP851448PI
CD40LG	Rhesus macaque	E.coli	N-terminal 6xHis-tagged	CSB-EP004937MOW
CD40LG	Rabbit	E.coli	N-terminal 6xHis-tagged	CSB-EP004937RB
CD40LG	Rabbit	Yeast	N-terminal 6xHis-tagged	CSB-YP004937RB
CD40LG	Human	E.coli	Tag-Free Active	CSB-AP005191HU
CD40LG	Human	Mammalian cell	N-terminal hFc-Flag-tagged Active	CSB-MP004937HU3
CD40LG	Human	Mammalian cell	N-terminal hFc-Avi-tagged	CSB-MP004937HU3-B
CD47	Human	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-AP005201HU
CD47	Human	Mammalian cell	C-terminal hFc-tagged Active	CSB-MP004940HU
CD52	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP004943HU
CD8A	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP004966HU
CD8A	Human	Mammalian cell	N-terminal 6xHis-Myc-tagged	CSB-MP004966HU
CD8A	Human	Yeast	N-terminal 6xHis-tagged	CSB-YP004966HU
Cd8a	Golden hamster	E.coli	N-terminal 6xHis-tagged and C-terminal Strep II-tagged	CSB-EP4831MRG1
CSF2	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP006045HU
CSF2	Human	Yeast	N-terminal 6xHis-tagged	CSB-YP006045HU
CSF2	Human	E.coli	Tag-Free Active	CSB-AP003581HU
CSF2	Human	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-AP003601HU
Csf2	Mouse	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-AP003651MO
Csf2	Mouse	E.coli	Tag-Free Active	CSB-AP003661MO
CSF2RA	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP006046HU
CSF2RA	Human	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-AP003591HU
Ctla4	Mouse	E.coli	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged	CSB-EP006163MO1
CTLA4	Guinea pig	Yeast	C-terminal 6xHis-tagged Active	CSB-AP005121GU

Target	Organism	Source	Tag Info	Code
Ctla4	Mouse	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-AP005301MO
CTLA4	Human	Mammalian cell	C-terminal hFc-tagged Active	CSB-MP006163HU1
DAG1	Rabbit	E.coli	N-terminal 6xHis-tagged	CSB-EP636735RB
DAG1	Human	Baculovirus	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-BP613486HU
DLG4	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP006938HU
Dmd	Rat	Yeast	C-terminal 6xHis-tagged	CSB-YP006963RA
EGFR	Human	E.coli	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged	CSB-EP007479HU3
EGFR	Human	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP007479HU
EGFR	Human	Mammalian cell	N-terminal 10xHis-tagged and C-terminal Myc-tagged Active	CSB-MP007479HU
EPOR	Dog	E.coli	C-terminal 6xHis-tagged	CSB-EP007744DO
Epor	Rat	E.coli	N-terminal 6xHis-SUMO-tagged and C-terminal 6xHis-tagged	CSB-EP007744RAa3
EPOR	Dog	E.coli	N-terminal 6xHis-GB1-tagged and C-terminal 6xHis-tagged	CSB-EP007744DOg9
Epor	Mouse	Baculovirus	C-terminal 10xHis-tagged	CSB-BP007744MO1
EPOR	Human	Mammalian cell	C-terminal 10xHis-tagged	CSB-MP007744HU1
ERBB2	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP007763HU
ERBB2	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP007763HU1
ERBB2	Human	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP007763HU3
ERBB2	Human	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-MP007763HU
F10	Human	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP007915HU
F10	Mouse	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP007915MO
FCGRT	Human	E.coli	N-terminal GST-tagged	CSB-EP008545HU
FCGRT	Human	Mammalian cell	N-terminal 6xHis-tagged	CSB-MP008545HU
FCGRT	Human	Yeast	N-terminal 6xHis-tagged	CSB-YP008545HU
FN1	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP008759HU1
FN1	Human	E.coli	Tag-Free Active	CSB-AP005751HU
FOXP3	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP861191HU1
FUS	Human	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP009069HU
FUS	Human	E.coli	N-terminal 10xHis-MBP-tagged	CSB-EP009069Hue7

Target	Organism	Source	Tag Info	Code
Fxn	Mouse	E.coli	N-terminal 6xHis-tagged	CSB-EP009086MOa0
Fxn	Rat	E.coli	Tag-Free	CSB-EP009086RAe1
FxN	Human	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP613687HU(A4)b1
FXN	Cynomolgus monkey	E.coli	N-terminal 6xHis-tagged	CSB-EP810224MOV
Fxn	Mouse	Yeast	N-terminal 6xHis-tagged	CSB-YP009086MO
gag	MoMLV	E.coli	N-terminal GST-tagged	CSB-EP356017MHKe0
gag	GLV	E.coli	N-terminal 6xHis-tagged	CSB-EP724163GBAS
gag	FIV	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP323409FAK
Galc	Mouse	in vitro E.coli expression system	Tag-Free	CSB-CF009196MO
GALC	Human	in vitro E.coli expression system	Tag-Free	CSB-CF009196HU(F1)e1
GDNF	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP009356HU
GDNF	Human	E.Coli	Tag-Free Active	CSB-AP002881HU
gfap	Zebrafish	E.coli	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged	CSB-EP009369DIL
GFAP	Mouse	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP009369MO
Gfap	Mouse	Baculovirus	Tag-Free	CSB-BP009369MOe1
GFAP	Pongo pygmaeus abelii	E.coli	N-terminal 10xHis-tagged	CSB-EP009369PYX
GHR	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP009411HU
GHR	Pongo pygmaeus abelii	E.coli	C-terminal 6xHis-tagged	CSB-EP009411MOW
GHR	Pig	E.coli	C-terminal 6xHis-tagged	CSB-EP009411PI
GHR	Rabbit	in vitro E.coli expression system	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-CF009411RB1
GHR	Human	Mammalian cell	C-terminal hFc-tagged Active	CSB-MP009411HU
GHR	Human	Mammalian cell	C-terminal mFc-Avi-tagged Active	CSB-MP009411HUj1-B
GHR	Dog	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP009411DO
GLP1R	Human	Baculovirus	N-terminal 10xHis-tagged	CSB-BP009514HU
GLP1R	Human	in vitro E.coli expression system	N-terminal 10xHis-tagged	CSB-CF009514HU(A4)
GLP1R	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP009514HU
Glp1r	Mouse	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP009514MO1
Glp1r	Rat	E.coli	N-terminal 6xHis-tagged	CSB-EP009514RA1

Target	Organism	Source	Tag Info	Code
GLP1R	Human	Mammalian cell	N-terminal 6xHis-tagged	CSB-MP009514HU
GLP1R	Human	Yeast	N-terminal 6xHis-tagged	CSB-YP009514HU
Glp1r	Rat	Yeast	N-terminal 6xHis-tagged	CSB-YP009514RA1
GLP1R	Human	Mammalian cell	N-terminal 10xHis-tagged and C-terminal Myc-tagged Active	CSB-MP009514HUb1
HEXB	Human	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP010316HU
IFNB1	Human	E.coli	N-terminal GST-tagged	CSB-EP011048HU
IFNB1	Human	E.coli	Tag-Free	CSB-EP011048Hue1
IFNG	Chicken	E.coli	N-terminal GST-tagged	CSB-EP011050CH
IFNG	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP011050HU
IFNG	Sheep	E.coli	N-terminal 6xHis-tagged	CSB-EP011050SH
Ifng	Guinea pig	E.coli	N-terminal 6xHis-tagged	CSB-EP1470GP1
IFNG	Woodchuck	Yeast	N-terminal 6xHis-tagged	CSB-YP011050MQG
IFNG	Sheep	Yeast	N-terminal 6xHis-tagged	CSB-YP011050SH
IFNG	Goat	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP011050GO
IFNG	Human	E.Coli	Tag-Free Active	CSB-AP002801HU
IFNG	Rhesus macaque	E.Coli	Tag-Free Active	CSB-AP003241MOW
IFNG	Human	E.coli	Tag-Free Active	CSB-AP004201HU
IFNG	Cynomolgus monkey	E.coli	C-terminal 6xHis-tagged	CSB-EP011050MOV
IFNG	Cynomolgus monkey	Yeast	C-terminal 6xHis-tagged	CSB-YP011050MOV
IFNG	Golden hamster	E.coli	N-terminal 6xHis-tagged and C-terminal Strep II-tagged	CSB-EP011050MRGm1
IFNG	Dog	E.coli	C-terminal 6xHis-tagged	CSB-EP011050DO
IGF1	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP356436HU
IGF1	Human	E.coli	N-terminal GST-tagged	CSB-EP356436Hue0
Igf1	Rat	E.coli	N-terminal GST-tagged	CSB-EP362300RA
IGF1	Human	E.coli	Tag-Free Active	CSB-AP003721HU
IGF1	Bovine	Yeast	N-terminal hFc-tagged	CSB-YP011086BO
IGF1R	Human	E.coli	N-terminal GST-tagged	CSB-EP011087HU
IGF1R	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP011087Hua0

Target	Organism	Source	Tag Info	Code
IL10	Rhesus macaque	Baculovirus	N-terminal 10xHis-tagged and C-terminal myc-tagged	CSB-BP011580MOW
IL10	Rhesus macaque	E.coli	N-terminal 6xHis-B2M-tagged	CSB-EP011580MOW
IL10	Rhesus macaque	Mammalian cell	C-terminal 10xHis-tagged	CSB-MP011580MOW
IL10	Rhesus macaque	Mammalian cell	C-terminal hFc-tagged	CSB-MP011580MOWd9
IL10	Human	E.Coli	Tag-Free Active	CSB-AP001791HU
IL10	Mouse	E.coli	Tag-Free Active	CSB-AP004841MO
IL10	Golden hamster	E.coli	N-terminal 6xHis-tagged and C-terminal Strep II-tagged	CSB-EP4842MRG
IL12B	Golden hamster	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP817228MRG
IL12B	Human	Yeast	C-terminal 6xHis-tagged	CSB-YP011587HU
IL17A	Rhesus macaque	E.coli	N-terminal GST-tagged	CSB-EP011597MOV
IL17A	Human	E.coli	Tag-Free Active	CSB-AP004331HU
IL17A	Human	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-AP004341HU
IL17a	Mouse	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-AP004781MO
IL17A	Human	Baculovirus	N-terminal 6xHis-tagged Active	CSB-BP624104HU(M)
IL17a	Rat	E.coli	N-terminal 6xHis-tagged	CSB-EP717213RA
IL17A	Pig	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP736792PI
IL17A	Pig	E.coli	C-terminal 6xHis-tagged	CSB-EP736792Pic7
IL1A	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP011613HU
IL1A	Rhesus macaque	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP011613MOW
IL1A	Human	Yeast	C-terminal 6xHis-tagged	CSB-YP011613Huc7
IL1A	Human	E.coli	N-terminal MBP-tagged and C-terminal 6xHis-Avi-tagged	CSB-EP011613HU-B
IL1A	Rhesus macaque	E.Coli	Tag-Free Active	CSB-AP003071MOW
IL1a	Mouse	E.Coli	Tag-Free Active	CSB-AP003251MO
IL1a	Golden hamster	E.coli	N-terminal 6xHis-tagged and C-terminal Strep II-tagged	CSB-EP4839MRG
IL1B	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP011614HU
IL1B	Human	E.coli	N-terminal GST-tagged	CSB-EP011614Hue0
IL1B	Sheep	E.coli	N-terminal 6xHis-tagged	CSB-EP011614SH
il1b	Zebrafish	E.coli	N-terminal 6xHis-tagged	CSB-EP2346DIL

Target	Organism	Source	Tag Info	Code
il1b	Zebrafish	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP2346DILb1
IL1B	Sheep	Yeast	N-terminal 6xHis-tagged	CSB-YP011614SH
IL1B	Rhesus macaque	E.Coli	Tag-Free Active	CSB-AP003081MOW
Il1b	Mouse	E.Coli	Tag-Free Active	CSB-AP003261MO
IL1B	Horse	E.Coli	Tag-Free Active	CSB-AP003481HO
IL1B	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP011614HUa0
Il1b	Golden hamster	E.coli	N-terminal 6xHis-tagged and C-terminal Strep II-tagged	CSB-EP4840MRG
IL2	Human	E.coli	C-terminal 6xHis-tagged	CSB-EP011629HUC7
IL2	Rhesus macaque	E.coli	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged	CSB-EP011629MOW
IL2	Rabbit	E.coli	N-terminal 6xHis-tagged	CSB-EP011629RB
IL2	Human	Yeast	N-terminal 6xHis-tagged	CSB-YP011629HU
IL2	Human	E.Coli	Tag-Free Active	CSB-AP001691HU
IL2	Human	E.coli	Tag-Free Active	CSB-AP004381HU
IL2	Human	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-AP004501HU
IL2	Human	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-AP004711HU
Il2	Mouse	E.coli	Tag-Free Active	CSB-AP004821MO
IL2	Llama	E.coli	N-terminal 6xHis-tagged	CSB-EP769693LBS
IL2	Human	Mammalian cell	N-terminal hFc-tagged Active	CSB-MP011629HU
IL2RA	Human	Mammalian cell	C-terminal 10xHis-tagged Active	CSB-MP011649HU3
IL4R	Pig	E.coli	N-terminal 6xHis-tagged	CSB-EP771258PI
IL4R	Human	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-AP004271HU
IL4R	Human	Mammalian cell	C-terminal Fc-tagged Active	CSB-AP004281HU
IL4R	Human	E.coli	C-terminal 6xHis-tagged	CSB-EP011661HU1
IL4R	Human	Mammalian cell	C-terminal 10xHis-tagged	CSB-MP011661HU1
IL4R	Human	Mammalian cell	C-terminal 10xHis-Avi-tagged Active	CSB-MP011661HU2-B
IL6	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP011664HU
IL6	Sheep	E.coli	N-terminal GST-tagged	CSB-EP011664SH
Il6	Rabbit	Yeast	N-terminal 6xHis-tagged	CSB-YP863071RB

Target	Organism	Source	Tag Info	Code
IL6	Rhesus macaque	E.Coli	Tag-Free Active	CSB-AP003101MOW
IL6	Human	E.coli	Tag-Free Active	CSB-AP004301HU
IL6	Llama	E.coli	C-terminal 6xHis-tagged	CSB-EP803040LBS
IL6	Human	Yeast	N-terminal 6xHis-tagged Active	CSB-YP011664HU
IL6R	Rabbit	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP5136RB1
INS	Bovine	E.coli	N-terminal 6xHis-KSI-tagged	CSB-EP011742BO
INSR	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP011753HU
LRRK2	Human	E.coli	C-terminal 6xHis-tagged	CSB-EP722493HU1
MAOB	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP013399HU
MAPT	Rhesus macaque	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP013481MOW
Mapt	Rat	E.coli	N-terminal 6xHis-tagged	CSB-EP013481RA
Mapt	Mouse	Yeast	N-terminal 6xHis-tagged	CSB-YP013481MO
Mapt	Rat	Yeast	N-terminal 6xHis-tagged	CSB-YP013481RA
MAPT	Human	E.coli	N-terminal 10xHis-tagged	CSB-EP013481HU(F8)
MAPT	Rhesus macaque	Mammalian cell	N-terminal 10xHis-tagged Active	CSB-MP013481MOW
Mapt	Mouse	Mammalian cell	N-terminal 10xHis-tagged Active	CSB-MP013481MO
Mapt	Rat	Mammalian cell	N-terminal 10xHis-tagged Active	CSB-MP013481RA
MAPT	Human	Mammalian cell	N-terminal 10xHis-tagged Active	CSB-MP013481HU(F8)
Mbp	Mouse	Yeast	N-terminal 6xHis-tagged	CSB-YP013551MO
MS4A1	Human	Yeast	N-terminal 6xHis-tagged	CSB-YP015007HU1
Ms4a1	Mouse	Yeast	N-terminal 6xHis-tagged	CSB-YP015007MO1
Ms4a1	Mouse	Yeast	N-terminal 6xHis-tagged	CSB-YP015007MO2
MS4A1	Human	Mammalian cell	C-terminal 10xHis-tagged (This tag can be tested only under denaturing conditions) Active	CSB-MP015007HU
MS4A1	Dog	Mammalian cell	C-terminal 10xHis-tagged (This tag can be tested only under denaturing conditions) Active	CSB-MP661636DO
MS4A1	Human	in vitro E.coli expression system	N-terminal 10xHis-tagged	CSB-CF015007HU
MS4A1	Dog	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP661636DO1
MS4A1	Cynomolgus monkey	Mammalian cell	C-terminal 10xHis-tagged (This tag can be tested only under denaturing conditions) Active	CSB-MP4516MOV
Ms4a1	Rat	Mammalian cell	C-terminal 10xHis-tagged (This tag can be tested only under denaturing conditions)	CSB-MP5193RA

Target	Organism	Source	Tag Info	Code
Ms4a1	Mouse	Mammalian cell	C-terminal 10xHis-tagged (This tag can be tested only under denaturing conditions)	CSB-MP015007MO
MSTN	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP015057HU
MSTN	Cat	E.coli	Tag-Free	CSB-EP3992CA
MTOR	Human	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP008968HU
NFKB1	Human	E.coli	C-terminal 6xHis-tagged	CSB-EP015759HU
NFKB1	Human	E.coli	C-terminal 6xHis-tagged	CSB-EP015759HU1
NGF	Pig	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP643662PI
Ngf	Rat	E.coli	N-terminal 6xHis-tagged	CSB-EP015779RA1
NGF	Human	E.coli	Tag-Free Active	CSB-AP003771HU
NGF	Human	Mammalian cell	Tag-Free Active	CSB-AP003951HU
Ngf	Mouse	E.coli	Tag-Free Active	CSB-AP004081MO
NGF	Human	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP015779HU
OPTN	Human	Baculovirus	C-terminal 9xHis-tagged	CSB-BP016363HU
OPTN	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP016363HU
OTC	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP017270HU
Otc	Mouse	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP017270MO
PAH	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP017396HU
PARK7	Human	E.coli	N-terminal GST-tagged	CSB-EP860342HU1
PARK7	Human	E.coli	Tag-Free	CSB-EP860342HU1e1
PDCD1	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP619964HU
PDCD1	Human	Mammalian cell	C-terminal 6xHis-tagged Active	CSB-MP619964HU1
PLG	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP018188HU1
PLG	Human	Baculovirus	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-BP018188HU(M1)
PLG	Human	Baculovirus	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-BP018188HU(M2)
PLG	Human	E.Coli	Tag-Free Active	CSB-AP002261HU
Ppt1	Mouse	Mammalian cell	C-terminal 10xHis-tagged	CSB-MP018587MO
PRKN	Human	Baculovirus	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-BP017451HU
PRKN	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP017451HU


Target	Organism	Source	Tag Info	Code
Prkn	Mouse	E.coli	N-terminal GST-tagged	CSB-EP893911MO
Prkn	Rat	Baculovirus	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-BP864348RA
Pycard	Mouse	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP861664MO
PYCARD	Human	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP890936HU
PYCARD	Human	Mammalian cell	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-MP890936HU
Pycard	Mouse	E.coli	Tag-Free	CSB-EP861664MOe1
RGMA	Human	Baculovirus	N-terminal MBP-tagged and C-terminal 6xHis-tagged	CSB-BP836186HUd8
RGMA	Human	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP836186HU
SCN1A	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP020834HU
Sema4d	Mouse	Mammalian cell	C-terminal mFc-tagged Active	CSB-MP020990MO
SEMA4D	Human	Mammalian cell	C-terminal 10xHis-tagged Active	CSB-MP835707HU
SEMA4D	Human	Mammalian cell	C-terminal hFc-tagged Active	CSB-MP835707HUd9
SEMA4D	Rhesus macaque	Mammalian cell	C-terminal 10xHis-tagged Active	CSB-MP3947MOW
Sirt3	Rat	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP848095RA
Sirt3	Mouse	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP848095MO
SIRT3	Human	E.coli	C-terminal 6xHis-tagged	CSB-EP882102HU1
SIRT3	Human	E.coli	C-terminal 6xHis-tagged	CSB-EP882102HU1(M)
SIRT3	Human	E.coli	C-terminal 6xHis-tagged	CSB-EP882102HU2(M)
Smn1	Mouse	Baculovirus	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-BP021838MO
Smn1	Mouse	Mammalian cell	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-MP021838MO
SNAP25	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP021873HU
SNCA	Human	E.coli	N-terminal GST-tagged	CSB-EP021912HU1
Snca	Mouse	E.coli	Tag-Free	CSB-EP021912MO
SNCA	Cynomolgus monkey	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP021912MOV
Snca	Rat	E.coli	N-terminal GST-tagged	CSB-EP021912RAe0
SNCA	Pig	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-EP666915PI
SNCA	Human	Yeast	N-terminal 6xHis-tagged	CSB-YP021912HU
SNCA	Cynomolgus monkey	Yeast	N-terminal 6xHis-tagged	CSB-YP021912MOV

Target	Organism	Source	Tag Info	Code
sod1	Zebrafish	E.coli	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged	CSB-EP022397DIL
Sod1	Mouse	E.coli	N-terminal 6xHis-tagged	CSB-EP022397MO
Sod1	Mouse	Yeast	N-terminal 6xHis-tagged	CSB-YP022397MO
SOD1	Human	E.Coli	N-terminal 10xHis-tagged Active	CSB-AP000121HU
SOD1	Human	E.Coli	Tag-Free Active	CSB-AP000131HU
Sod1	Rat	Yeast	N-terminal 6xHis-tagged	CSB-YP022397RA
Sort1	Rat	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP022412RA
Sort1	Rat	E.coli	N-terminal 6xHis-tagged	CSB-EP022412RAa0
TARDBP	Human	Yeast	N-terminal 6xHis-tagged	CSB-YP023129HU
TARDBP	Human	E.coli	Tag-Free	CSB-EP023129HU(A4)
Tfrc	Rat	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP858708RA1
TFRC	Human	Mammalian cell	N-terminal 10xHis-tagged Active	CSB-MP3648HU
TGFB1	Bovine	E.coli	N-terminal 6xHis-tagged	CSB-EP023446BO
TGFB1	Chicken	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP023446CH
Tgfb1	Rat	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP023446RA
Tgfb1	Rat	E.coli	N-terminal 6xHis-tagged	CSB-EP023446RAa0
Tgfb1	Rat	E.coli	Tag-Free	CSB-EP023446RAe1
Tgfb1	Rat	Mammalian cell	N-terminal 6xHis-tagged	CSB-MP023446RA
TGFB1	Human	Mammalian cell	Tag-Free Active	CSB-AP003861HU
Tgfb1	Mouse	Mammalian cell	Tag-Free Active	CSB-AP004111MO
TLR2	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP023601HU1
TLR4	Human	E.coli	N-terminal 10xHis-tagged	CSB-EP023603HU1
TNF	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP023955HU
TNF	Human	E.coli	Tag-Free Active	CSB-AP004861HU
TNF	Human	E.coli	C-terminal 6xHis-tagged Active	CSB-AP004871HU
TNF	Human	E.coli	N-terminal 6xHis-tagged Active	CSB-AP004881HU
Tnf	Mouse	E.coli	Tag-Free Active	CSB-AP005031MO
TNF	Rabbit	E.coli	Tag-Free Active	CSB-AP005051RB


Target	Organism	Source	Tag Info	Code
TNF	Human	Yeast	N-terminal 6xHis-tagged	CSB-EP022397DIL
TNF	Human	E.coli	C-terminal 6xHis-tagged	CSB-EP022397MO
TNFSF9	Human	E.Coli	Tag-Free	CSB-YP022397MO
TNFSF9	Human	Mammalian cell	N-terminal hFc-Myc-tagged	CSB-AP000121HU
TNFSF9	Human	Mammalian cell	C-terminal hFc-Avi-tagged	CSB-AP000131HU
TP53	Human	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-YP022397RA
TP53	Mouse	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP022412RA
TP53	Rat	E.coli	N-terminal 6xHis-tagged	CSB-EP022412RAa0
TP53	Human	Baculovirus	C-terminal 6xHis-tagged	CSB-YP023129HU
TP53	Human	Baculovirus	C-terminal 6xHis-tagged	CSB-EP023129HU(A4)
TP53	Human	Baculovirus	C-terminal 6xHis-tagged	CSB-EP858708RA1
TREM2	Human	E.coli	N-terminal 6xHis-tagged	CSB-MP3648HU
Trem2	Mouse	E.coli	N-terminal 6xHis-tagged	CSB-EP023446BO
Trem2	Mouse	Mammalian cell	N-terminal 6xHis-Myc-tagged	CSB-EP023446CH
TREM2	Human	E.coli	N-terminal 10xHis-GST-tagged and C-terminal Myc-tagged	CSB-EP023446RA
TTR	Human	E.coli	N-terminal 6xHis-tagged	CSB-EP023446RAa0
Ttr	Mouse	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-EP023446RAe1
Ttr	Mouse	E.coli	N-terminal GST-tagged	CSB-MP023446RA
Ttr	Mouse	E.coli	Tag-Free	CSB-AP003861HU
TTR	Cynomolgus monkey	E.coli	N-terminal 6xHis-tagged	CSB-AP004111MO
TTR	Cynomolgus monkey	Mammalian cell	N-terminal 6xHis-Myc-tagged	CSB-EP023601HU1
Ttr	Mouse	Yeast	N-terminal 6xHis-tagged	CSB-EP023603HU1
TTR	Human	Mammalian cell	C-terminal hFc-Myc-tagged Active	CSB-EP023955HU
Ttr	Mouse	Mammalian cell	C-terminal 10xHis-tagged Active	CSB-AP004861HU
VEGFA	Pig	E.coli	N-terminal 6xHis-tagged	CSB-AP004871HU
VEGFA	Rabbit	E.coli	N-terminal 6xHis-SUMO-tagged	CSB-AP004881HU
VEGFA	Rabbit	Yeast	N-terminal 6xHis-tagged	CSB-AP005031MO
Vwf	Mouse	E.coli	N-terminal 10xHis-tagged and C-terminal Myc-tagged	CSB-AP005051RB

● Antibodies


Phospho-AKT1 (Ser473) Recombinant Monoclonal Antibody
CSB-RA001553A473phHU




EGFR Recombinant Monoclonal Antibody
CSB-RA159341A0HU



GFAP Recombinant Monoclonal Antibody
CSB-RA009369MA1HU



VCP Recombinant Monoclonal Antibody
CSB-RA182227A0HU



Recombinant Antibody

Product Name	Code	Target	Species	Application
ACVR2B Recombinant Monoclonal Antibody	CSB-RA260702A0HU	ACVR2B	Human	ELISA, IHC
Phospho-AKT1 (T450) Recombinant Monoclonal Antibody	CSB-RA001553A450phHU	AKT1	Human	ELISA, WB
Phospho-AKT1 (Ser473) Recombinant Monoclonal Antibody	CSB-RA001553A473phHU	AKT1	Human	ELISA, WB, IHC, IF, IP
AKT1 Recombinant Monoclonal Antibody	CSB-RA917625A0HU	AKT1	Human, Mouse, Rat	ELISA, WB, IHC, IF
ALB Recombinant Monoclonal Antibody	CSB-RA264109A0HU	ALB	Human	ELISA, FC
ANGPT2 Recombinant Monoclonal Antibody	CSB-RA191985A0HU	ANGPT2	Human	ELISA, IHC
ANGPT2 Recombinant Monoclonal Antibody	CSB-RA001707MA01HU	ANGPT2	Human	ELISA
APP Recombinant Monoclonal Antibody	CSB-RA994273A0HU	APP	Human, Mouse, Rat	ELISA, WB, IF
BACE1 Recombinant Monoclonal Antibody	CSB-RA177574A0HU	BACE1	Human, Mouse, Rat	ELISA, WB
BCHE Recombinant Monoclonal Antibody	CSB-RA252650A0HU	BCHE	Human	ELISA, WB
BDNF Recombinant Monoclonal Antibody	CSB-RA989910A0HU	BDNF	Human, Mouse	ELISA, WB
C3 Recombinant Monoclonal Antibody	CSB-RA303909A0HU	C3	Human	ELISA, WB, IHC, FC
C3 Recombinant Monoclonal Antibody	CSB-RA699863A0HU	C3	Human	ELISA, IHC
C5 Recombinant Monoclonal Antibody	CSB-RA003995MA1HU	C5	Human	ELISA, FC
CD19 Recombinant Monoclonal Antibody	CSB-RA004888A0HU	CD19	Human	ELISA
CD19 Recombinant Monoclonal Antibody	CSB-RA780821A0HU	CD19	Human	ELISA, WB, IHC, IP
CD19 Recombinant Monoclonal Antibody	CSB-RA004888MA1HU	CD19	Human, Mouse	ELISA, WB
CD247 Recombinant Monoclonal Antibody	CSB-RA244537A0HU	CD247	Human	ELISA, FC
CD274 Recombinant Monoclonal Antibody	CSB-RA977797A0HU	CD274	Human	ELISA, WB, IHC
CD274 Recombinant Monoclonal Antibody	CSB-RA878942MA1HU	CD274	Human	ELISA, IHC, FC
CD38 Recombinant Monoclonal Antibody	CSB-RA796695A0HU	CD38	Human	ELISA, WB, IHC
CD3E Recombinant Monoclonal Antibody	CSB-RA004931MA1HU	CD3E	Human	ELISA, FC
CD4 Recombinant Monoclonal Antibody	CSB-RA004935A0HU	CD4	Human	ELISA, WB, IHC
CD40 Recombinant Monoclonal Antibody	CSB-RA004936MA1HU	CD40	Human	ELISA, IF, FC
CD47 Recombinant Monoclonal Antibody	CSB-RA802124A0HU	CD47	Human	ELISA, WB, IHC, IF
CD80 Recombinant Monoclonal Antibody	CSB-RA246383A0HU	CD80	Human	ELISA, WB, IHC, IP
CD8A Recombinant Monoclonal Antibody	CSB-RA934881A0HU	CD8A	Human	ELISA, WB, IHC, IP
CTLA4 Recombinant Monoclonal Antibody	CSB-RA213310A0HU	CTLA4	Human	ELISA, IHC

Product Name	Code	Target	Species	Application
CTLA4 Recombinant Monoclonal Antibody	CSB-RA006163MA1HU	CTLA4	Human, Mouse	ELISA, IHC
DLG4 Recombinant Monoclonal Antibody	CSB-RA255792A0HU	DLG4	Human, Mouse, Rat	ELISA, WB
DMD Recombinant Monoclonal Antibody	CSB-RA271553A0HU	DMD	Human	ELISA, WB, IHC, IF, IP
Phospho-EGFR (Y1092) Recombinant Monoclonal Antibody	CSB-RA007479A1092phHU	EGFR	Human	ELISA, WB, IHC, IF
Phospho-EGFR (Y1068) Recombinant Monoclonal Antibody	CSB-RA007479A1068phHU	EGFR	Human	ELISA, FC
EGFR Recombinant Monoclonal Antibody	CSB-RA159341A0HU	EGFR	Human	ELISA, IHC
EGFR Recombinant Monoclonal Antibody	CSB-RA794061A0HU	EGFR	Human	ELISA
EGFR Recombinant Monoclonal Antibody	CSB-RA159341MA1HU	EGFR	Human	ELISA, WB, IF
ERBB2 Recombinant Monoclonal Antibody	CSB-RA260392A0HU	ERBB2	Human	ELISA, WB
ERBB2 Recombinant Monoclonal Antibody	CSB-RA007763MA1HU	ERBB2	Human	ELISA, WB
Phospho-ERBB2 (Y1139) Recombinant Monoclonal Antibody	CSB-RA588766A0HU	ERBB2	Human	ELISA, WB
ERBB2 Recombinant Monoclonal Antibody	CSB-RA007763MA2HU	ERBB2	Human	ELISA, WB, IHC, FC
F10 Recombinant Monoclonal Antibody	CSB-RA599469A0HU	F10	Human	ELISA, IHC
FN1 Recombinant Monoclonal Antibody	CSB-RA921689A0HU	FN1	Human	ELISA, FC
FN1 Recombinant Monoclonal Antibody	CSB-RA919940A0HU	FN1	Human	ELISA
FN1 Recombinant Monoclonal Antibody	CSB-RA008759MA1HU	FN1	Human	ELISA, WB, IHC, IP
GAA Recombinant Monoclonal Antibody	CSB-RA566370A0HU	GAA	Human	ELISA, WB
GDNF Recombinant Monoclonal Antibody	CSB-RA009356A0HU	GDNF	Human	ELISA, FC
GFAP Recombinant Monoclonal Antibody	CSB-RA241995A0HU	GFAP	Human	ELISA, WB, IHC
GFAP Recombinant Monoclonal Antibody	CSB-RA009369MA1HU	GFAP	Human, Mouse, Rat	ELISA, IHC, FC
GHR Recombinant Monoclonal Antibody	CSB-RA978909A0HU	GHR	Human	ELISA, WB, IHC
GJB2 Recombinant Monoclonal Antibody	CSB-RA009452MA1HU	GJB2	Human	ELISA, FC
GLP1R Recombinant Monoclonal Antibody	CSB-RA009514MA1HU	GLP1R	Human	ELISA, WB, IHC
HTT Recombinant Monoclonal Antibody	CSB-RA802195A0HU	HTT	Human	ELISA, IF, FC
IFNAR1 Recombinant Monoclonal Antibody	CSB-RA011046A0HU	IFNAR1	Human	ELISA, WB, IHC, IF
IGF1R Recombinant Monoclonal Antibody	CSB-RA011087MA1HU	IGF1R	Human	ELISA, WB, IHC, IP
IL10 Recombinant Monoclonal Antibody	CSB-RA011580A0HU	IL10	Human	ELISA, WB, IHC, IP
IL12B Recombinant Monoclonal Antibody	CSB-RA895985A0HU	IL12B	Human	ELISA, IHC

Product Name	Code	Target	Species	Application
IL17A Recombinant Monoclonal Antibody	CSB-RA624104MA1HU	IL17A	Human	ELISA
IL1R1 Recombinant Monoclonal Antibody	CSB-RA011621MA1HU	IL1R1	Human	ELISA, FC
IL2RA Recombinant Monoclonal Antibody	CSB-RA011649MA1HU	IL2RA	Human	ELISA, IF
IL6 Recombinant Monoclonal Antibody	CSB-RA011664MA1HU	IL6	Human	ELISA
INS Recombinant Monoclonal Antibody	CSB-RA584163A0HU	INS	Human	ELISA, IHC
INSR Recombinant Monoclonal Antibody	CSB-RA155156A0HU	INSR	Human	ELISA, IF
KRT20 Recombinant Monoclonal Antibody	CSB-RA176322A0HU	KRT20	Human	ELISA, IHC, FC
LDLR Recombinant Monoclonal Antibody	CSB-RA575353A0HU	LDLR	Human	ELISA, IHC, IF, FC
LRRK2 Recombinant Monoclonal Antibody	CSB-RA797150A0HU	LRRK2	Human	ELISA, IHC
MAOB Recombinant Monoclonal Antibody	CSB-RA999377A0HU	MAOB	Human	ELISA, WB, IHC
Phospho-MAPT (T231) Recombinant Monoclonal Antibody	CSB-RA013481A231phHU	MAPT	Human	ELISA, WB
Phospho-MAPT (S396) Recombinant Monoclonal Antibody	CSB-RA013481A396phHU	MAPT	Human	ELISA, WB
Phospho-MAPT (S324) Recombinant Monoclonal Antibody	CSB-RA013481A324phHU	MAPT	Human	ELISA, IF
MAPT Recombinant Monoclonal Antibody	CSB-RA246354A0HU	MAPT	Human	ELISA, IHC
MAPT Recombinant Monoclonal Antibody	CSB-RA013481A1HU	MAPT	Mouse, Macaca mulatta	ELISA
Phospho-MAPT (S199) Recombinant Monoclonal Antibody	CSB-RA051594A0HU	MAPT	Human	ELISA, IHC, FC
Phospho-MAPT (S404) Recombinant Monoclonal Antibody	CSB-RA901354A0HU	MAPT	Human	ELISA, IHC
Phospho-MAPT (S214) Recombinant Monoclonal Antibody	CSB-RA050476A0HU	MAPT	Human	ELISA, IHC
Phospho-MAPT (S202) Recombinant Monoclonal Antibody	CSB-RA082445A0HU	MAPT	Human	ELISA, FC
MAPT Recombinant Monoclonal Antibody	CSB-RA013481MA1HU	MAPT	Human	ELISA, FC
CD20 Recombinant Monoclonal Antibody	CSB-RA015007A0HU	MS4A1	Human	ELISA
MS4A1 Recombinant Monoclonal Antibody	CSB-RA015007A1HU	MS4A1	Human	ELISA
MS4A1 Recombinant Monoclonal Antibody	CSB-RA015007MA3HU	MS4A1	Human, Dog, Macaca fascicularis	ELISA
MS4A1 Recombinant Monoclonal Antibody	CSB-RA015007MA4HU	MS4A1	Human	ELISA, IHC, IF
Phospho-MTOR (S2481) Recombinant Monoclonal Antibody	CSB-RA008968A2481phHU	MTOR	Human	ELISA, WB, IF
Phospho-MTOR (S2448) Recombinant Monoclonal Antibody	CSB-RA008968A2448phHU	MTOR	Human	ELISA, WB, IHC, IF
NGF Recombinant Monoclonal Antibody	CSB-RA288635A0HU	NGF	Human, Rat	ELISA, WB, IHC
NTRK2 Recombinant Monoclonal Antibody	CSB-RA619082MA1HU	NTRK2	Human	ELISA, IF

Product Name	Code	Target	Species	Application
PAH Recombinant Monoclonal Antibody	CSB-RA169000A0HU	PAH	Human	ELISA, IHC, IF
PCSK9 Recombinant Monoclonal Antibody	CSB-RA259644A0HU	PCSK9	Human	ELISA, WB, IHC
PDCD1 Recombinant Monoclonal Antibody	CSB-RA240597A0HU	PDCD1	Human	ELISA, WB, IHC
SEMA4D Recombinant Monoclonal Antibody	CSB-RA835707A1HU	SEMA4D	Human	ELISA, FC
SEMA4D Recombinant Monoclonal Antibody	CSB-RA835707A2HU	SEMA4D	Human, Mouse	ELISA, FC
SMN1 Recombinant Monoclonal Antibody	CSB-RA567382A0HU	SMN1	Human	ELISA, IHC
Phospho-SNCA (S129) Recombinant Monoclonal Antibody	CSB-RA021912A129phHU	SNCA	Human	ELISA, IHC, IF
SNCA Recombinant Monoclonal Antibody	CSB-RA585595A0HU	SNCA	Human	ELISA, IHC
SNCA Recombinant Monoclonal Antibody	CSB-RA021912MA1HU	SNCA	Human	ELISA, IHC, IF, FC
SOD1 Recombinant Monoclonal Antibody	CSB-RA829583A0HU	SOD1	Human	ELISA, WB, FC
SOD1 Recombinant Monoclonal Antibody	CSB-RA022397MA1HU	SOD1	Human	ELISA, IF, FC
TBK1 Recombinant Monoclonal Antibody	CSB-RA932866A0HU	TBK1	Human	ELISA, IHC, FC
TFRC Recombinant Monoclonal Antibody	CSB-RA023441MA1HU	TFRC	Human	ELISA
Phospho-TP53 (T55) Recombinant Monoclonal Antibody	CSB-RA024077A55phHU	TP53	Human	ELISA, WB, IF
Phospho-TP53 (S392) Recombinant Monoclonal Antibody	CSB-RA024077A392phHU	TP53	Human	ELISA, WB
Phospho-TP53 (S33) Recombinant Monoclonal Antibody	CSB-RA024077A33phHU	TP53	Human	ELISA, IHC, IF
Phospho-TP53 (S9) Recombinant Monoclonal Antibody	CSB-RA024077A09phHU	TP53	Human	ELISA, IF
TP53 Recombinant Monoclonal Antibody	CSB-RA236796A0HU	TP53	Human	ELISA, WB, IHC, IF
TP53 Recombinant Monoclonal Antibody	CSB-RA825742A0HU	TP53	Human	ELISA, WB, IHC, IF
TP53 Recombinant Monoclonal Antibody	CSB-RA214807A0HU	TP53	Human	ELISA, IHC
TP53 Recombinant Monoclonal Antibody	CSB-RA592348A0HU	TP53	Human	ELISA, IF
TP53 Recombinant Monoclonal Antibody	CSB-RA024077MA1HU	TP53	Human	ELISA, WB, IF, FC
TTR Recombinant Monoclonal Antibody	CSB-RA437590A0HU	TTR	Human	ELISA, IHC
TTR Recombinant Monoclonal Antibody	CSB-RA588381A0HU	TTR	Human	ELISA, IHC
VCP Recombinant Monoclonal Antibody	CSB-RA182227A0HU	VCP	Human, Rat	ELISA, WB, IHC, IF, IP

Monoclonal Antibody

Product Name	Code	Target	Species	Application
AKT1 Monoclonal Antibody	CSB-MA011037	AKT1	Human,Mouse,Pig,Rat	ELISA, WB
ALB Monoclonal Antibody	CSB-MA076752	ALB	Human	ELISA, WB, IHC
APP Monoclonal Antibody	CSB-MA0019501A0m	APP	Human, Mouse	ELISA, WB, IHC, IF
CD19 Monoclonal Antibody	CSB-MA004888A0m	CD19	Human	ELISA, WB, IHC, IF, FC
CD19 Monoclonal Antibody	CSB-MA004888A1m	CD19	Human	ELISA, WB, IHC, FC
CD19 Monoclonal Antibody,Purified	CSB-MA798261	CD19	Human	ELISA, FC
CD38 Monoclonal Antibody	CSB-MA785760	CD38	Human	ELISA, IHC
CD3E Monoclonal Antibody,Purified	CSB-MA947479	CD3E	Human	ELISA, IF, FC
CD4 Monoclonal Antibody	CSB-MA000228	CD4	Human,Mouse,Rat	ELISA, IHC
CD4 Monoclonal Antibody,Purified	CSB-MA783233	CD4	Human	ELISA, FC
CD8A Monoclonal Antibody,Purified	CSB-MA927984	CD8A	Human	ELISA, FC
CTLA4 Monoclonal Antibody	CSB-MA006163A0m	CTLA4	Human	ELISA, WB, IHC
EGFR Monoclonal Antibody	CSB-MA000199	EGFR	Human	ELISA, WB, IHC, IF
EGFR Monoclonal Antibody	CSB-MA000293	EGFR	Human	ELISA, WB
ERBB2 Monoclonal Antibody	CSB-MA000200	ERBB2	Human,Mouse,Rat	ELISA, WB, IHC
ERBB2 Monoclonal Antibody	CSB-MA000237	ERBB2	Human	ELISA, WB
ERBB2 Monoclonal Antibody	CSB-MA172796	ERBB2	Human	ELISA, WB, IHC
FN1 Monoclonal Antibody	CSB-MA000296	FN1	Human,Mouse,Rat	ELISA, WB
FN1 Monoclonal Antibody	CSB-MA000318	FN1	Mouse,Rat,Human	ELISA, WB
FOXP3 Monoclonal Antibody	CSB-MA545218	FOXP3	Human	ELISA, WB
FOXP3 Monoclonal Antibody	CSB-MA925351	FOXP3	Human	ELISA, WB, ICC
GFAP Monoclonal Antibody	CSB-MA009369A0m	GFAP	Human, Mouse, Rat	ELISA, WB, IHC, IF, FC
GFAP Monoclonal Antibody	CSB-MA000216	GFAP	Rat,Mouse	ELISA, WB
GFAP Monoclonal Antibody	CSB-MA293354	GFAP	Human	ELISA, IHC
IFNG Monoclonal Antibody	CSB-MA065241A0m	IFNG	Human	ELISA, IHC
IL1B Monoclonal Antibody	CSB-MA065741A0m	IL1B	Human	ELISA
IL1B Monoclonal Antibody	CSB-MA010884	IL1B	Human	ELISA, WB, IHC, ICC
IL2RA Monoclonal Antibody	CSB-MA000266	IL2RA	Human	ELISA, IHC

Product Name	Code	Target	Species	Application
IL6 Monoclonal Antibody	CSB-MA067571A0m	IL6	Human	ELISA, WB, IHC
MBP Monoclonal Antibody	CSB-MA000061M0m	MBP	All	ELISA, WB
MS4A1 Monoclonal Antibody	CSB-MA000204	MS4A1	Human,Mouse,Rat	ELISA, IHC
MS4A1 Monoclonal Antibody	CSB-MA996413	MS4A1	Human	ELISA, IHC
SEMA4D Monoclonal Antibody	CSB-MA835707A0m	SEMA4D	Human	ELISA, WB, FC
TNF Monoclonal Antibody	CSB-MA084771A0m	TNF	Bovine	ELISA
TNF Monoclonal Antibody	CSB-MA080271	TNF	Human,Mouse,Rat	ELISA, WB, IHC
TNF Monoclonal Antibody	CSB-MA080272	TNF	Human,Mouse,Rat	ELISA, WB, IHC
TP53 Monoclonal Antibody	CSB-MA0240771A0m	TP53	Human	ELISA, WB, IHC
TP53 Monoclonal Antibody	CSB-MA000208	TP53	Human	ELISA, WB, IHC, IF
TP53 Monoclonal Antibody	CSB-MA989534	TP53	Human	ELISA, WB
TP53 Monoclonal Antibody	CSB-MA189580	TP53	Human,Dog,Mouse,Rat	ELISA, WB
TTR Monoclonal Antibody	CSB-MA263346	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA633769	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA116768	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA287792	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA111203	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA382290	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA306092	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA110724	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA351513	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA126741	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA681670	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA362195	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA118621	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA831325	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA336957	TTR	Human	ELISA, WB, IHC
TTR Monoclonal Antibody	CSB-MA447738	TTR	Human	ELISA, WB, IHC

Polyclonal Antibody

Product Name	Code	Target	Species	Application
ABCD1 Antibody	CSB-PA001068LA01HU	ABCD1	Human	ELISA, IHC
ABCD1 Antibody	CSB-PA008163	ABCD1	Human	ELISA, WB
ABCD1 Antibody	CSB-PA001068GA01HU	ABCD1	Human, Mouse, Rat	ELISA, WB, IHC
ABCD1 Antibody	CSB-PA214656	ABCD1	Human	ELISA, WB
ACHE Antibody	CSB-PA5619A0Rb	ACHE	Human	ELISA, IHC, IF
ache Antibody	CSB-PA863598LA01DIL	ache	Zebrafish	ELISA, WB
ACHE Antibody	CSB-PA010202	ACHE	Human, Mouse, Rat	ELISA, WB, IHC
ACHE Antibody	CSB-PA010203	ACHE	Human, Mouse, Rat	ELISA, WB
ACHE Antibody	CSB-PA001154GA01HU	ACHE	Human, Mouse, Rat	ELISA, WB
ACHE Antibody	CSB-PA904467	ACHE	Human, Mouse, Rat	ELISA, WB
ACHE Antibody	CSB-PA114463	ACHE	Human, Mouse, Rat	ELISA, WB, IHC
ACHE Antibody	CSB-PA547258	ACHE	Human	ELISA, WB
ACVR2B Antibody	CSB-PA623829ESR1HU	ACVR2B	Human	ELISA, IHC
ACVR2B Antibody	CSB-PA292270	ACVR2B	Human, Mouse, Rat	ELISA, IHC
ACVR2B Antibody	CSB-PA206692	ACVR2B	Human, Mouse, Rat	ELISA, WB, IHC
ADCYAP1R1 Antibody	CSB-PA207627	ADCYAP1R1	Human, Mouse	ELISA, WB, IHC
AGER Antibody	CSB-PA14209A0Rb	AGER	Human	ELISA, IHC, IF
AGER Antibody	CSB-PA007810	AGER	Human	ELISA, WB, IHC
AGER Antibody	CSB-PA066696	AGER	Human	ELISA, WB, IHC
AGER Antibody	CSB-PA874417	AGER	Human	ELISA, WB, IHC
AKT1 Antibody	CSB-PA15905A0Rb	AKT1	Human	ELISA, WB, IHC, IF, IP
AKT1 Antibody	CSB-PA000852	AKT1	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-AKT1 (S124) Antibody	CSB-PA008118	AKT1	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-AKT1 (T72) Antibody	CSB-PA008120	AKT1	Human, Mouse, Rat	ELISA, WB, IHC
AKT1 Antibody	CSB-PA008122	AKT1	Human, Mouse, Rat	ELISA, IHC
AKT1 Antibody	CSB-PA080256	AKT1	Human, Mouse, Rat	WB
Phospho-AKT1 (S129) Antibody	CSB-PA000645	AKT1	Human, Mouse, Rat	ELISA, WB
Phospho-AKT1 (T308) Antibody	CSB-PA000657	AKT1	Human, Mouse, Rat	ELISA, WB

Product Name	Code	Target	Species	Application
Phospho-AKT1 (S473) Antibody	CSB-PA000732	AKT1	Human, Mouse, Rat	ELISA, WB, IHC, IF
AKT1 Antibody	CSB-PA000851	AKT1	Human, Mouse, Rat, Monkey	ELISA, WB, IHC, IF
Phospho-AKT1 (S246) Antibody	CSB-PA000465	AKT1	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-AKT1 (T450) Antibody	CSB-PA000468	AKT1	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-AKT1 (Y326) Antibody	CSB-PA000649	AKT1	Human, Mouse, Rat	ELISA, WB
Phospho-AKT1 (Y474) Antibody	CSB-PA000704	AKT1	Human, Mouse, Rat, Monkey	ELISA, WB, IHC
AKT1 Antibody	CSB-PA000855	AKT1	Human, Mouse, Rat	ELISA, WB, IHC
AKT1 Antibody	CSB-PA000849	AKT1	Human, Mouse, Rat	ELISA, WB, IHC
AKT1 Antibody	CSB-PA001553GA01HU	AKT1	Human	ELISA, WB
AKT1 Antibody	CSB-PA098004	AKT1	Human, Mouse, Rat	ELISA, IHC
AKT1 Antibody	CSB-PA721992	AKT1	Human, Mouse, Rat	ELISA, WB
AKT1 Antibody	CSB-PA942120	AKT1	Human, Mouse, Rat	ELISA, WB
AKT1 (Ab-124) Antibody	CSB-PA919481	AKT1	Human, Mouse, Rat	ELISA, WB
AKT1 Antibody	CSB-PA149203	AKT1	Human, Mouse, Rat	ELISA, WB
AKT1 (Ab-473) Antibody	CSB-PA063905	AKT1	Human, Mouse, Rat	ELISA, WB, IHC
AKT1 (Ab-308) Antibody	CSB-PA288499	AKT1	Human, Mouse, Rat	ELISA, WB, IHC
AKT1 Antibody	CSB-PA204523	AKT1	Human, Mouse, Rat	ELISA, IHC
AKT1 (Ab-450) Antibody	CSB-PA969941	AKT1	Human, Mouse, Rat	ELISA, WB, IHC, IF
AKT1 (Ab-129) Antibody	CSB-PA996941	AKT1	Human, Mouse, Rat	ELISA, WB, IHC, IF
AKT1 (Ab-326) Antibody	CSB-PA927050	AKT1	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-AKT1 (Ser129) Antibody	CSB-PA115404	AKT1	Human, Mouse, Rat	ELISA, WB
Phospho-AKT1 (Tyr326) Antibody	CSB-PA063650	AKT1	Human, Mouse, Rat	ELISA, WB
Phospho-AKT1 (Thr450) Antibody	CSB-PA237749	AKT1	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-AKT1 (Ser473) Antibody	CSB-PA573922	AKT1	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-AKT1 (Thr308) Antibody	CSB-PA010227	AKT1	Human, Mouse, Rat	ELISA, WB, IHC, IF
ALB Antibody	CSB-PA00260E1Rb	ALB	Canine	ELISA
ALB Antibody	CSB-PA00950E1Rb	ALB	Bovine	ELISA
ALB Antibody	CSB-PA001561LA01HU	ALB	Human, Mouse	ELISA, WB, IHC, IF

Product Name	Code	Target	Species	Application
Alb Antibody	CSB-PA001561LA01RA	ALB	Rat	ELISA, WB
ALB Antibody	CSB-PA001561GA01HU	ALB	Human, Mouse, Rat	ELISA, IHC
ALB Antibody	CSB-PA079344	ALB	Human	ELISA, WB, IHC
BSA Antibody	CSB-PA592874	ALB	Bovine	ELISA, WB
BSA Antibody	CSB-PA272395	ALB	Bovine	ELISA, WB
ANGPT2 Antibody	CSB-PA05749A0Rb	ANGPT2	Human, Mouse, Rat	ELISA, WB, IHC
Angpt2 Antibody	CSB-PA16959A0Rb	Angpt2	Mouse	ELISA, WB, IHC
ANGPT2 Antibody	CSB-PA083784	ANGPT2	Human, Mouse, Rat	ELISA, WB, IHC
ANGPT2 Antibody	CSB-PA694467	ANGPT2	Human, Mouse, Rat	ELISA, IHC
APOE Antibody	CSB-PA07599A0Rb	APOE	Human	ELISA, IHC, IF
APOE Antibody	CSB-PA009485	APOE	Human, Mouse, Rat	ELISA, WB, IHC
APOE Antibody	CSB-PA009486	APOE	Human	ELISA, WB
APOE Antibody	CSB-PA227407	APOE	Human	ELISA, WB
APP Antibody	CSB-PA07174A0Rb	APP	Human	ELISA, IHC
APP Antibody	CSB-PA07179A0Rb	APP	Human, Mouse	ELISA, WB, IHC, IF
Phospho-APP (T743) Antibody	CSB-PA000724	APP	Human, Mouse, Rat	ELISA, WB, IHC, IF
APP Antibody	CSB-PA000875	APP	Human, Mouse, Rat	ELISA, WB, IHC, IF
APP Antibody	CSB-PA001950GA01HU	APP	Human, Mouse, Rat	ELISA, WB
APP Antibody	CSB-PA891395	APP	Human	ELISA, IHC
APP Antibody	CSB-PA123968	APP	Human, Mouse, Rat	ELISA, WB, IHC
APP (Ab-668) Antibody	CSB-PA907645	APP	Human, Mouse, Rat	ELISA, WB, IF
Phospho-APP (Thr668) Antibody	CSB-PA296382	APP	Human, Mouse, Rat	ELISA, WB
Phospho-APP (Thr743/668) Antibody	CSB-PA212801	APP	Human, Mouse, Rat	ELISA, WB, IHC, IF
ARSA Antibody	CSB-PA002141ESR1HU	ARSA	Human	ELISA
ARSA Antibody	CSB-PA009975	ARSA	Human, Mouse, Rat	ELISA, WB, IHC
ARSA Antibody	CSB-PA002141GA01HU	ARSA	Human, Mouse, Rat	ELISA, WB, IHC
ARSA Antibody	CSB-PA514644	ARSA	Human, Mouse	ELISA, IHC
ARSA Antibody	CSB-PA901669	ARSA	Human, Mouse, Rat	ELISA, WB, IHC

Product Name	Code	Target	Species	Application
ATXN2 Antibody	CSB-PA859520LA01HU	ATXN2	Human	ELISA, WB, IF
ATXN2 Antibody	CSB-PA000944	ATXN2	Human	ELISA, WB, IHC
ATXN2 Antibody	CSB-PA002441GA01HU	ATXN2	Human, Mouse, Rat	ELISA, WB
BACE1 Antibody	CSB-PA002524LA01HU	BACE1	Human	ELISA, IHC, IF
BACE1 Antibody	CSB-PA008418	BACE1	Human, Mouse, Rat	ELISA, IHC, IF
BACE1 Antibody	CSB-PA002524GA01HU	BACE1	Human, Mouse, Rat	ELISA, WB
BACE1 Antibody	CSB-PA173556	BACE1	Human, Mouse, Rat	ELISA, IHC
BACE1 Antibody	CSB-PA150463	BACE1	Human, Mouse, Rat	ELISA, WB
BCHE Antibody	CSB-PA002604LA01DO	BCHE	Dog	ELISA, IHC
BCHE Antibody	CSB-PA002604LA01HU	BCHE	Human	ELISA, WB, IF
BDNF Antibody	CSB-PA002655LA01HU	BDNF	Human	ELISA, IHC
BDNF Antibody	CSB-PA05775A0Rb	BDNF	Human, Mouse	ELISA, WB, IHC
Bdnf Antibody	CSB-PA16639A0Rb	Bdnf	Rat, Human	ELISA, WB
BDNF Antibody	CSB-PA002655GA01HU	BDNF	Human, Mouse, Rat	ELISA
BDNF Antibody	CSB-PA616184	BDNF	Human, Mouse, Rat	ELISA, WB, IHC
BDNF Antibody	CSB-PA988388	BDNF	Human, Mouse, Rat	ELISA, IHC
C1R Antibody	CSB-PA003655LA01HU	C1R	Human, Rat	ELISA, WB, IHC, IF
Cleaved-C1R (I464) Antibody	CSB-PA000016	C1R	Human	ELISA, WB
Cleaved-C1R (R463) Antibody	CSB-PA050282	C1R	Human	ELISA, WB
C1S Antibody	CSB-PA003657LA01HU	C1S	Human	ELISA, IHC
Cleaved-C1S (R437) Antibody	CSB-PA000017	C1S	Human	ELISA, WB
C1S Antibody	CSB-PA060002	C1S	Human	ELISA, WB, IHC
C1S Antibody	CSB-PA003657GA01HU	C1S	Human, Mouse, Rat	ELISA, WB, IHC
C1S Antibody	CSB-PA145542	C1S	Human	ELISA, WB, IHC
C3 Antibody	CSB-PA10599A0Rb	C3	Human	ELISA, IHC, IF
C3 Antibody	CSB-PA365507LA01MO	C3	Mouse, Human	ELISA, WB
C3 Antibody	CSB-PA360475HA01RA	C3	Rat	ELISA
C5 Antibody	CSB-PA003995LA01HU	C5	Human	ELISA, IHC

Product Name	Code	Target	Species	Application
C5 Antibody	CSB-PA003995YA01PI	C5	Pig	ELISA, IHC
C9orf72 Antibody	CSB-PA842710LA01HU	C9orf72	Human	ELISA, IHC, IF
CD19 Antibody	CSB-PA001431	CD19	Human, Mouse, Monkey	ELISA, WB
Phospho-CD19 (Y531) Antibody	CSB-PA070286	CD19	Human, Mouse, Monkey	ELISA, WB
CD19 Antibody	CSB-PA006159	CD19	Human	ELISA, WB
CD19 (Ab-531) Antibody	CSB-PA253616	CD19	Human, Mouse	ELISA, WB
Phospho-CD19 (Tyr531) Antibody	CSB-PA058389	CD19	Human, Mouse	ELISA, WB
CD247 Antibody	CSB-PA10499A0Rb	CD247	Human	ELISA, WB, IHC, IF, IP
Phospho-CD247 (Y142) Antibody	CSB-PA070253	CD247	Human, Mouse, Rat, Monkey	ELISA, WB, IHC
CD247 Antibody	CSB-PA070254	CD247	Human, Mouse, Rat, Monkey	ELISA, WB, IHC
CD247 Antibody	CSB-PA273375	CD247	Human	ELISA, IHC
CD247 Antibody	CSB-PA004904GA01HU	CD247	Human, Mouse, Rat	ELISA, WB, IHC
CD274 Antibody	CSB-PA878942ESR1HU	CD274	Human	ELISA, IHC
CD274 Antibody	CSB-PA004911GA01HU	CD274	Human, Mouse, Rat	ELISA, WB, IHC
CD276 Antibody	CSB-PA745871ESR1HU	CD276	Human	ELISA, WB, IHC
CD276 Antibody	CSB-PA745871ESR2HU	CD276	Human	ELISA, IHC
CD276 Antibody	CSB-PA005320	CD276	Human, Mouse, Rat	ELISA, WB
CD276 Antibody	CSB-PA004912GA01HU	CD276	Human, Mouse, Rat	ELISA, WB, IHC
CD33 Antibody	CSB-PA004925ESR1HU	CD33	Human	ELISA, WB, IHC
CD33 Antibody	CSB-PA004925ESR2HU	CD33	Human	ELISA, IHC
CD33 Antibody	CSB-PA004925LA01HU	CD33	Human	ELISA, WB, IHC
CD33 Antibody	CSB-PA006175	CD33	Human	ELISA, WB
CD33 Antibody	CSB-PA004925GA01HU	CD33	Human	ELISA, WB, IHC, IF
CD33 Antibody	CSB-PA070366	CD33	Human, Mouse	ELISA, WB, IHC
CD33 Antibody	CSB-PA243594	CD33	Human	ELISA, WB, IHC
CD38 Antibody	CSB-PA004929ESR1HU	CD38	Human	ELISA, IHC
CD38 Antibody	CSB-PA005628	CD38	Human	ELISA, WB
CD38 Antibody	CSB-PA134528	Cd38	Human	ELISA, WB

Product Name	Code	Target	Species	Application
CD38 Antibody	CSB-PA926996	CD38	Human	ELISA, WB, IHC
CD3E Antibody	CSB-PA10484A0Rb	CD3E	Human	ELISA, WB, IHC, IF
CD3E Antibody	CSB-PA004931HA01HU	CD3E	Human	ELISA, IHC, IF
CD3E Antibody	CSB-PA005833	CD3E	Human	ELISA, WB
CD3E Antibody	CSB-PA070244	CD3E	Human	ELISA, WB, IHC
CD3E Antibody	CSB-PA004931GA01HU	CD3E	Human	ELISA, WB, IHC
CD3E Antibody	CSB-PA949010	CD3E	Human, Mouse, Rat	ELISA, WB
CD4 Antibody	CSB-PA004935ESR2HU	CD4	Human	ELISA, IHC
Phospho-CD4 (S433) Antibody	CSB-PA070257	CD4	Human, Mouse	ELISA, IHC, IF
CD4 Antibody	CSB-PA001453	CD4	Human, Mouse	ELISA, WB, IHC, IF
CD4 Antibody	CSB-PA004935GA01HU	CD4	Human, Mouse, Rat	ELISA, WB
CD40 Antibody	CSB-PA004936LA01HU	CD40	Human, Mouse	ELISA, WB, IF
CD40 Antibody	CSB-PA001455	CD40	Human, Mouse, Monkey	ELISA, WB, IHC
Phospho-CD40 (T254) Antibody	CSB-PA080031	CD40	Human, Mouse	ELISA, WB
CD40 Antibody	CSB-PA004936GA01HU	CD40	Human	ELISA, WB
CD40 Antibody	CSB-PA215588	CD40	Human, Mouse	ELISA, IHC
Phospho-CD40 (Thr254) Antibody	CSB-PA176134	CD40	Human, Mouse	ELISA, WB
CD40LG Antibody	CSB-PA06005A0Rb	CD40LG	Human	ELISA
CD40LG Antibody	CSB-PA004937GA01HU	CD40LG	Human, Mouse, Rat	ELISA, WB, IHC
CD40LG Antibody	CSB-PA067336	CD40LG	Human, Mouse	ELISA, WB
CD47 Antibody	CSB-PA004940LA01HU	CD47	Human	ELISA, IHC, IF
CD47 Antibody	CSB-PA005993	CD47	Human	ELISA, WB
CD47 Antibody	CSB-PA004940GA01HU	CD47	Human	ELISA, WB
CD47 Antibody	CSB-PA445679	CD47	Human	ELISA, WB, IHC
CD47 Antibody	CSB-PA911004	CD47	Human	ELISA, IHC
CD52 Antibody	CSB-PA338154	CD52	Human	ELISA, IHC
CD80 Antibody	CSB-PA129040	CD80	Human	ELISA, IHC
CD80 Antibody	CSB-PA004959GA01HU	Cd80	Human	ELISA, WB

Product Name	Code	Target	Species	Application
CD80 Antibody	CSB-PA277650	CD80	Human	ELISA, IHC
CD80 Antibody	CSB-PA130815	CD80	Human	ELISA, IHC
CD8A Antibody	CSB-PA004966LA01HU	CD8A	Human	ELISA, IHC
CNGA3 Antibody	CSB-PA613670DSR1HU	CNGA3	Human	ELISA, IHC
CNGA3 Antibody	CSB-PA005644GA01HU	CNGA3	Human, Mouse, Rat	ELISA, WB
CNGA3 Antibody	CSB-PA292433	CNGA3	Human, Mouse	ELISA, WB
CNGA3 Antibody	CSB-PA185605	CNGA3	Human, Mouse	ELISA, WB
CNTFR Antibody	CSB-PA005684EA01HU	CNTFR	Human	ELISA, WB, IHC, IF
CNTFR Antibody	CSB-PA005684GA01HU	CNTFR	Human, Mouse, Rat	ELISA, WB
CNTFR Antibody	CSB-PA276708	CNTFR	Human, Mouse, Rat	ELISA, WB
CNTFR Antibody	CSB-PA204475	CNTFR	Human, Mouse, Rat	ELISA, WB
CRHR2 Antibody	CSB-PA623797LA01HU	CRHR2	Human, Rat, Mouse	ELISA, WB, IHC, IF
CRHR2 Antibody	CSB-PA007281	CRHR2	Human, Mouse, Rat	ELISA, WB, IHC, IF
CRHR2 Antibody	CSB-PA001788	CRHR2	Human, Mouse, Rat	ELISA, WB, IF
CRHR2 Antibody	CSB-PA350000	CRHR2	Human, Rat	ELISA, WB
CRHR2 Antibody	CSB-PA977299	CRHR2	Human, Mouse, Rat	ELISA, IHC
CRHR2 Antibody	CSB-PA557863	CRHR2	Human, Mouse, Rat	ELISA, WB, IHC
CSF2 Antibody	CSB-PA06475A0Rb	CSF2	Human	ELISA, IHC
CSF2 Antibody	CSB-PA006043GA01HU	CSF2	Human	ELISA, WB
CSF2 Antibody	CSB-PA006045GA01HU	CSF2	Human	ELISA, WB, IHC
CSF2 Antibody	CSB-PA573559	CSF2	Human, Mouse	ELISA, WB
CSF2 Antibody	CSB-PA268599	CSF2	Human, Mouse	ELISA, WB, IHC
CSF2RA Antibody	CSB-PA006046LA01HU	CSF2RA	Human	ELISA, IHC, IF
CSF2RA Antibody	CSB-PA005243	CSF2RA	Human	ELISA, WB
CSF2RA Antibody	CSB-PA006046GA01HU	CSF2RA	Human	ELISA, WB, IHC
CSF2RA Antibody	CSB-PA644750	CSF2RA	Human	ELISA, WB
CSF2RA Antibody	CSB-PA441042	CSF2RA	Human, Mouse	ELISA, WB
CSF2RA Antibody	CSB-PA193987	CSF2RA	Human	ELISA, WB, IHC

Product Name	Code	Target	Species	Application
CSF2RA Antibody	CSB-PA046967	CSF2RA	Human	ELISA, WB
CSF3R Antibody	CSB-PA860321LA01HU	CSF3R	Human	ELISA, IHC, IF
CSF3R Antibody	CSB-PA004693	CSF3R	Human	ELISA, WB
CSF3R Antibody	CSB-PA006057	CSF3R	Human	ELISA, WB
CSF3R Antibody	CSB-PA300322	CSF3R	Human	ELISA, IHC
CSF3R Antibody	CSB-PA094529	CSF3R	Human	ELISA, IHC
CTLA4 Antibody	CSB-PA006163LA01HU	CTLA4	Human	ELISA, IF
CTLA4 Antibody	CSB-PA252417	CTLA4	Human	ELISA, IHC
CTLA4 Antibody	CSB-PA006163KA01HU	CTLA4	Human	ELISA, WB, IHC
CTLA4 Antibody	CSB-PA224600	CTLA4	Human, Mouse, Rat	ELISA, IHC
DAG1 Antibody	CSB-PA613486NA01HU	DAG1	Human	ELISA, IHC
DAG1 Antibody	CSB-PA006488GA01HU	DAG1	Human, Mouse, Rat	ELISA, WB, IHC
DHPS Antibody	CSB-PA006854LA01HU	DHPS	Human	ELISA, IHC
DHPS Antibody	CSB-PA002145	DHPS	Human	ELISA, WB
DLG4 Antibody	CSB-PA006938LA01HU	DLG4	Human	ELISA, IHC
DLG4 Antibody	CSB-PA007805	DLG4	Human, Mouse, Rat	ELISA, WB
DLG4 Antibody	CSB-PA284504	DLG4	Human, Mouse, Rat	ELISA, IHC
DLG4 Antibody	CSB-PA575557	DLG4	Human, Mouse, Rat	ELISA, IHC
DMD Antibody	CSB-PA006963LA01HU	DMD	Human	ELISA, IHC
DMD Antibody	CSB-PA006963GA01HU	DMD	Human, Mouse, Rat	ELISA, WB
DMD Antibody	CSB-PA503789	DMD	Human, Mouse	ELISA, IHC
DMPK Antibody	CSB-PA006968LA01HU	DMPK	Human	ELISA, IHC, IF
DMPK Antibody	CSB-PA002164	DMPK	Human	ELISA, WB
DYSF Antibody	CSB-PA007314LA01HU	DYSF	Human	ELISA, IHC, IF
DYSF Antibody	CSB-PA060232	DYSF	Human, Mouse	ELISA, WB, IF
DYSF Antibody	CSB-PA184651	DYSF	Human, Mouse	ELISA, WB
EGFR Antibody	CSB-PA10279A0Rb	EGFR	Human	ELISA, IHC, IF
EGFR Antibody	CSB-PA005571LA01HU	EGFR	Human	ELISA, WB, IF

Product Name	Code	Target	Species	Application
Phospho-EGFR (T693) Antibody	CSB-PA000514	EGFR	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-EGFR (Y1092) Antibody	CSB-PA000515	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-EGFR (Y1110) Antibody	CSB-PA000516	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (T678) Antibody	CSB-PA000664	EGFR	Human, Mouse, Rat, Monkey	ELISA, WB, IF
EGFR Antibody	CSB-PA002264	EGFR	Human, Mouse, Rat	ELISA, WB, IHC, IF
EGFR Antibody	CSB-PA002273	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
EGFR Antibody	CSB-PA002275	EGFR	Human, Mouse	ELISA, WB
Phospho-EGFR (Y1172) Antibody	CSB-PA007951	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (Y869) Antibody	CSB-PA007952	EGFR	Human, Mouse, Rat	ELISA, WB, IF
Phospho-EGFR (Y1069) Antibody	CSB-PA007953	EGFR	Human, Mouse, Rat, Monkey	ELISA, WB
Phospho-EGFR (S695) Antibody	CSB-PA007954	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (Y1016) Antibody	CSB-PA007969	EGFR	Human, Mouse, Rat	ELISA, WB, IHC, IF
EGFR Antibody	CSB-PA007973	EGFR	Human, Mouse, Rat	ELISA, WB
EGFR Antibody	CSB-PA007976	EGFR	Human, Mouse, Rat	ELISA, WB
EGFR Antibody	CSB-PA007977	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-EGFR (S1070) Antibody	CSB-PA007950	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-EGFR (S1026) Antibody	CSB-PA007966	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (S1071) Antibody	CSB-PA007968	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (Y1197) Antibody	CSB-PA000517	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
EGFR Antibody	CSB-PA002267	EGFR	Human, Mouse, Rat	ELISA, WB, IHC, IF
EGFR Antibody	CSB-PA002281	EGFR	Human, Mouse, Rat	ELISA, WB, IHC, IF
EGFR Antibody	CSB-PA007971	EGFR	Human, Mouse, Rat	ELISA, WB
EGFR Antibody	CSB-PA007978	EGFR	Human, Mouse, Rat	ELISA, WB
EGFR Antibody	CSB-PA007479GA01HU	EGFR	Human	ELISA, WB
EGFR Antibody	CSB-PA162597	EGFR	Human, Mouse	ELISA, IHC
EGFR (Ab-678) Antibody	CSB-PA929356	EGFR	Human, Mouse, Rat	ELISA, WB
EGFR (Ab-693) Antibody	CSB-PA179555	EGFR	Human, Mouse, Rat	ELISA, WB
EGFR (Ab-1172) Antibody	CSB-PA033139	EGFR	Human, Mouse, Rat	ELISA, WB

Product Name	Code	Target	Species	Application
EGFR (Ab-1197) Antibody	CSB-PA911813	EGFR	Human, Mouse, Rat	ELISA, WB
EGFR (Ab-1197) Goat Polyclonal Antibody	CSB-PA956575	EGFR	Human, Mouse, Rat	ELISA, WB
EGFR (Ab-1071) Antibody	CSB-PA294390	EGFR	Human	ELISA, WB
EGFR (Ab-1092) Antibody	CSB-PA953009	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
EGFR (Ab-869) Antibody	CSB-PA877124	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
EGFR (Ab-1110) Antibody	CSB-PA111241	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
EGFR (Ab-1026) Antibody	CSB-PA079221	EGFR	Human	ELISA, WB
EGFR (Ab-1070) Antibody	CSB-PA174357	EGFR	Human, Mouse, Rat	ELISA, WB, IHC, IF
EGFR (Ab-1172) Antibody	CSB-PA579851	EGFR	Human, Mouse	ELISA, WB, IHC, IF
Phospho-EGFR (Tyr869) Antibody	CSB-PA184845	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (Tyr1110) Antibody	CSB-PA938612	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (Tyr1197) Antibody	CSB-PA183869	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (Tyr1016) Antibody	CSB-PA559519	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (Ser695) Antibody	CSB-PA189998	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (Ser1071) Antibody	CSB-PA119181	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (Ser1026) Antibody	CSB-PA957122	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (Ser1070) Antibody	CSB-PA209840	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-EGFR (Tyr1092) Antibody	CSB-PA069395	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-EGFR (Thr678) Antibody	CSB-PA561961	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-EGFR (Thr693) Antibody	CSB-PA979756	EGFR	Human, Mouse, Rat	ELISA, WB
Phospho-EGFR (Tyr1172) Antibody	CSB-PA277676	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-EGFR (Tyr1197) Antibody	CSB-PA049634	EGFR	Human, Mouse, Rat	ELISA, WB, IHC
EPHA4 Antibody	CSB-PA007724LA01HU	EPHA4	Human	ELISA, IHC, IF
Phospho-EPHA4 (Y596) Antibody	CSB-PA008062	EPHA4	Human, Mouse	ELISA, WB
EPHA4 Antibody	CSB-PA830643	EPHA4	Human, Mouse	ELISA, IHC
Phospho-EPHA4 (Tyr596) Antibody	CSB-PA959570	EPHA4	Human, Mouse	ELISA, WB
EPOR Antibody	CSB-PA007744LA01HU	EPOR	Human	ELISA, WB, IHC, IF
Phospho-EPOR (Y426) Antibody	CSB-PA000757	EPOR	Human, Mouse, Rat	ELISA, WB

Product Name	Code	Target	Species	Application
Phospho-EPOR (Y368) Antibody	CSB-PA008072	EPOR	Human, Mouse, Rat	ELISA, WB, IF
EPOR Antibody	CSB-PA008073	EPOR	Human, Mouse, Rat, Monkey	ELISA, WB, IF
EPOR Antibody	CSB-PA008074	EPOR	Human, Mouse, Rat	ELISA, WB
EPOR Antibody	CSB-PA914047	EPOR	Human, Mouse, Rat	ELISA, WB, IHC
EPOR Antibody	CSB-PA980670	EPOR	Human, Mouse, Rat	ELISA, IHC
EPOR (Ab-426) Antibody	CSB-PA826447	EPOR	Human, Mouse, Rat	ELISA, WB
Phospho-EPOR (Tyr368) Antibody	CSB-PA598485	EPOR	Human, Mouse, Rat	ELISA, WB
ERBB2 Antibody	CSB-PA09759A0Rb	ERBB2	Human	ELISA, IHC
Phospho-ERBB2 (Y1221/Y1222) Antibody	CSB-PA000577	ERBB2	Human, Mouse, Rat	ELISA, WB, IHC, IP
Phospho-ERBB2 (Y877) Antibody	CSB-PA000578	ERBB2	Human, Mouse, Rat	ELISA, WB, IHC, IP
Phospho-ERBB2 (Y1221/Y1222) Antibody	CSB-PA000733	ERBB2	Human, Mouse, Rat, Monkey	ELISA, WB, IHC, IF
ERBB2 Antibody	CSB-PA003407	ERBB2	Human, Mouse, Rat	ELISA, WB
Phospho-ERBB2 (Y1112) Antibody	CSB-PA008078	ERBB2	Human, Mouse, Rat	ELISA, WB
ERBB2 Antibody	CSB-PA008081	ERBB2	Human, Mouse, Rat	ELISA, WB
Phospho-ERBB2 (T686) Antibody	CSB-PA008079	ERBB2	Human, Mouse, Rat	ELISA, IF
ERBB2 Antibody	CSB-PA008088	ERBB2	Human, Mouse, Rat	ELISA, IF
ERBB2 Antibody	CSB-PA008089	ERBB2	Human, Mouse, Rat	ELISA, IF
Phospho-ERBB2 (Y1248) Antibody	CSB-PA008077	ERBB2	Human, Mouse, Rat	ELISA, WB, IHC
ERBB2 Antibody	CSB-PA008092	ERBB2	Human, Mouse, Rat	ELISA, WB, IHC
ERBB2 Antibody	CSB-PA003405	ERBB2	Human, Mouse, Rat	ELISA, WB, IHC
ERBB2 Antibody	CSB-PA007763GA01HU	ERBB2	Human	ELISA, WB, IHC
ERBB2 Antibody	CSB-PA972166	ERBB2	Human	ELISA, WB
ERBB2 Antibody	CSB-PA007763KA01HU	ERBB2	Human, Mouse, Rat	ELISA, WB, IHC
ERBB2 (Ab-877) Antibody	CSB-PA200841	ERBB2	Human, Mouse, Rat	ELISA, WB, IHC
ERBB2 (Ab-1221/1222) Antibody	CSB-PA899547	ERBB2	Human	ELISA, WB, IHC, IF
ERBB2 (Ab-1248) Antibody	CSB-PA115264	ERBB2	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-ERBB2 (Tyr1112) Antibody	CSB-PA761682	ERBB2	Human, Mouse, Rat	ELISA, WB
Phospho-ERBB2 (Tyr877) Antibody	CSB-PA846722	ERBB2	Human, Mouse, Rat	ELISA, WB, IHC, IF

Product Name	Code	Target	Species	Application
Phospho-ERBB2 (Tyr1221/Tyr1222) Antibody	CSB-PA281422	ERBB2	Human	ELISA, WB, IF
Phospho-ERBB2 (Tyr1248) Antibody	CSB-PA997735	ERBB2	Human, Mouse, Rat	ELISA, WB, IHC, IF
Cleaved-F10 (A41) Antibody	CSB-PA000067	F10	Human	ELISA, WB
Cleaved-F10 (I235) Antibody	CSB-PA008167	F10	Human, Mouse, Rat	ELISA, WB
F10 Antibody	CSB-PA829543	F10	Human, Mouse, Rat	ELISA, WB
F10 Antibody	CSB-PA245689	F10	Human, Mouse, Rat	ELISA, WB
FCGRT Antibody	CSB-PA008545EA01HU	FCGRT	Human, Mouse	ELISA, WB, IHC, IF
FCGRT Antibody	CSB-PA851443LA01MOV	FCGRT	Macaca fascicularis	ELISA
FCGRT Antibody	CSB-PA008545GA01HU	FCGRT	Human	ELISA, WB
FKRP Antibody	CSB-PA862050LA01HU	FKRP	Human	ELISA, IHC
FKRP Antibody	CSB-PA060135	FKRP	Human, Mouse, Rat	ELISA, WB, IF
FKRP Antibody	CSB-PA035204	FKRP	Human, Mouse, Rat	ELISA, WB
FN1 Antibody	CSB-PA13187A0Rb	FN1	Human	ELISA
Fn1 Antibody	CSB-PA008759HA01RA	Fn1	Rat	ELISA, IHC
FN1 Antibody	CSB-PA002552	FN1	Human, Mouse, Rat	ELISA, WB, IHC
FN1 Antibody	CSB-PA558006	FN1	Human	ELISA, IHC
FN1 Antibody	CSB-PA955751	FN1	Human, Mouse, Rat	ELISA, IHC
FN1 Antibody	CSB-PA266111	FN1	Human, Mouse, Rat	ELISA, IHC
FOXP3 Antibody	CSB-PA15229A0Rb	FOXP3	Human, Mouse	ELISA, WB, IHC
FOXP3 Antibody	CSB-PA005784	FOXP3	Human, Mouse, Rat	ELISA, WB
FOXP3 Antibody	CSB-PA999916	FOXP3	Human, Mouse	ELISA, WB
FUS Antibody	CSB-PA02704A0Rb	FUS	Human	ELISA, WB, IHC, IF
FUS Antibody	CSB-PA009069GA01HU	FUS	Human, Mouse, Rat	ELISA, WB, IHC, IF
FUS Antibody	CSB-PA229936	FUS	Human, Mouse	ELISA, IHC
FUS Antibody	CSB-PA815873	FUS	Human, Mouse	ELISA, IHC
Fxn Antibody	CSB-PA009086LA01MO	Fxn	Mouse, Human	ELISA, WB
FXN Antibody	CSB-PA009086GA01HU	FXN	Human, Mouse, Rat	ELISA, WB
GAA Antibody	CSB-PA009125LA01HU	GAA	Human, Mouse	ELISA, WB, IHC, IF

Product Name	Code	Target	Species	Application
GAA Antibody	CSB-PA009125GA01HU	GAA	Human, Mouse, Rat	ELISA, WB, IHC
gag Antibody	CSB-PA15019A0Rb	gag	Human immunodeficiency virus type 1 group M subtype B	ELISA
gag Antibody	CSB-PA356017DA01MHK	gag	Moloney murine leukemia virus	ELISA
gag Antibody	CSB-PA356017EA01MHK	gag	Moloney murine leukemia virus	ELISA
gag Antibody	CSB-PA356299XA01HKF	gag	Human immunodeficiency virus type 1 group M subtype A (isolate MAL) (HIV-1)	ELISA, WB
gag Antibody	CSB-PA328493XA01HKG	gag	Human immunodeficiency virus type 1 group M subtype A (isolate U455) (HIV-1)	ELISA, WB
GAG Antibody	CSB-PA315187XA01ASN	GAG	Avian leukosis virus RSA (RSV-SRA) (Rous sarcoma virus (strain Schmidt-Ruppin A))	ELISA, WB
gag Antibody	CSB-PA315187XA11ASN	gag	Avian leukosis virus RSA	ELISA, WB
GALC Antibody	CSB-PA009196GA01HU	GALC	Human, Mouse, Rat	ELISA, WB
GALC Antibody	CSB-PA277712	GALC	Human, Mouse	ELISA, WB
GALC Antibody	CSB-PA931313	GALC	Human, Mouse	ELISA, WB
GFAP Antibody	CSB-PA009369LA01HU	GFAP	Human	ELISA, WB, IHC, IF
Gfap Antibody	CSB-PA009369LA01RA	Gfap	Rat, Mouse	ELISA, WB
gfap Antibody	CSB-PA009369LA01DIL	gfap	Zebrafish, Human	ELISA, WB
Phospho-GFAP (S38) Antibody	CSB-PA008785	GFAP	Human	ELISA, WB, IHC, IF
GFAP Antibody	CSB-PA008787	GFAP	Human	ELISA, WB, IHC, IF
GFAP Antibody	CSB-PA009369GA01HU	GFAP	Human, Mouse, Rat	ELISA, WB, IHC
GFAP Antibody	CSB-PA289722	GFAP	Human, Mouse, Rat	ELISA, IHC
GFAP Antibody	CSB-PA047116	GFAP	Human, Rat	ELISA, WB, IF
Phospho-GFAP (Ser38) Antibody	CSB-PA148282	GFAP	Human	ELISA, WB
GHR Antibody	CSB-PA009411LA01HU	GHR	Human	ELISA, IHC
GHR Antibody	CSB-PA006152	GHR	Human, Mouse, Rat	ELISA, WB
GHR Antibody	CSB-PA278805	GHR	Human, Mouse, Rat	ELISA, WB
GHR Antibody	CSB-PA381442	GHR	Human, Mouse, Rat	ELISA, WB
GJB2 Antibody	CSB-PA009452LA01HU	GJB2	Human, Mouse	ELISA, WB, IHC, IF
GJB2 Antibody	CSB-PA001757	GJB2	Human, Mouse, Rat	ELISA, WB
GJB2 Antibody	CSB-PA009452GA01HU	GJB2	Human, Mouse, Rat	ELISA, WB
GJB2 Antibody	CSB-PA549840	GJB2	Human, Rat	ELISA, IHC

Product Name	Code	Target	Species	Application
GLP1R Antibody	CSB-PA009514LA01HU	GLP1R	Human	ELISA, WB, IHC, IF
GLP1R Antibody	CSB-PA009514YA01HU	GLP1R	Human, Mouse	ELISA, WB, IF
Glp1r Antibody	CSB-PA009514HA01MO	Glp1r	Mouse	ELISA, WB
GLP1R Antibody	CSB-PA002699	GLP1R	Human, Mouse	ELISA, WB, IF
GLP1R Antibody	CSB-PA183226	GLP1R	Human, Mouse, Rat	ELISA, IHC
GRN Antibody	CSB-PA009939ESR2HU	GRN	Human	ELISA, IHC
GRN Antibody	CSB-PA009939GA01HU	GRN	Human, Mouse, Rat	ELISA, WB, IHC
HEXA Antibody	CSB-PA12259A0Rb	HEXA	Human, Rat	ELISA, WB, IHC, IF
HEXA Antibody	CSB-PA005178	HEXA	Human, Mouse, Rat	ELISA, WB
HEXA Antibody	CSB-PA010315GA01HU	HEXA	Human, Mouse, Rat	ELISA, WB, IHC
HEXB Antibody	CSB-PA010316LA01HU	HEXB	Human	ELISA, IHC
HEXB Antibody	CSB-PA009172	HEXB	Human	ELISA, WB, IHC, IF
HEXB Antibody	CSB-PA010316GA01HU	HEXB	Human, Mouse, Rat	ELISA, WB, IHC
HEXB Antibody	CSB-PA224713	HEXB	Human	ELISA, WB, IHC
Phospho-HTT (S421) Antibody	CSB-PA009169	HTT	Human, Mouse, Rat	ELISA, IHC
HTT Antibody	CSB-PA009170	HTT	Human, Mouse, Rat	ELISA, IHC
IER3 Antibody	CSB-PA011002LA01HU	IER3	Human	ELISA, IF
Phospho-IFNAR1 (Y466) Antibody	CSB-PA009471	IFNAR1	Human	ELISA, WB, IHC
IFNAR1 Antibody	CSB-PA009473	IFNAR1	Human	ELISA, IHC, IF
IFNAR1 Antibody	CSB-PA107635	IFNAR1	Human	ELISA, IHC
IFNAR1 Antibody	CSB-PA011046GA01HU	IFNAR1	Human	ELISA, WB
Phospho-IFNAR1 (Ser535/Ser539) Antibody	CSB-PA041971	IFNAR1	Human	ELISA, WB
IFNB1 Antibody	CSB-PA06515A0Rb	IFNB1	Human	ELISA, IHC
IFNG Antibody	CSB-PA011050HA01HU	IFNG	Human	ELISA, WB
IFNG Antibody	CSB-PA06524A0Rb	IFNG	Human	ELISA
IFNG Antibody	CSB-PA06525A0Rb	IFNG	Human	ELISA, WB
IFNG Antibody	CSB-PA14709A0Rb	IFNG	Guinea Pig	ELISA, WB
Ifng Antibody	CSB-PA15574A0Rb	Ifng	Rat	ELISA

Product Name	Code	Target	Species	Application
Ifng Antibody	CSB-PA15587A0Rb	Ifng	Mouse	ELISA, IHC
IFNG Antibody	CSB-PA15594A0Rb	IFNG	Chicken	ELISA
IFNG Antibody	CSB-PA011050YA01BO	IFNG	Bovine	ELISA, WB
IFNG Antibody	CSB-PA003000	IFNG	Human	ELISA, WB, IHC
IFNG Antibody	CSB-PA011050GA01HU	IFNG	Human, Mouse, Rat	ELISA, WB, IHC
IFNG Antibody	CSB-PA994634	IFNG	Human	ELISA, IHC
IFNG Antibody	CSB-PA995781	IFNG	Human	ELISA, WB, IHC
IGF1 Antibody	CSB-PA08399A0Rb	IGF1	Human	ELISA, IHC
IGF1 Antibody	CSB-PA009484	IGF1	Human, Mouse, Rat	ELISA, IHC
IGF1 Antibody	CSB-PA070545	IGF1	Human	ELISA, IHC
IGF1 Antibody	CSB-PA011086KA01HU	IGF1	Human, Mouse	ELISA, WB, IHC
IGF1R Antibody	CSB-PA08395A0Rb	IGF1R	Human	ELISA, IF
IGF1R Antibody	CSB-PA009488	IGF1R	Human, Mouse, Rat	ELISA, WB
IGF1R Antibody	CSB-PA009490	IGF1R	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-IGF1R (Y1346) Antibody	CSB-PA009487	IGF1R	Human, Mouse, Rat	ELISA, IHC
IGF1R Antibody	CSB-PA009489	IGF1R	Human, Mouse, Rat	ELISA, WB
Phospho-IGF1R (Y1165/Y1166) Antibody	CSB-PA000549	IGF1R	Human, Mouse, Rat, Monkey	ELISA, WB
IGF1R Antibody	CSB-PA011087GA01HU	IGF1R	Human, Mouse, Rat	ELISA, WB, IP
IGF1R (Ab-1165/1166) Antibody	CSB-PA142373	IGF1R	Human, Mouse, Rat	ELISA, WB
IGF1R (Ab-1280) Antibody	CSB-PA443662	IGF1R	Human, Mouse, Rat	ELISA, WB, IHC
IGF1R (Ab-1346) Antibody	CSB-PA277594	IGF1R	Human, Mouse, Rat	ELISA, WB, IHC
IGF1R (Ab-1161) Antibody	CSB-PA067618	IGF1R	Human, Mouse, Rat	ELISA, WB, IHC, IF
IGF1R (Ab-1346) Antibody	CSB-PA942675	IGF1R	Human, Mouse	ELISA, WB, IHC
Phospho-IGF1R (Tyr1280) Antibody	CSB-PA449888	IGF1R	Human	ELISA, WB
Phospho-IGF1R (Tyr1346) Antibody	CSB-PA945418	IGF1R	Human, Mouse, Rat	ELISA, IHC
Phospho-IGF1R (Tyr1165/Tyr1166) Antibody	CSB-PA597654	IGF1R	Human, Mouse, Rat	ELISA, WB, IF
Phospho-IGF1R (Tyr1161) Antibody	CSB-PA120059	IGF1R	Human, Mouse, Rat	ELISA, WB, IHC, IF
IGF2R Antibody	CSB-PA011093LA01HU	IGF2R	Human	ELISA, IHC, IF

Product Name	Code	Target	Species	Application
IGF2R Antibody	CSB-PA005283	IGF2R	Human	ELISA, WB
Phospho-IGF2R (S2409) Antibody	CSB-PA009491	IGF2R	Human, Mouse, Monkey	ELISA, WB, IHC
IGF2R Antibody	CSB-PA009494	IGF2R	Human, Mouse	ELISA, IHC, IF
IGF2R Antibody	CSB-PA009497	IGF2R	Human, Mouse, Rat	ELISA, IHC, IF
IGF2R Antibody	CSB-PA011093GA01HU	IGF2R	Human	ELISA, WB
IGF2R Antibody	CSB-PA632749	IGF2R	Human, Mouse	ELISA, WB, IHC
IGF2R Antibody	CSB-PA016408	IGF2R	Human, Mouse	ELISA, IHC
Phospho-IGF2R (Ser2409) Antibody	CSB-PA973142	IGF2R	Human, Mouse	ELISA, WB, IHC
IL10 Antibody	CSB-PA011580LA01BO	IL10	Bovine	ELISA
IL10 Antibody	CSB-PA011580LA01MOW	IL10	Macaca mulatta	ELISA
IL10 Antibody	CSB-PA004870	IL10	Human	ELISA, WB
IL10 Antibody	CSB-PA011580KA01HU	IL10	Human, Mouse, Rat	ELISA, WB, IHC
IL10RA Antibody	CSB-PA621688ESR2HU	IL10RA	Human, Mouse	ELISA, WB, IHC
IL10RA Antibody	CSB-PA621688EA01HU	IL10RA	Human	ELISA, WB, IF
IL10RA Antibody	CSB-PA003017	IL10RA	Human	ELISA, WB, IF
IL10RA Antibody	CSB-PA006290	IL10RA	Human	ELISA, WB
Phospho-IL10RA (Y496) Antibody	CSB-PA009575	IL10RA	Human	ELISA, WB
IL10RA Antibody	CSB-PA011581GA01HU	IL10RA	Human	ELISA, WB
IL10RA Antibody	CSB-PA062287	IL10RA	Human	ELISA, IHC
IL10RA Antibody	CSB-PA153938	IL10RA	Human	ELISA, IHC
Phospho-IL10RA (Tyr496) Antibody	CSB-PA051687	IL10RA	Human	ELISA, WB
Phospho-IL10RA (Ser319/323) Antibody	CSB-PA949178	IL10RA	Human	ELISA, WB
IL12A Antibody	CSB-PA16819A0Rb	IL12A	Human	ELISA, IHC
IL12A Antibody	CSB-PA004736	IL12A	Human	ELISA, WB
IL12A Antibody	CSB-PA005776	IL12A	Human	ELISA, WB
IL12A Antibody	CSB-PA009581	IL12A	Human, Mouse, Rat	ELISA, WB, IHC, ICC
IL12A Antibody	CSB-PA877335	IL12A	Human	ELISA, IHC
IL12B Antibody	CSB-PA011587LA01SH	IL12B	Sheep	ELISA

Product Name	Code	Target	Species	Application
IL12B Antibody	CSB-PA011587YA01SH	IL12B	Sheep	ELISA, WB, IHC
IL12B Antibody	CSB-PA06615A0Rb	IL12B	Human, Mouse	ELISA, WB, IHC
IL12B Antibody	CSB-PA009590	IL12B	Human	ELISA, WB
IL12B Antibody	CSB-PA011587GA01HU	IL12B	Human	ELISA, WB
IL13RA2 Antibody	CSB-PA618901ESR1HU	IL13RA2	Human	ELISA, WB, IHC
IL13RA2 Antibody	CSB-PA618901ESR2HU	IL13RA2	Human	ELISA, IHC
IL13RA2 Antibody	CSB-PA009597	IL13RA2	Human	ELISA, WB
IL13RA2 Antibody	CSB-PA006276	IL13RA2	Human	ELISA, WB
IL13RA2 Antibody	CSB-PA011592GA01HU	IL13RA2	Human, Mouse, Rat	ELISA, WB, IHC
IL13RA2 Antibody	CSB-PA239289	IL13RA2	Human	ELISA, IHC
IL13RA2 Antibody	CSB-PA569374	IL13RA2	Human	ELISA, IHC
IL17A Antibody	CSB-PA06659A0Rb	IL17A	Human	ELISA, IHC, IF
IL17A Antibody	CSB-PA011597GA01HU	IL17A	Human	ELISA, IHC
IL1A Antibody	CSB-PA011613LA01MOW	IL1A	Macaca mulatta, Human	ELISA, WB
IL1A Antibody	CSB-PA06565A0Rb	IL1A	Human	ELISA
IL1A Antibody	CSB-PA009537	IL1A	Human, Mouse, Rat	ELISA, IHC
IL1A Antibody	CSB-PA011613GA01HU	IL1A	Human	ELISA, WB
IL1A Antibody	CSB-PA107052	IL1A	Human	ELISA, IHC
IL1A Antibody	CSB-PA584603	IL1A	Human	ELISA, WB, IHC
IL1B Antibody	CSB-PA011614HA01SH	IL1B	Sheep	ELISA, WB
IL1B Antibody	CSB-PA011614YA01SH	IL1B	Sheep	ELISA, WB
IL1B Antibody	CSB-PA06574A0Rb	IL1B	Human	ELISA
Il1b Antibody	CSB-PA17347A0Rb	Il1b	Mouse	ELISA
IL1B Antibody	CSB-PA005122	IL1B	Human, Mouse, Rat	ELISA, WB
IL1B Antibody	CSB-PA003023	IL1B	Human, Mouse, Rat	ELISA, WB, IHC
IL1B Antibody	CSB-PA011614GA01HU	IL1B	Human, Mouse, Rat	ELISA, WB
IL1B Antibody	CSB-PA829263	IL1B	Human, Mouse	ELISA, WB
IL1B Antibody	CSB-PA825267	IL1B	Human, Mouse	ELISA, WB, IHC

Product Name	Code	Target	Species	Application
IL1R1 Antibody	CSB-PA011621ESR1HU	IL1R1	Human, Mouse	ELISA, WB, IHC
IL1R1 Antibody	CSB-PA011621ESR2HU	IL1R1	Human	ELISA, IHC
IL1R1 Antibody	CSB-PA011621EA01HU	IL1R1	Human	ELISA, WB, IF
IL1R1 Antibody	CSB-PA005251	IL1R1	Human	ELISA, WB
Phospho-IL1R1 (Y496) Antibody	CSB-PA009540	IL1R1	Human, Mouse, Rat	ELISA, WB, IHC
IL1R1 Antibody	CSB-PA727952	IL1R1	Human, Mouse	ELISA, IHC
IL1R1 Antibody	CSB-PA072110	IL1R1	Human, Mouse	ELISA, IHC
IL2 Antibody	CSB-PA011629LA01HU	IL2	Human	ELISA, IHC
IL2 Antibody	CSB-PA003024	IL2	Human, Rat	ELISA, WB
IL2RA Antibody	CSB-PA011649LA01HU	IL2RA	Human, Rat	ELISA, WB, IHC, IF
IL2RA Antibody	CSB-PA005198	IL2RA	Human	ELISA, WB
Phospho-IL2RA (S268) Antibody	CSB-PA009542	IL2RA	Human	ELISA, WB, IHC, IF
IL2RA Antibody	CSB-PA009543	IL2RA	Human	ELISA, WB, IF
IL2RA Antibody	CSB-PA904921	IL2RA	Human, Mouse, Rat	ELISA, WB
IL2RA Antibody	CSB-PA013952	IL2RA	Human, Mouse, Rat	ELISA, IHC
Phospho-IL4R (Y497) Antibody	CSB-PA009554	IL4R	Human, Mouse	ELISA, WB, IHC, IF
IL4R Antibody	CSB-PA009556	IL4R	Human, Mouse	ELISA, IHC
IL4R Antibody	CSB-PA003031	IL4R	Human, Mouse	ELISA, WB, IF
IL4R Antibody	CSB-PA512009	IL4R	Human	ELISA, WB
IL6 Antibody	CSB-PA06757A0Rb	IL6	Human, Mouse	ELISA, WB, IHC
IL6 Antibody	CSB-PA005546	IL6	Human	ELISA, WB
IL6R Antibody	CSB-PA006108	IL6R	Human	ELISA, WB
IL6ST Antibody	CSB-PA009560	IL6ST	Human, Mouse	ELISA, WB, IHC, IF
IL7R Antibody	CSB-PA011670ESR1HU	IL7R	Human	ELISA, WB, IHC
IL7R Antibody	CSB-PA011670ESR2HU	IL7R	Human	ELISA, IHC
IL7R Antibody	CSB-PA011670EA01HU	IL7R	Human	ELISA, IHC, IF
Phospho-IL7R (Y449) Antibody	CSB-PA009561	IL7R	Human, Mouse	ELISA, WB, IF
IL7R Antibody	CSB-PA003034	IL7R	Human, Mouse, Monkey	ELISA, WB

Product Name	Code	Target	Species	Application
IL7R Antibody	CSB-PA006261	IL7R	Human	ELISA, WB
IL7R Antibody	CSB-PA011670GA01HU	IL7R	Human, Mouse	ELISA, WB
INS Antibody	CSB-PA12089A0Rb	INS	Human	ELISA, IHC
INS Antibody	CSB-PA011222	INS	Human, Mouse, Rat	ELISA, IHC
INS Antibody	CSB-PA011742GA01HU	INS	Human, Mouse, Rat	ELISA, IHC
INS Antibody	CSB-PA285156	INS	Human	ELISA, IHC
INSR Antibody	CSB-PA011753LA01HU	INSR	Human	ELISA, IHC, IF
INSR Antibody	CSB-PA003043	INSR	Human, Mouse, Rat, Monkey	ELISA, WB
Phospho-INSR (Y1355) Antibody	CSB-PA011225	INSR	Human	ELISA, WB, IHC
Phospho-INSR (Y1361) Antibody	CSB-PA011226	INSR	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-INSR (T1375) Antibody	CSB-PA011224	INSR	Human, Mouse, Rat	ELISA, WB, IHC
INSR Antibody	CSB-PA011227	INSR	Human, Mouse, Rat	ELISA, IF
INSR Antibody	CSB-PA011228	INSR	Human, Mouse, Rat	ELISA, WB, IHC
INSR Antibody	CSB-PA011229	INSR	Human	ELISA, IHC, IF
INSR Antibody	CSB-PA410086	INSR	Human	ELISA, IHC
INSR Antibody	CSB-PA011753GA01HU	INSR	Human, Rat	ELISA, WB
INSR (Ab-1375) Antibody	CSB-PA945193	INSR	Human, Mouse	ELISA, WB
Phospho-INSR (Tyr1361) Antibody	CSB-PA202534	INSR	Human, Mouse, Rat	ELISA, WB
Phospho-INSR (Tyr1355) Antibody	CSB-PA203919	INSR	Human	ELISA, WB
Phospho-INSR (Thr1375) Antibody	CSB-PA994563	INSR	Human, Mouse	ELISA, WB
KRT20 Antibody	CSB-PA012536LA01HU	KRT20	Human	ELISA, WB, IHC, IF
KRT20 Antibody	CSB-PA002048	KRT20	Human, Mouse, Rat	ELISA, WB, IHC, IF
KRT20 Antibody	CSB-PA012536GA01HU	KRT20	Human, Mouse, Rat	ELISA, WB, IHC
KRT20 Antibody	CSB-PA285300	KRT20	Human, Mouse, Rat	ELISA, IHC
KRT20 Antibody	CSB-PA032909	KRT20	Human, Mouse, Rat	ELISA, WB, IHC
LDLR Antibody	CSB-PA09949A0Rb	LDLR	Human	ELISA, IF
LDLR Antibody	CSB-PA012846GA01HU	LDLR	Human, Mouse, Rat	ELISA, WB
LRRK2 Antibody	CSB-PA722493LA01HU	LRRK2	Human	ELISA, IHC, IF

Product Name	Code	Target	Species	Application
LRRK2 Antibody	CSB-PA444926	LRRK2	Human	ELISA, IHC
LRRK2 Antibody	CSB-PA902150	LRRK2	Human	ELISA, IHC
MAOB Antibody	CSB-PA013399LA01HU	MAOB	Human, Mouse	ELISA, WB, IHC
MAOB Antibody	CSB-PA013399GA01HU	MAOB	Human, Mouse, Rat	ELISA, WB, IHC
Mapt Antibody	CSB-PA013481YA01RA	Mapt	Rat	ELISA
MAPT Antibody	CSB-PA11169A0Rb	MAPT	Human, Mouse	ELISA, WB, IHC, IF
Phospho-MAPT (S235) Antibody	CSB-PA000627	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (S404) Antibody	CSB-PA000629	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT Antibody	CSB-PA004234	MAPT	Human	ELISA, WB
Phospho-MAPT (S356) Antibody	CSB-PA010012	MAPT	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-MAPT (T181) Antibody	CSB-PA010016	MAPT	Human	ELISA, WB
Phospho-MAPT (T205) Antibody	CSB-PA010017	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (S214) Antibody	CSB-PA010009	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (S262) Antibody	CSB-PA010010	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (T231) Antibody	CSB-PA010020	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (S519) Antibody	CSB-PA010022	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (T534) Antibody	CSB-PA010023	MAPT	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-MAPT (S516/199) Antibody	CSB-PA010024	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (S422) Antibody	CSB-PA010025	MAPT	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-MAPT (T212) Antibody	CSB-PA010026	MAPT	Human, Mouse, Rat	ELISA, IHC
MAPT Antibody	CSB-PA010039	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT Antibody	CSB-PA010043	MAPT	Human, Mouse, Rat	ELISA, WB, IHC
MAPT Antibody	CSB-PA010060	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT Antibody	CSB-PA010061	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT Antibody	CSB-PA010063	MAPT	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-MAPT (S396) Antibody	CSB-PA000628	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT Antibody	CSB-PA004235	MAPT	Human, Mouse, Rat, Monkey	ELISA, WB
MAPT Antibody	CSB-PA004230	MAPT	Human, Mouse, Rat	ELISA, WB

Product Name	Code	Target	Species	Application
MAPT Antibody	CSB-PA013481GA01HU	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT Antibody	CSB-PA776085	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT Antibody	CSB-PA272079	MAPT	Human, Mouse, Rat	ELISA, IHC
MAPT Antibody	CSB-PA994185	MAPT	Human, Mouse	ELISA, WB, IHC
MAPT Antibody	CSB-PA029620	MAPT	Human, Mouse	ELISA, WB, IHC
MAPT (Ab-356) Antibody	CSB-PA782905	MAPT	Human, Mouse, Rat	ELISA, IF
MAPT (Ab-235) Antibody	CSB-PA028826	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT (Ab-181) Antibody	CSB-PA036348	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT (Ab-205) Antibody	CSB-PA247791	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT (Ab-214) Antibody	CSB-PA265670	MAPT	Human, Mouse, Rat	ELISA, IF
MAPT (Ab-262) Antibody	CSB-PA103274	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT (Ab-404) Antibody	CSB-PA214875	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT (Ab-212) Antibody	CSB-PA549925	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT (Ab-422) Antibody	CSB-PA087062	MAPT	Human, Mouse, Rat	ELISA, WB
Mapt Antibody	CSB-PA945056	Mapt	Human, Mouse, Rat	ELISA, WB
MAPT (Ab-235) Antibody	CSB-PA795646	MAPT	Human, Mouse, Rat	ELISA, WB
MAPT (Ab-534/217) Antibody	CSB-PA786032	MAPT	Human, Mouse	ELISA, WB
MAPT (Ab-717/400) Antibody	CSB-PA112612	MAPT	Human, Mouse	ELISA, WB
MAPT (Ab-396) Antibody	CSB-PA727104	MAPT	Human, Mouse, Rat	ELISA, WB, IHC
MAPT (Ab-231) Antibody	CSB-PA235644	MAPT	Human, Mouse, Rat	ELISA, WB, IHC
MAPT (Ab-396) Antibody	CSB-PA904096	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (Ser235) Antibody	CSB-PA076235	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (Ser214) Antibody	CSB-PA135965	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (Thr212) Antibody	CSB-PA237372	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (Ser422) Antibody	CSB-PA277096	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (Ser416) Antibody	CSB-PA050509	MAPT	Human, Rat, Mouse	ELISA, WB
Phospho-MAPT (Ser202) Antibody	CSB-PA257727	MAPT	Human, Mouse, Rat	ELISA, WB
Phospho-MAPT (Thr217) Antibody	CSB-PA141564	MAPT	Human, Mouse, Rat	ELISA, WB

Product Name	Code		Target	Species	Application
Phospho-MAPT (Ser516/199) Antibody	CSB-PA176217	MAPT	Human, Mouse, Rat		ELISA, WB
Phospho-MAPT (Ser356) Antibody	CSB-PA252810	MAPT	Human, Mouse, Rat		ELISA, WB, IF
Phospho-MAPT (Ser396) Antibody	CSB-PA563392	MAPT	Human, Mouse, Rat		ELISA, WB, IHC
Phospho-MAPT (Thr181) Antibody	CSB-PA145569	MAPT	Human, Mouse, Rat		ELISA, WB, IHC
Phospho-MAPT (Thr205) Antibody	CSB-PA271293	MAPT	Human, Mouse, Rat		ELISA, WB, IHC
Phospho-MAPT (Thr231) Antibody	CSB-PA832791	MAPT	Human, Mouse, Rat		ELISA, WB, IHC
Phospho-MAPT (Ser262) Antibody	CSB-PA045644	MAPT	Human, Mouse, Rat		ELISA, WB, IHC, IF
Phospho-MAPT (Ser404) Antibody	CSB-PA963344	MAPT	Human, Mouse, Rat		ELISA, WB, IHC, IF
MBP Antibody	CSB-PA20287A0Rb	MBP	Macaca fascicularis		ELISA
MBP Antibody	CSB-PA013551LA01HU	MBP	Human, Mouse, Rat		ELISA, WB, IHC, IF
MBP Antibody	CSB-PA013551GA01HU	MBP	Human, Mouse, Rat		ELISA, WB, IHC
MBP Antibody	CSB-PA261629	MBP	Human, Mouse, Rat		ELISA, WB
MECP2 Antibody	CSB-PA013638ESR2HU	MECP2	Human		ELISA, IHC
MECP2 Antibody	CSB-PA000357	MECP2	Human		WB
MECP2 Antibody	CSB-PA013638GA01HU	MECP2	Human, Mouse, Rat		ELISA, WB, IHC
MECP2 Antibody	CSB-PA440474	MECP2	Human, Mouse, Rat		ELISA, WB, IHC
MECP2 Antibody	CSB-PA912786	MECP2	Human		ELISA, WB
MOG Antibody	CSB-PA619083ESR2HU	MOG	Human		ELISA, IHC
MOG Antibody	CSB-PA014709GA01HU	MOG	Human, Mouse, Rat		ELISA, WB
Ms4a1 Antibody	CSB-PA015007LA01MO	Ms4a1	Mouse, Human		ELISA, WB
MS4A1 Antibody	CSB-PA015007DSR2HU	MS4A1	Human		ELISA, IHC
MS4A1 Antibody	CSB-PA001433	MS4A1	Human		ELISA, WB
MS4A1 Antibody	CSB-PA006161	MS4A1	Human		ELISA, WB
MSTN Antibody	CSB-PA11869A0Rb	MSTN	Human		ELISA, WB, IHC
MSTN Antibody	CSB-PA008776	MSTN	Human, Mouse		ELISA, IHC
MSTN Antibody	CSB-PA015057GA01HU	MSTN	Human, Mouse, Rat		ELISA, WB, IHC
MSTN Antibody	CSB-PA293957	MSTN	Human, Mouse, Rat		ELISA, IHC
MTM1 Antibody	CSB-PA614418ESR1HU	MTM1	Human		ELISA, WB, IHC

Product Name	Code	Target	Species	Application
MTM1 Antibody	CSB-PA614418ESR2HU	MTM1	Human	ELISA, IHC
MTM1 Antibody	CSB-PA003350	MTM1	Human, Mouse	ELISA, WB
MTM1 Antibody	CSB-PA015173GA01HU	MTM1	Human, Mouse, Rat	ELISA, WB, IHC
MTOR Antibody	CSB-PA008968HA01HU	MTOR	Human	ELISA, IHC, IF
Phospho-MTOR (T2446) Antibody	CSB-PA000652	MTOR	Human, Mouse, Rat	ELISA, WB
MTOR Antibody	CSB-PA008481	MTOR	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-MTOR (S2481) Antibody	CSB-PA008479	MTOR	Human, Mouse, Rat	ELISA, IHC, IF
Phospho-MTOR (S2448) Antibody	CSB-PA000576	MTOR	Human, Mouse, Rat	ELISA, WB, IHC, IF
MTOR Antibody	CSB-PA003333	MTOR	Human, Mouse, Rat	ELISA, WB, IHC, IF
MTOR Antibody	CSB-PA008968GA01HU	MTOR	Human, Mouse, Rat	ELISA, WB, IHC
MTOR Antibody	CSB-PA876444	MTOR	Human, Mouse, Rat	ELISA, IHC
MTOR Antibody	CSB-PA273227	MTOR	Human, Mouse, Rat	ELISA, IHC
MTOR (Ab-2481) Antibody	CSB-PA054657	MTOR	Human, Mouse, Rat	ELISA, IHC
MTOR (Ab-2448) Antibody	CSB-PA208208	MTOR	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-MTOR (Ser2448) Antibody	CSB-PA271384	MTOR	Human, Mouse, Rat	ELISA, WB, IHC, IF
NFKB1 Antibody	CSB-PA10354A0Rb	NFKB1	Human	ELISA, WB, IHC, IF, ChIP
Phospho-NFKB1 (S893) Antibody	CSB-PA000579	NFKB1	Human	ELISA, WB, IHC
Phospho-NFKB1 (S907) Antibody	CSB-PA000580	NFKB1	Human	ELISA, WB, IHC, IP
Phospho-NFKB1 (S932) Antibody	CSB-PA000581	NFKB1	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-NFKB1 (S337) Antibody	CSB-PA000683	NFKB1	Human, Mouse	ELISA, WB, IHC
Phospho-NFKB1 (S927) Antibody	CSB-PA000730	NFKB1	Human, Mouse	ELISA, WB, IHC, IF
NFKB1 Antibody	CSB-PA003431	NFKB1	Human, Mouse	ELISA, WB, IHC, IF
NFKB1 Antibody	CSB-PA004763	NFKB1	Human	ELISA, WB
Phospho-NFKB1 (S337) Antibody	CSB-PA020084	NFKB1	Human, Mouse, Rat	ELISA, WB, IHC, IF
NFKB1 Antibody	CSB-PA020085	NFKB1	Human, Mouse, Rat	ELISA, WB, IHC
NFKB1 Antibody	CSB-PA020086	NFKB1	Human, Mouse, Rat	ELISA, IHC
NFKB1 Antibody	CSB-PA020090	NFKB1	Human	ELISA, WB, IHC
NFKB1 Antibody	CSB-PA020096	NFKB1	Human, Mouse, Rat	ELISA, WB, IHC

Product Name	Code	Target	Species	Application
NFKB1 Antibody	CSB-PA003430	NFKB1	Human	ELISA, WB, IHC, IP
NFKB1 Antibody	CSB-PA003434	NFKB1	Human, Mouse	ELISA, WB, IHC, IF
NFKB1 Antibody	CSB-PA015759GA01HU	NFKB1	Human, Mouse, Rat	ELISA, WB
NFKB1 Antibody	CSB-PA259078	NFKB1	Human, Mouse, Rat	ELISA, IHC
NFKB1 Antibody	CSB-PA546304	NFKB1	Human, Mouse, Rat	ELISA, IHC
NFKB1 Antibody	CSB-PA441503	NFKB1	Human, Mouse, Rat	ELISA, WB
NFKB1 (Ab-337) Antibody	CSB-PA186775	NFKB1	Human, Mouse, Rat	ELISA, WB, IHC
NFKB1 (Ab-893) Antibody	CSB-PA108426	NFKB1	Human	ELISA, WB, IHC
NFKB1 (Ab-907) Antibody	CSB-PA696291	NFKB1	Human	ELISA, WB, IHC
NFKB1 (Ab-932) Antibody	CSB-PA190132	NFKB1	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-NFKB1 (Ser932) Antibody	CSB-PA142672	NFKB1	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-NFKB1 (Ser893) Antibody	CSB-PA258273	NFKB1	Human	ELISA, IHC, IF
Phospho-NFKB1 (Ser907) Antibody	CSB-PA925652	NFKB1	Human	ELISA, WB, IHC
Phospho-NFKB1 (Ser927) Antibody	CSB-PA283590	NFKB1	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-NFKB1 (Ser337) Antibody	CSB-PA632717	NFKB1	Human, Mouse, Rat	ELISA, WB, IHC, IF
Ngf Antibody	CSB-PA16859A0Rb	Ngf	Rat	ELISA, IP
NGF Antibody	CSB-PA643662LA01PI	NGF	Pig	ELISA
NGF Antibody	CSB-PA643662LA01HU	NGF	Human	ELISA, IHC
NGF Antibody	CSB-PA006132	NGF	Human	ELISA, WB
NGF Antibody	CSB-PA003446	NGF	Human, Mouse, Rat	ELISA, WB, IHC
NGF Antibody	CSB-PA015779GA01HU	NGF	Human	ELISA, IHC
NGF Antibody	CSB-PA942456	NGF	Human, Mouse, Rat	ELISA, IHC
NGF Antibody	CSB-PA558560	NGF	Human, Mouse, Rat	ELISA, IHC
NGFR Antibody	CSB-PA10479A0Rb	NGFR	Human, Mouse	ELISA, WB, IHC
NGFR Antibody	CSB-PA003447	NGFR	Human, Mouse, Rat	ELISA, WB
NGFR Antibody	CSB-PA003448	NGFR	Human, Mouse, Rat	ELISA, WB, IF
NGFR Antibody	CSB-PA015780GA01HU	NGFR	Human, Mouse, Rat	ELISA, WB, IHC
NGFR Antibody	CSB-PA319979	NGFR	Human, Mouse, Rat	ELISA, WB, IHC

Product Name	Code	Target	Species	Application
NGFR Antibody	CSB-PA084270	NGFR	Human, Mouse, Rat	ELISA, WB, IHC
NGFR Antibody	CSB-PA038076	NGFR	Human, Mouse, Rat	ELISA, WB, IF
NTRK2 Antibody	CSB-PA619082ESR1HU	NTRK2	Human	ELISA, IHC
Phospho-NTRK2 (Y516) Antibody	CSB-PA020171	NTRK2	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-NTRK2 (Y706) Antibody	CSB-PA020172	NTRK2	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-NTRK2 (Y706/Y707) Antibody	CSB-PA020173	NTRK2	Human, Mouse, Rat	ELISA, IHC
NTRK2 Antibody	CSB-PA004332	NTRK2	Human, Mouse, Rat	ELISA, WB, IHC
NTRK2 Antibody	CSB-PA020174	NTRK2	Human, Mouse, Rat	ELISA, WB, IHC, IF
NTRK2 Antibody	CSB-PA016134GA01HU	NTRK2	Human, Mouse, Rat	ELISA, WB
NTRK2 (Ab-705) Antibody	CSB-PA292268	NTRK2	Human, Mouse, Rat	ELISA, WB, IF
NTRK2 Antibody	CSB-PA936032	NTRK2	Human, Mouse	ELISA, WB, IHC
Phospho-NTRK2 (Tyr515) Antibody	CSB-PA197704	NTRK2	Human, Mouse, Rat	ELISA, IF
Phospho-NTRK2 (Tyr705) Antibody	CSB-PA108948	NTRK2	Human, Mouse, Rat	ELISA, WB
Phospho-NTRK2 (Tyr706/Tyr707) Antibody	CSB-PA789970	NTRK2	Human, Mouse, Rat	ELISA, IHC
OPTN Antibody	CSB-PA016363GA01HU	OPTN	Human, Mouse, Rat	ELISA, WB, IHC, IF
OPTN Antibody	CSB-PA548038	OPTN	Human, Mouse, Rat	ELISA, IHC
OPTN Antibody	CSB-PA924414	OPTN	Human, Mouse, Rat	ELISA, IHC
OTC Antibody	CSB-PA017270LA01HU	OTC	Human	ELISA, IHC
OTC Antibody	CSB-PA189809	OTC	Human, Mouse, Rat	ELISA, WB, IHC
OTC Antibody	CSB-PA033931	OTC	Human, Mouse, Rat	ELISA, WB, IHC
PAH Antibody	CSB-PA017396LA01HU	PAH	Human, Mouse	ELISA, WB, IHC, IF
PAH Antibody	CSB-PA003698	PAH	Human, Mouse, Rat	ELISA, WB, IHC
PARK7 Antibody	CSB-PA03724A0Rb	PARK7	Human, Mouse, Rat	ELISA, WB, IHC
PARK7 Antibody	CSB-PA03729A0Rb	PARK7	Human	ELISA, WB, IHC, IF
PARK7 Antibody	CSB-PA002152	PARK7	Human, Mouse, Rat	ELISA, WB, IHC
PARK7 Antibody	CSB-PA003707	PARK7	Human, Mouse	ELISA, WB, IHC, IF
PARK7 Antibody	CSB-PA017454GA01HU	PARK7	Human, Mouse, Rat	ELISA, WB, IHC
PARK7 Antibody	CSB-PA948689	PARK7	Human, Mouse, Rat	ELISA, WB, IHC

Product Name	Code	Target	Species	Application
PARK7 Antibody	CSB-PA261296	PARK7	Human, Mouse, Rat	ELISA, WB
PARK7 Antibody	CSB-PA108078	PARK7	Human, Mouse, Rat	ELISA, WB, IHC, IF
PCSK9 Antibody	CSB-PA07039A0Rb	PCSK9	Human	ELISA, IHC, IF
PCSK9 Antibody	CSB-PA017647LA01HU	PCSK9	Human	ELISA, WB, IHC, IF
PCSK9 Antibody	CSB-PA017647GA01HU	PCSK9	Human, Mouse, Rat	ELISA, WB
PCSK9 Antibody	CSB-PA638311	PCSK9	Human	ELISA, WB
PCSK9 Antibody	CSB-PA108805	PCSK9	Human	ELISA, WB
PDCD1 Antibody	CSB-PA619964LA01HU	PDCD1	Human, Mouse	ELISA, WB, IHC, IF
PDCD1 Antibody	CSB-PA017664GA01HU	PDCD1	Human	ELISA, WB, IHC
PDCD1 Antibody	CSB-PA483633	PDCD1	Human	ELISA, IHC
PINK1 Antibody	CSB-PA863144LA01HU	PINK1	Human	ELISA, IHC
PINK1 Antibody	CSB-PA447054	PINK1	Human, Mouse	ELISA
PLG Antibody	CSB-PA018188LA01HU	PLG	Human	ELISA, IHC
Cleaved-PLG (V98) Antibody	CSB-PA000083	PLG	Human, Mouse, Rat	ELISA, WB
PLG Antibody	CSB-PA030085	PLG	Human	ELISA, WB
PLG Antibody	CSB-PA018188GA01HU	PLG	Human, Mouse, Rat	ELISA, WB
PLG Antibody	CSB-PA050445	PLG	Human, Rat	ELISA, WB
PLG Antibody	CSB-PA157482	PLG	Human, Rat	ELISA, WB, IHC
PLG Antibody	CSB-PA970985	PLG	Human	ELISA, WB
PPP1R1B Antibody	CSB-PA880070DSR1HU	PPP1R1B	Human	ELISA, IHC
Phospho-PPP1R1B (T75) Antibody	CSB-PA000510	PPP1R1B	Human, Mouse, Rat, Monkey	ELISA, WB, IHC
Phospho-PPP1R1B (T34) Antibody	CSB-PA000742	PPP1R1B	Human, Mouse, Rat	ELISA, WB, IHC, IF
PPP1R1B Antibody	CSB-PA002073	PPP1R1B	Human, Mouse, Rat, Monkey	ELISA, WB, IHC
PPP1R1B Antibody	CSB-PA002075	PPP1R1B	Human, Mouse, Rat	ELISA, WB, IHC
PPP1R1B Antibody	CSB-PA018535GA01HU	PPP1R1B	Human, Mouse, Rat	ELISA, WB, IHC
PPP1R1B Antibody	CSB-PA961944	PPP1R1B	Human, Rat	ELISA, IHC
PPP1R1B Antibody	CSB-PA981476	PPP1R1B	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-PPP1R1B (Thr75) Antibody	CSB-PA778974	PPP1R1B	Human, Mouse, Rat	ELISA, IF

Product Name	Code	Target	Species	Application
PPT1 Antibody	CSB-PA018587LA01HU	PPT1	Human	ELISA, WB, IHC
PPT1 Antibody	CSB-PA030196	PPT1	Human, Mouse, Rat	ELISA, IHC
PPT1 Antibody	CSB-PA018587GA01HU	PPT1	Human, Mouse, Rat	ELISA, IHC
PRKN Antibody	CSB-PA017451LA01HU	PRKN	Human	ELISA, WB, IHC, IF
PARK2 Antibody	CSB-PA003708	PRKN	Human	ELISA, WB, IHC
Phospho-PARK2 (S131) Antibody	CSB-PA020235	PRKN	Human	ELISA, WB, IHC
PARK2 Antibody	CSB-PA020236	PRKN	Human, Mouse, Rat	ELISA, WB, IHC, IF
PARK2 Antibody	CSB-PA076233	PRKN	Human, Mouse, Rat	ELISA, WB, IHC
PSEN1 Antibody	CSB-PA018846ESR2HU	PSEN1	Human	ELISA, WB, IHC
PSEN1 Antibody	CSB-PA018846EA01HU	PSEN1	Human	ELISA, WB, IHC
PSEN1 Antibody	CSB-PA040015	PSEN1	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-PSEN1 (S357) Antibody	CSB-PA040014	PSEN1	Human, Mouse, Rat	ELISA, WB, IHC
PSEN1 Antibody	CSB-PA040016	PSEN1	Human	ELISA, WB, IHC
PSEN1 Antibody	CSB-PA018846GA01HU	PSEN1	Human, Mouse, Rat	ELISA, WB
PSEN1 Antibody	CSB-PA695720	PSEN1	Human	ELISA, WB, IHC
PYCARD Antibody	CSB-PA890936LA01HU	PYCARD	Human, Mouse	ELISA, WB, IHC
PYCARD Antibody	CSB-PA000940	PYCARD	Human, Mouse, Rat	ELISA, WB, IHC, IF
PYCARD Antibody	CSB-PA019114GA01HU	PYCARD	Human	ELISA, WB, IHC
RGMA Antibody	CSB-PA836186LA01HU	RGMA	Human, Mouse	ELISA, WB, IHC
RGMA Antibody	CSB-PA019628GA01HU	RGMA	Human, Mouse, Rat	ELISA, WB
RTN4 Antibody	CSB-PA003461	RTN4	Human, Mouse, Rat	ELISA, WB, IHC
RTN4 Antibody	CSB-PA020572GA01HU	RTN4	Human, Mouse, Rat	ELISA, WB, IHC
RTN4 Antibody	CSB-PA931699	RTN4	Human, Mouse, Rat	ELISA, IHC
RTN4 Antibody	CSB-PA054695	RTN4	Human, Mouse, Rat	ELISA, WB, IHC
SCN1A Antibody	CSB-PA020834ESR1HU	SCN1A	Human	ELISA
SCN1A Antibody	CSB-PA409848	SCN1A	Human, Rat	ELISA, IHC
SCN1A Antibody	CSB-PA044197	SCN1A	Human, Rat	ELISA, IHC
SCN2A Antibody	CSB-PA040291	SCN2A	Human, Mouse, Rat	ELISA, IHC

Product Name	Code	Target	Species	Application
SCN2A Antibody	CSB-PA040292	SCN2A	Human, Mouse, Rat	ELISA, IHC
SCN2A Antibody	CSB-PA937791	SCN2A	Human, Rat	ELISA, IHC
SCN2A Antibody	CSB-PA902560	SCN2A	Human, Rat	ELISA, IHC
SELP Antibody	CSB-PA020978LA01HU	SELP	Human	ELISA, IF
SELP Antibody	CSB-PA006202	SELP	Human, Mouse, Rat	ELISA, WB
SEMA4D Antibody	CSB-PA835707LA01HU	SEMA4D	Human	ELISA, IF
SEMA4D Antibody	CSB-PA005238	SEMA4D	Human, Mouse, Rat	ELISA, WB
SIRT3 Antibody	CSB-PA882102ESR1HU	SIRT3	Human	ELISA, WB, IHC
SIRT3 Antibody	CSB-PA882102ESR2HU	SIRT3	Human	ELISA, WB, IHC
SIRT3 Antibody	CSB-PA882102LA01HU	SIRT3	Human, Mouse	ELISA, WB, IF
SIRT3 Antibody	CSB-PA008365	SIRT3	Human	ELISA, WB
SIRT3 Antibody	CSB-PA100396	SIRT3	Human, Mouse	ELISA, IHC
SIRT3 Antibody	CSB-PA937428	SIRT3	Human, Mouse	ELISA, WB, IHC
SIRT3 Antibody	CSB-PA229200	SIRT3	Human	ELISA, WB
SLC6A3 Antibody	CSB-PA021701LA01HU	SLC6A3	Human, Mouse	ELISA, WB, IF
SLC6A3 Antibody	CSB-PA198930	SLC6A3	Human, Mouse, Rat	ELISA, IHC
SLC6A3 Antibody	CSB-PA791171	SLC6A3	Human, Mouse, Rat	ELISA, IHC
SMN1 Antibody	CSB-PA614982LA01HU	SMN1	Human	ELISA, WB, IHC, IF
SMN1 Antibody	CSB-PA135133	SMN1	Human	ELISA, WB
SMN1 Antibody	CSB-PA678622	SMN1	Human	ELISA, WB
SNAP25 Antibody	CSB-PA021873ESR1HU	SNAP25	Human	ELISA, IHC
SNAP25 Antibody	CSB-PA021873ESR2HU	SNAP25	Human	ELISA, IHC
SNAP25 Antibody	CSB-PA050112	SNAP25	Human, Mouse, Rat	ELISA, WB, IHC
SNAP25 Antibody	CSB-PA021873GA01HU	SNAP25	Human, Mouse, Rat	ELISA, WB
SNAP25 Antibody	CSB-PA564264	SNAP25	Human, Mouse, Rat	ELISA, WB, IHC
SNAP25 Antibody	CSB-PA585168	SNAP25	Human, Mouse, Rat	ELISA, WB, IHC
Snca Antibody	CSB-PA021912DA01RA	Snca	Rat	ELISA
SNCA Antibody	CSB-PA06994A0Rb	SNCA	Human	ELISA, WB, IHC

Product Name	Code	Target	Species	Application
SNCA Antibody	CSB-PA06997A0Rb	SNCA	Human	ELISA, IHC, IF
SNCA Antibody	CSB-PA021912DSR1HU	SNCA	Human	ELISA, IHC
SNCA Antibody	CSB-PA021912OA01HU	SNCA	Human	ELISA, IHC, IF
Phospho-SNCA (Y125) Antibody	CSB-PA000748	SNCA	Human, Mouse, Rat	ELISA, WB, IHC
SNCA Antibody	CSB-PA004209	SNCA	Human, Mouse, Rat	ELISA, WB, IHC
SNCA Antibody	CSB-PA004210	SNCA	Human, Mouse, Rat	ELISA, WB, IHC
SNCA Antibody	CSB-PA004212	SNCA	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-SNCA (S129) Antibody	CSB-PA050113	SNCA	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-SNCA (Y133) Antibody	CSB-PA050114	SNCA	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-SNCA (Y136) Antibody	CSB-PA050115	SNCA	Human, Mouse, Rat	ELISA, WB
SNCA Antibody	CSB-PA050116	SNCA	Human, Mouse, Rat	ELISA, WB, IHC
SNCA Antibody	CSB-PA021912GA01HU	SNCA	Human, Mouse, Rat	ELISA, WB, IHC
SNCA Antibody	CSB-PA592134	SNCA	Human, Mouse, Rat	ELISA, IHC
SNCA Antibody	CSB-PA169656	SNCA	Human, Mouse, Rat	ELISA, IHC
SNCA (Ab-129) Antibody	CSB-PA238735	SNCA	Human, Mouse, Rat	ELISA, WB
SNCA (Ab-125) Antibody	CSB-PA578342	SNCA	Human, Mouse, Rat	ELISA, WB
SNCA (Ab-133) Antibody	CSB-PA783926	SNCA	Human, Mouse, Rat	ELISA, WB, IF
SNCA Antibody	CSB-PA036896	SNCA	Human, Mouse, Rat	ELISA, WB, IHC
Phospho-SNCA (Ser129) Antibody	CSB-PA585265	SNCA	Human, Mouse, Rat	ELISA, WB
Phospho-SNCA (Tyr125) Antibody	CSB-PA145405	SNCA	Human	ELISA, WB, IF
Phospho-SNCA (Tyr133) Antibody	CSB-PA053564	SNCA	Human, Mouse, Rat	ELISA, WB, IF
Sod1 Antibody	CSB-PA022397LA01MO	Sod1	Mouse, Rat, Human	ELISA, WB, IHC, IP
Sod1 Antibody	CSB-PA022397YA01MO	Sod1	Mouse	ELISA
SOD1 Antibody	CSB-PA02864A0Rb	SOD1	Human, Mouse	ELISA, WB, IHC, IF
SOD1 Antibody	CSB-PA004131	SOD1	Human, Mouse, Rat	ELISA, WB
SOD1 Antibody	CSB-PA080200	SOD1	Human, Mouse, Rat	WB, IHC
SOD1 Antibody	CSB-PA022397GA01HU	SOD1	Human, Mouse, Rat	ELISA, WB, IHC
SOD1 Antibody	CSB-PA238446	SOD1	Human, Mouse, Rat	ELISA, WB

Product Name	Code	Target	Species	Application
sod1 Antibody	CSB-PA022397XA01DIL	sod1	Danio rerio (Zebrafish) (Brachydanio rerio)	ELISA, WB
Sort1 Antibody	CSB-PA022412LA01RA	Sort1	Rat	ELISA
SORT1 Antibody	CSB-PA859102LA01HU	SORT1	Human	ELISA, IHC
SORT1 Antibody	CSB-PA022412GA01HU	SORT1	Human, Mouse, Rat	ELISA, WB, IHC
SORT1 Antibody	CSB-PA448728	SORT1	Human, Mouse, Rat	ELISA, WB, IHC
SORT1 Antibody	CSB-PA194090	SORT1	Human, Mouse, Rat	ELISA, IHC
SQSTM1 Antibody	CSB-PA615696LA01HU	SQSTM1	Human	ELISA, IHC, IF, IP
SQSTM1 Antibody	CSB-PA022647GA01HU	SQSTM1	Human	ELISA, WB, IF
TARDBP Antibody	CSB-PA621652ESR1HU	TARDBP	Human, Mouse	ELISA, WB, IHC, IP
TARDBP Antibody	CSB-PA023129EA01HU	TARDBP	Human	ELISA, WB, IHC, IF
TARDBP Antibody	CSB-PA023129GA01HU	TARDBP	Human, Mouse, Rat, Zebrafish	ELISA, WB, IHC, IF, IP
TBK1 Antibody	CSB-PA023232LA01HU	TBK1	Human, Mouse	ELISA, WB, IHC, IF
TBK1 Antibody	CSB-PA890690NA01HU	TBK1	Human	ELISA, WB, IF
TBK1 (Ab-172) Antibody	CSB-PA890690NA172nphHU	TBK1	Human	ELISA, IHC
TBK1 Antibody	CSB-PA009019	TBK1	Human, Mouse	ELISA, IHC
TFRC Antibody	CSB-PA07219A0Rb	TFRC	Human	ELISA, IHC, IF
TFRC Antibody	CSB-PA001477	TFRC	Human, Mouse	ELISA, WB, IHC
TFRC Antibody	CSB-PA005607	TFRC	Human	ELISA, WB
Phospho-TFRC (S24) Antibody	CSB-PA010092	TFRC	Human, Mouse	ELISA, WB, IHC
TFRC Antibody	CSB-PA023441GA01HU	TFRC	Human, Mouse, Rat	ELISA, WB, IHC
TFRC Antibody	CSB-PA104885	TFRC	Human, Mouse, Rat	ELISA, WB, IHC
TFRC Antibody	CSB-PA834425	TFRC	Human, Mouse, Rat	ELISA, WB, IHC
TGFB1 Antibody	CSB-PA023446LA01BO	TGFB1	Bovine	ELISA
TGFB1 Antibody	CSB-PA023446LA01CH	TGFB1	Chicken	ELISA
tgfb1 Antibody	CSB-PA023446LA01OEI	tgfb1	Oncorhynchus mykiss, Human	ELISA, WB
TGFB1 Antibody	CSB-PA07309A0Rb	TGFB1	Human	ELISA, IHC
Tgfb1 Antibody	CSB-PA023446LA01MO	Tgfb1	Mouse, Human	ELISA, WB, IP
TGFB1 Antibody	CSB-PA080237	TGFB1	Human, Mouse, Rat	WB

Product Name	Code	Target	Species	Application
TGFB1 Antibody	CSB-PA080233	TGFB1	Human, Mouse, Rat	WB, IHC
TGFB1 Antibody	CSB-PA004279	TGFB1	Human, Mouse, Rat	ELISA, WB, IHC, IF
TGFB1 Antibody	CSB-PA023446GA01HU	TGFB1	Human, Mouse, Rat	ELISA, IHC
TGFB1 Antibody	CSB-PA697937	TGFB1	Human, Mouse, Rat	ELISA, WB, IHC
TGFB1 Antibody	CSB-PA060587	TGFB1	Human, Mouse, Rat	ELISA, IHC
TGFB1 Antibody	CSB-PA786987	TGFB1	Human, Mouse, Rat	ELISA, WB, IHC
TGFB1 Antibody	CSB-PA286763	TGFB1	Human, Mouse, Rat	ELISA, WB
TH Antibody	CSB-PA023470LA01HU	TH	Human, Mouse, Rat	ELISA, WB, IHC
th Antibody	CSB-PA24419A0Rb	th	Zebrafish	ELISA, WB
TH Antibody	CSB-PA004284	TH	Human, Mouse, Rat	ELISA, WB, IHC, IF
TH Antibody	CSB-PA004286	TH	Rat	ELISA, WB, IHC
Phospho-TH (S8) Antibody	CSB-PA006345	TH	Rat	ELISA, WB, IHC
Phospho-TH (S62) Antibody	CSB-PA050254	TH	Human, Mouse, Rat, Monkey	ELISA, WB, IHC
TH Antibody	CSB-PA050255	TH	Human, Mouse, Rat, Monkey	ELISA, WB, IHC
Phospho-TH (S19) Antibody	CSB-PA000731	TH	Human, Mouse, Rat	ELISA, WB, IHC, IF
TH Antibody	CSB-PA004285	TH	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-TH (S71) Antibody	CSB-PA011030	TH	Human, Mouse, Rat	ELISA, WB, IHC
TH Antibody	CSB-PA195276	TH	Human	ELISA, IHC
TH Antibody	CSB-PA903853	TH	Human	ELISA, IHC
TH (Ab-31) Antibody	CSB-PA695636	TH	Mouse, Rat	ELISA, WB
TH (Ab-19) Antibody	CSB-PA219605	TH	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-TH (Ser40) Antibody	CSB-PA154990	TH	Human, Mouse, Rat	ELISA, IF
Phospho-TH (Ser19) Antibody	CSB-PA439504	TH	Rat	ELISA, IF
TLR2 Antibody	CSB-PA023601LA01HU	TLR2	Human	ELISA
TLR2 Antibody	CSB-PA080185	TLR2	Human, Mouse, Rat	ELISA, IHC
TLR2 Antibody	CSB-PA205006	TLR2	Human, Mouse	ELISA, IHC
TLR2 Antibody	CSB-PA183593	TLR2	Human, Mouse	ELISA, IHC
TLR4 Antibody	CSB-PA080181	TLR4	Human, Mouse, Rat	ELISA, IHC

Product Name	Code	Target	Species	Application
TLR4 Antibody	CSB-PA001434	TLR4	Human, Mouse	ELISA, WB, IHC
TLR4 Antibody	CSB-PA023603GA01HU	TLR4	Human, Mouse, Rat, Bovin	ELISA, WB
TLR4 Antibody	CSB-PA237733	TLR4	Human, Mouse, Rat	ELISA, IHC
TLR4 Antibody	CSB-PA902590	TLR4	Human, Mouse, Rat	ELISA, WB, IHC
TLR4 Antibody	CSB-PA776633	TLR4	Human	ELISA, IHC
TLR4 Antibody	CSB-PA585936	TLR4	Human	ELISA, IHC
TNF Antibody	CSB-PA023955DA01SH	TNF	Sheep	ELISA, WB
TNF Antibody	CSB-PA07427A0Rb	TNF	Human	ELISA, WB, IHC, IF
TNF Antibody	CSB-PA08474A0Rb	TNF	Bovine	ELISA, WB
TNF Antibody	CSB-PA523593ZA01MQG	TNF	Marmota monax	ELISA, WB
TNF Antibody	CSB-PA004310	TNF	Human, Mouse, Rat	ELISA, WB, IHC, IF
TNF Antibody	CSB-PA004963	TNF	Human	ELISA, WB
TNF Antibody	CSB-PA080267	TNF	Human, Mouse, Rat	WB, IHC
TNF Antibody	CSB-PA080266	TNF	Human, Mouse, Rat	WB
TNF Antibody	CSB-PA954619	TNF	Human	ELISA, IHC
TNF Antibody	CSB-PA023955KA01HU	TNF	Human, Mouse, Rat	ELISA, WB, IHC
TNF Antibody	CSB-PA147864	TNF	Human, Mouse	ELISA, IHC, IF
TNFSF13 Antibody	CSB-PA023989ESR1HU	TNFSF13	Human	ELISA, WB, IHC
TNFSF13 Antibody	CSB-PA023989ESR2HU	TNFSF13	Human	ELISA, IHC
TNFSF13 Antibody	CSB-PA000909	TNFSF13	Human, Mouse, Rat	ELISA, WB, IHC, IF
TNFSF13 Antibody	CSB-PA023989GA01HU	TNFSF13	Human, Mouse, Rat	ELISA, WB
TNFSF13 Antibody	CSB-PA877550	TNFSF13	Human, Mouse	ELISA, IHC
TNFSF13 Antibody	CSB-PA935662	TNFSF13	Human, Mouse	ELISA, IHC
TNFSF13 Antibody	CSB-PA599200	TNFSF13	Human, Mouse	ELISA, IHC
TNFSF13 Antibody	CSB-PA587516	TNFSF13	Human, Mouse	ELISA, IHC
TNFSF13 Antibody	CSB-PA244663	TNFSF13	Human, Mouse, Rat	ELISA, WB, IHC, IF
TNFSF9 Antibody	CSB-PA08047A0Rb	TNFSF9	Human	ELISA, IHC, IF
TNFSF9 Antibody	CSB-PA001424	TNFSF9	Human	ELISA, IF

Product Name	Code	Target	Species	Application
TNFSF9 Antibody	CSB-PA203832	TNFSF9	Human	ELISA, WB, IF
TP53 Antibody	CSB-PA024077LA01MO	TP53	Mouse	ELISA
TP53 Antibody	CSB-PA07889A0Rb	TP53	Human, Mouse, Rat	ELISA, WB, IHC, ChIP
TP53 Antibody	CSB-PA15509A0Rb	TP53	Rat, Mouse	ELISA, WB
Di-Methyl-TP53 (K370) Antibody	CSB-PA000096	TP53	Human	ELISA, WB
Phospho-TP53 (S20) Antibody	CSB-PA000596	TP53	Human, Mouse, Rat, Monkey	ELISA, WB, IHC
Phospho-TP53 (S315) Antibody	CSB-PA000597	TP53	Human	ELISA, WB, IHC, IP
Phospho-TP53 (S33) Antibody	CSB-PA000598	TP53	Human	ELISA, WB
Phospho-TP53 (S46) Antibody	CSB-PA000599	TP53	Human, Monkey	ELISA, WB
Phospho-TP53 (S6) Antibody	CSB-PA000600	TP53	Human, Monkey	ELISA, WB, IHC, IP
Phospho-TP53 (S9) Antibody	CSB-PA000602	TP53	Human, Monkey	ELISA, WB, IHC
Phospho-TP53 (S37) Antibody	CSB-PA000632	TP53	Human	ELISA, WB, IHC
TP53 Antibody	CSB-PA003668	TP53	Human, Monkey	ELISA, WB, IHC, IF
TP53 Antibody	CSB-PA003672	TP53	Human, Monkey	ELISA, WB, IHC, IF
TP53 Antibody	CSB-PA003675	TP53	Human, Mouse, Rat	ELISA, WB
TP53 Antibody	CSB-PA003680	TP53	Human	ELISA, WB, IHC, IF
TP53 Antibody	CSB-PA005967	TP53	Human, Mouse, Rat	ELISA, WB
Acetyl-TP53 (K381) Antibody	CSB-PA050292	TP53	Human, Mouse, Rat	ELISA, WB, IHC, IF
Phospho-TP53 (S366) Antibody	CSB-PA050296	TP53	Human	ELISA, WB, IHC
Acetyl-TP53 (K317) Antibody	CSB-PA050291	TP53	Human	ELISA, WB, IHC, IF
Acetyl-TP53 (K386) Antibody	CSB-PA050293	TP53	Human	ELISA, IHC, IF
Phospho-TP53 (S378) Antibody	CSB-PA050295	TP53	Human, Mouse, Rat	ELISA, WB
Phospho-TP53 (S392) Antibody	CSB-PA050297	TP53	Human, Monkey	ELISA, WB, IHC, IF
Phospho-TP53 (T81) Antibody	CSB-PA050298	TP53	Human	ELISA, IHC
Phospho-TP53 (T387) Antibody	CSB-PA050299	TP53	Human	ELISA, IHC
Phospho-TP53 (S376) Antibody	CSB-PA050300	TP53	Human, Mouse, Rat	ELISA, IHC
TP53 Antibody	CSB-PA050301	TP53	Human	ELISA, WB, IHC
TP53 Antibody	CSB-PA050302	TP53	Human, Mouse, Rat	ELISA, WB, IHC

Product Name	Code	Target	Species	Application
TP53 Antibody	CSB-PA050303	TP53	Human	ELISA, WB
TP53 Antibody	CSB-PA050304	TP53	Human	ELISA, IHC
TP53 Antibody	CSB-PA050305	TP53	Human, Mouse, Rat, Monkey	ELISA, WB, IHC, IF
Acetyl-TP53 (K370) Antibody	CSB-PA000123	TP53	Human	ELISA, WB
Acetyl-TP53 (K372) Antibody	CSB-PA000124	TP53	Human	ELISA, WB
Acetyl-TP53 (K373) Antibody	CSB-PA000135	TP53	Human	ELISA, WB
Acetyl-TP53 (K382) Antibody	CSB-PA000136	TP53	Human	ELISA, WB
Acetyl-TP53 (K320) Antibody	CSB-PA000140	TP53	Human	ELISA, WB
Phospho-TP53 (S15) Antibody	CSB-PA000595	TP53	Human, Rat	ELISA, WB, IHC, IF, IP
Phospho-TP53 (T18) Antibody	CSB-PA000603	TP53	Human, Mouse, Rat	ELISA, WB
TP53 Antibody	CSB-PA003667	TP53	Human, Mouse, Rat, Monkey	ELISA, WB
TP53 Antibody	CSB-PA003670	TP53	Human, Mouse, Rat	ELISA, WB, IHC, IF
TP53 Antibody	CSB-PA003677	TP53	Human, Mouse, Rat	ELISA, WB, IHC
TP53 Antibody	CSB-PA003678	TP53	Human, Rat, Monkey	ELISA, WB, IHC, IF
TP53 Antibody	CSB-PA050306	TP53	Human, Mouse, Rat, Monkey	ELISA, WB, IHC, IF
TP53 Antibody	CSB-PA080174	TP53	Mouse, Rat	WB, IHC
TP53 Antibody	CSB-PA003666	TP53	Human, Monkey	ELISA, WB
TP53 Antibody	CSB-PA003676	TP53	Human	ELISA, WB
Phospho-TP53 (Ser366) Antibody	CSB-PA847762	TP53	Human	ELISA, WB
TP53 Antibody	CSB-PA024077GA01HU	TP53	Human, Mouse, Rat	ELISA, WB, IHC
TP53 Antibody	CSB-PA150453	TP53	Human	ELISA, WB
TP53 (Ab-6) Antibody	CSB-PA277007	TP53	Human, Mouse	ELISA, WB
TP53 (Ab-9) Antibody	CSB-PA599282	TP53	Human	ELISA, WB
TP53 (Ab-18) Antibody	CSB-PA992550	TP53	Human	ELISA, WB
TP53 (Ab-33) Antibody	CSB-PA593512	TP53	Human	ELISA, WB
TP53 (Ab-37) Antibody	CSB-PA923986	TP53	Human	ELISA, WB
TP53 (Ab-46) Antibody	CSB-PA445446	TP53	Human	ELISA, WB
TP53 (Ab-315) Antibody	CSB-PA964395	TP53	Human	ELISA, WB

Product Name	Code	Target	Species	Application
TP53 (Ab-378) Antibody	CSB-PA870744	TP53	Human, Mouse	ELISA, WB
TP53 (Ab-15) Antibody	CSB-PA063047	TP53	Human	ELISA, WB, IHC, IF
Phospho-TP53 (Thr18) Antibody	CSB-PA109326	TP53	Human	ELISA, WB
Phospho-TP53 (Ser20) Antibody	CSB-PA226168	TP53	Human, Mouse, Rat	ELISA, WB
Phospho-TP53 (Ser366) Antibody	CSB-PA449705	TP53	Human, Rat	ELISA, WB
Phospho-TP53 (Ser392) Antibody	CSB-PA052612	TP53	Human	ELISA, WB
Phospho-TP53 (Ser6) Antibody	CSB-PA572207	TP53	Human, Mouse	ELISA, WB, IHC
Phospho-TP53 (Ser37) Antibody	CSB-PA445853	TP53	Human, Mouse, Rat	ELISA, WB, IF
Phospho-TP53 (Ser33) Antibody	CSB-PA029276	TP53	Human	ELISA, WB, IF
Phospho-TP53 (Ser46) Antibody	CSB-PA945066	TP53	Human	ELISA, WB, IF
Phospho-TP53 (Ser315) Antibody	CSB-PA826876	TP53	Human	ELISA, WB, IHC, IF
Phospho-TP53 (Ser15) Antibody	CSB-PA157242	TP53	Human	ELISA, WB, IHC, IF
TPP1 Antibody	CSB-PA024113LA01HU	TPP1	Human	ELISA, IHC
TPP1 Antibody	CSB-PA024113GA01HU	TPP1	Human, Mouse, Rat	ELISA, WB, IHC
Trem2 Antibody	CSB-PA024405LA01MO	Trem2	Mouse, Human	ELISA, WB
TREM2 Antibody	CSB-PA024405LA01HU	TREM2	Human	ELISA, IHC, IF
TREM2 Antibody	CSB-PA436852	TREM2	Human, Mouse	ELISA, WB
TTR Antibody	CSB-PA00340E0Gp	TTR	Human	ELISA
Ttr Antibody	CSB-PA025270DA01MO	Ttr	Mouse	ELISA, IHC
Ttr Antibody	CSB-PA025270HA01MO	Ttr	Mouse	ELISA
TTR Antibody	CSB-PA025270HA01MOV	TTR	Macaca fascicularis	ELISA
TTR Antibody	CSB-PA12077A0Rb	TTR	Human	ELISA, IHC
TTR Antibody	CSB-PA025270GA01HU	TTR	Human	ELISA, WB, IHC
TTR Antibody	CSB-PA143118	TTR	Human, Mouse, Rat	ELISA, IHC
TTR Antibody	CSB-PA041819	TTR	Human, Mouse, Rat	ELISA, WB, IHC
UBE3A Antibody	CSB-PA025488LA01HU	UBE3A	Human	ELISA, IHC, IF
UBE3A Antibody	CSB-PA025488GA01HU	UBE3A	Human, Mouse, Rat	ELISA, WB, IF
VCP Antibody	CSB-PA025813LA01HU	VCP	Human	ELISA, WB, IHC

Product Name	Code	Target	Species	Application
vcp Antibody	CSB-PA772050LA01DIL	vcp	Zebrafish	ELISA, WB
Phospho-VCP (S352) Antibody	CSB-PA060059	VCP	Human, Mouse, Rat	ELISA, WB
VCP Antibody	CSB-PA060060	VCP	Human, Mouse, Rat	ELISA, WB, IHC
VCP Antibody	CSB-PA025813GA01HU	VCP	Human, Mouse, Rat	ELISA, WB, IHC
VEGFA Antibody	CSB-PA025833LA01PI	VEGFA	Pig	ELISA
VEGFA Antibody	CSB-PA07529A0Rb	VEGFA	Human	ELISA, WB
Vegfa Antibody	CSB-PA08249A0Rb	Vegfa	Rat	ELISA, WB, IHC
VEGFA Antibody	CSB-PA004762	VEGFA	Human	ELISA, WB
VEGFA Antibody	CSB-PA060064	VEGFA	Human, Mouse, Rat	ELISA, IHC
VEGFA Antibody	CSB-PA025833GA01HU	VEGFA	Human, Mouse, Rat	ELISA, WB, IHC
VEGFA Antibody	CSB-PA379558	VEGFA	Human	ELISA, WB, IHC
VEGFA Antibody	CSB-PA793155	VEGFA	Human	ELISA, WB
VWF Antibody	CSB-PA12779A0Rb	VWF	Human	ELISA, IHC
VWF Antibody	CSB-PA025960GA01HU	VWF	Human, Mouse, Rat	ELISA, WB, IHC
VWF Antibody	CSB-PA292051	VWF	Human	ELISA, IHC
VWF Antibody	CSB-PA195869	VWF	Human	ELISA, IHC

CUSABIO Featured Products Related to Neuroscience

● Elisa Kits

Product Name	Code	Target	Detection Range	Sensitivity
Human Acetylcholinesterase,AChE ELISA Kit	CSB-E09670h	ACHE	0.94 mU/mL-60 mU/mL	0.23 mU/mL
Rat Acetylcholinesterase,AChE ELISA Kit	CSB-E11304r	ACHE	7.8 pg/mL-500 pg/mL	1.95 pg/mL
Fish Acetylcholinesterase(AChE)ELISA Kit	CSB-E17001Fh	AChE	2.5 ng/mL-40 ng/mL	2.5 ng/mL
Mouse Acetylcholinesterase(AChE)ELISA Kit	CSB-E17521m	AChE	0.312 mU/mL-20 mU/mL	0.078 mU/mL
Human Receptor for advanced glycation end products, RAGE/AGER ELISA Kit	CSB-E09354h	AGER	78 pg/mL-5000 pg/mL	19.5 pg/mL
Rat Advanced glycosylation end product-specific receptor(AGER) ELISA kit	CSB-EL001441RA	AGER	31.25 pg/mL-2000 pg/mL	7.81 pg/mL
Mouse Advanced glycosylation end product-specific receptor(AGER) ELISA kit	CSB-EL001441MO	AGER	0.312 ng/ml-20 ng/ml	0.078 ng/ml
Human RAC-alpha serine/threonine-protein kinase(AKT1) ELISA kit	CSB-EL001553HU	AKT1	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Bovine Albumin(Alb) ELISA kit	CSB-E08664b	ALB	0.313 µg/mL-80 µg/mL	0.148 µg/mL
Mouse Albumin (Alb)ELISA Kit	CSB-E13878m	ALB	0.781 µg/mL-200 µg/mL	0.110 µg/mL
Dog Albumin(Alb) ELISA Kit	CSB-E15757c	ALB	0.273 µg/mL-35 µg/mL	0.266 µg/mL
Sheep Albumin (Alb)ELISA Kit	CSB-E15839Sh	ALB	350 µg/mL-21.88 µg/mL	5.5 µg/mL
Horse Albumin(Alb) ELISA Kit	CSB-E16205Hs	ALB	1.25 µg/mL-80 µg/mL	1.247 µg/mL
Pig Albumin (Alb) ELISA Kit	CSB-E16207p	ALB	0.03 µg/mL-7000 µg/mL	0.03 µg/mL
Goat Albumin (Alb)ELISA Kit	CSB-E15758G	ALB	0.17 µg/mL-350 µg/mL	0.17 µg/mL
Guinea Pig Albumin(Alb) ELISA kit	CSB-E15765Gp	ALB	0.156 µg/mL-10 µg/mL	0.09 µg/mL
Rabbit Albumin(Alb) ELISA Kit	CSB-E16204Rb	ALB	0.273 µg/mL-70 µg/mL	0.117 µg/mL
Human Angiotensin II,ANG-2 ELISA Kit	CSB-E04500h	ANGPT2	39 pg/mL-2500 pg/mL	9.75 pg/mL
Rat Angiotensin II,ANG-2 ELISA Kit	CSB-E07304r	ANGPT2	Request Information	Request Information
Mouse Angiotensin II,ANG-2 ELISA Kit	CSB-E07305m	ANGPT2	0.156 pg/mL-10 pg/mL	0.039 pg/mL
Pig Angiotensin II,ANG-2 ELISA Kit	CSB-E08460p	ANGPT2	0.94 ng/mL-60 ng/mL	0.235 ng/mL
Human apolipoprotein E (Apo-E) ELISA Kit	CSB-E09748h	APOE	15.62 ng/mL-1000 ng/mL	3.9 ng/mL
Rat apolipoprotein E (Apo-E) ELISA Kit	CSB-E09749r	APOE	15.6 ng/mL-1000 ng/mL	15.6 ng/mL
Mouse apolipoprotein E (Apo-E) ELISA Kit	CSB-E09750m	APOE	6.25 ng/mL-400 ng/mL	1.56 ng/mL
Pig apolipoprotein E (Apo-E) ELISA Kit	CSB-E17889p	APOE	31.25 ng/mL-2000 ng/mL	7.81 ng/mL
Monkey Apolipoprotein E(APOE) ELISA kit	CSB-EL001936RH	APOE	6.25 ng/mL-400 ng/mL	1.56 ng/mL
Rat Amyloid beta A4 protein(APP) ELISA kit	CSB-EL001950RA	APP	15.6 ng/mL-1000 ng/mL	3.9 ng/mL
Human β-site APP-Cleaving Enzyme 1,BACE1 ELISA kit	CSB-E09824h	BACE1	3.12 pg/mL-200 pg/mL	0.78 pg/mL

Product Name	Code	Target	Detection Range	Sensitivity
Human Brain derived neurotrophic factor,BDNF ELISA Kit	CSB-E04501h	BDNF	0.3125 ng/mL-20 ng/mL	0.063 ng/mL
Rat Brain derived neurotrophic facor,BDNF ELISA Kit	CSB-E04504r	BDNF	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Mouse Brain derived neurotrophic facor,BDNF ELISA Kit	CSB-E04505m	BDNF	15.6 pg/mL-1000 pg/mL	1.95 pg/mL
Pig Brain-derived neurotrophic factor(BDNF) ELISA kit	CSB-EL002655PI	BDNF	9.375 ng/mL-600 ng/mL	2.35 ng/mL
Monkey Brain-derived neurotrophic factor(BDNF) ELISA kit	CSB-EL002655RH	BDNF	0.78 ng/mL-50 ng/mL	0.195 ng/mL
Human Complement C1s subcomponent(C1S) ELISA kit	CSB-EL003657HU	C1S	25 ng/mL-1600 ng/mL	6.25 ng/mL
Human Complement 3,C3 ELISA Kit	CSB-E08665h	C3	4.7 ng/mL-300 ng/mL	4.031 ng/mL
Pig Complement 3,C3 ELISA Kit	CSB-E06920p	C3	0.625 µg/mL-10 µg/mL	0.625 µg/mL
Rat Complement 3,C3 ELISA Kit	CSB-E08666r	C3	0.625 µg/mL-40 µg/mL	0.156 µg/mL
Mouse Complement 3,C3 ELISA Kit	CSB-E08667m	C3	93.75 ng/mL-6000 ng/mL	23.44 ng/mL
Rabbit Complement 3,C3 ELISA Kit	CSB-E08668Rb	C3	27.3 ng/mL-7000 ng/mL	27.3 ng/mL
Fish Complement 3,C3 ELISA Kit	CSB-E09727s	C3	12.5 µg/mL-500 µg/mL	12.5 µg/mL
Monkey Complement 3(C3)ELISA Kit	CSB-E16528Mk	C3	4.68 ng/mL-300 ng/mL	1.17 ng/mL
Dog Complement 3(C3)ELISA Kit	CSB-EL003888DO	C3	0.78 µg/mL-50 µg/mL	0.195 µg/mL
Mouse Complement 5,C5 ELISA Kit	CSB-E13802m	C5	0.25 ng/mL-16 ng/mL	0.06 ng/mL
Rat Complement fragment 5,C5 ELISA Kit	CSB-E14363r	C5	9.38 ng/mL-600 ng/mL	2.34 ng/mL
Monkey complement C5 (C5) ELISA kit	CSB-EL003995MK	C5	0.625 ng/mL-40 ng/mL	0.156 ng/mL
Human Uncharacterized protein C9orf72(C9orf72) ELISA kit	CSB-EL004345HU	C9orf72	31.25 pg/mL-2000 pg/mL	7.8 pg/mL
Human T-cell surface glycoprotein CD3 zeta chain (CD247) ELISA kit	CSB-EL004904HU	CD247	0.156 ng/mL-10 ng/mL	0.039 ng/mL
Human Programmed Death Ligand-1(PD-L1/CD274) ELISA Kit	CSB-E13644h	CD274	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Mouse Programmed cell death 1 ligand 1(CD274) ELISA kit	CSB-EL004911MO	CD274	31.25 pg/mL-2000 pg/mL	7.81 pg/mL
Human Soluble CD276,sCD276/sB7-H3 ELISA Kit	CSB-E14285h	CD276	3.12 ng/mL-200 ng/mL	0.78 ng/mL
Human Myeloid cell surface antigen CD33(CD33) ELISA kit	CSB-EL004925HU	CD33	125 pg/mL-8000 pg/mL	31.25 pg/mL
Human Soluble Cluster of differentiation 38,sCD38 ELISA Kit	CSB-E09884h	CD38	0.9 ng/mL-60 ng/mL	0.225 ng/mL
Mouse ADP-ribosyl cyclase 1(CD38) ELISA kit	CSB-EL004929MO	CD38	0.625 ng/mL-40 ng/mL	0.156 ng/mL
Human cluster of differentiation 4,CD4 ELISA Kit	CSB-E08956h	CD4	0.625 ng/mL-40 ng/mL	0.156 ng/mL
Rat cluster of differentiation 4,CD4 ELISA Kit	CSB-E12928r	CD4	0.78 U/mL-50 U/mL	0.195 U/mL
Chicken Cluster of Differentiation 4(CD4) ELISA Kit	CSB-E13114C	Cd4	3.125 ng/mL-800 ng/mL	3.125 ng/mL

Product Name	Code	Target	Detection Range	Sensitivity
Mouse Tumor necrosis factor receptor superfamily member 5(CD40) ELISA kit	CSB-EL004936MO	CD40	Request Information	Request Information
Human Soluble Cluster of differentiation 40 ligand, sCD40L ELISA Kit	CSB-E04716h	CD40LG	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Mouse Soluble Cluster of differentiation 40 ligand, sCD40L ELISA Kit	CSB-E04717m	CD40LG	0.274 ng/ml-200 ng/ml.	0.274 ng/ml.
Rat Soluble Cluster of differentiation 40 ligand, sCD40L ELISA Kit	CSB-E07395r	CD40LG	1.87 ng/mL-120 ng/mL	0.47 ng/mL
Human Leukocyte surface antigen CD47(CD47) ELISA kit	CSB-EL004940HU	CD47	3.12 pg/mL-200 pg/mL	0.78 pg/mL
Mouse Leukocyte surface antigen CD47(CD47) ELISA kit	CSB-EL004940MO	CD47	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human CAMPATH-1 antigen(CD52) ELISA kit	CSB-EL004943HU	CD52	78 ng/mL-5000 ng/mL	19.5 ng/mL
Mouse T-lymphocyte activation antigen CD80(CD80) ELISA kit	CSB-EL004959MO	CD80	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human Soluble Cluster of Differentiation 80 (sCD80) ELISA Kit	CSB-E15768h	CD80	125 pg/mL-8000 pg/mL	31.25 pg/mL
Human Cluster of differentiation 8,CD8 ELISA Kit	CSB-E12707h	CD8A	62.5 U/mL-4000 U/mL	15.6 U/mL
Human calcitonin gene related peptide,CGRP ELISA Kit	CSB-E08210h	CGRP	1.56 pg/mL-100 pg/mL	0.39 pg/mL
Rat calcitonin gene related peptide,CGRP ELISA Kit	CSB-E08211r	CGRP	6.25 pg/mL-400 pg/mL	1.56 pg/mL
Pig calcitonin gene related peptide(CGRP)ELISA Kit	CSB-E16438p	CGRP	31.25 pg/mL-2000 pg/mL	7.81 pg/mL
Mouse calcitonin gene related peptide (CGRP) ELISA kit	CSB-EQ027706MO	CGRP	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human Granulocyte-Macrophage Colony Stimulating Factor,GM-CSF ELISA Kit	CSB-E04568h	CSF2	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Mouse Granulocyte-Macrophage Colony Stimulating Factor,GM-CSF ELISA Kit	CSB-E04569m	CSF2	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Rat Granulocyte-Macrophage Colony Stimulating Factor, GM-CSF ELISA Kit	CSB-E04570r	CSF2	15.6 pg/ml - 1000 pg/ml	3.9 pg/ml
Human cytotoxic T lymphocyte associated antigen 4, CTLA-4 ELISA Kit	CSB-E09171h	CTLA4	125 pg/mL-8000 pg/mL	31.25 pg/mL
Mouse cytotoxic T lymphocyte associated antigen 4, CTLA-4 ELISA Kit	CSB-E13637m	CTLA4	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Human dystroglycan 1 (dystrophin-associated glycoprotein 1) (DAG1) ELISA kit	CSB-EL006488HU	DAG1	125 pg/mL-8000 pg/mL	31.25 pg/mL
Human epidermal growth factor receptor,EGFR ELISA Kit	CSB-E12124h	EGFR	0.312 ng/mL-20 ng/mL	0.019 ng/mL
Human Ephrin type-A receptor 4(EPHA4) ELISA kit	CSB-EL007724HU	EPHA4	25 pg/mL-1600 pg/mL	6.25 pg/mL
Human Erythropoietin receptor,EPOR ELISA Kit	CSB-E04537h	EPOR	0.156 ng/ml - 10 ng/ml	0.039 ng/ml
Human epidermal growth factor receptor 2 (sp185/HER2) ELISA Kit	CSB-E11161h	ERBB2	0.156 ng/mL-10 ng/mL	0.04 ng/mL
Human coagulation factor X ,F X ELISA Kit	CSB-E08440h	F10	0.625 ng/mL-40 ng/mL	0.156 ng/mL
Rat coagulation factor X ,F X ELISA Kit	CSB-E08441r	F10	1.56 ng/mL-100 ng/mL	0.39 ng/mL
Mouse coagulation factor X ,F X ELISA Kit	CSB-E08442m	F10	0.312 µg/mL-20 µg/mL	0.078 µg/mL
Pig coagulation factor X (F10) ELISA kit	CSB-EL007915PI	F10	0.312 µg/mL-20 µg/mL	0.027 µg/mL

Product Name	Code	Target	Detection Range	Sensitivity
Human IgG receptor FcRn large subunit p51(FCGRT) ELISA kit	CSB-EL008545HU	FCGRT	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human Fibronectin,FN ELISA Kit	CSB-E04551h	FN1	12.5 ng/mL-800 ng/mL	3.12 ng/mL
Mouse Fibronectin,FN ELISA Kit	CSB-E04552m	FN1	2.7 ng/ml - 2000 ng/ml	2.7 ng/ml
Rat Fibronectin,FN ELISA Kit	CSB-E04553r	FN1	0.78 ng/mL-50 ng/mL	0.195 ng/mL
Bovine Fibronectin,FN ELISA Kit	CSB-E13781B	FN1	0.137 µg/mL-100 µg/mL	0.137 µg/mL
Rabbit Fibronectin(FN1) ELISA kit	CSB-EL008759RB	FN1	1.25 ng/mL-20 ng/mL	1.25 ng/mL
Human RNA-binding protein FUS (FUS/TLS) ELISA kit	CSB-E17376h	FUS	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human Lysosomal alpha-glucosidase(GAA) ELISA kit	CSB-EL009125HU	GAA	25 pg/mL-1600 pg/mL	6.25 pg/mL
Human glial cell line-derived neurotrophic factor (GDNF) ELISA kit	CSB-E04565h	GDNF	78 pg/mL-5000 pg/mL	19.5 pg/mL
Rat glial cell line-derived neurotrophic factor,GDNF ELISA KIT	CSB-E04566r	GDNF	31.25 pg/mL-2000 pg/mL	7.81 pg/mL
Mouse glial cell line-derived neurotrophic factor, GDNF ELISA KIT	CSB-E07341m	GDNF	10.9 pg/ml - 8000 pg/ml	10.9 pg/ml
Human glial fibrillary acidic protein,GFAP ELISA Kit	CSB-E08601h	GFAP	0.625 ng/mL-40 ng/mL	0.156 ng/mL
Rat glial fibrillary acidic protein,GFAP ELISA Kit	CSB-E08602r	GFAP	78 pg/mL-5000 pg/mL	19.5 pg/mL
Mouse glial fibrillary acidic protein,GFAP ELISA Kit	CSB-E08603m	GFAP	3.12 pg/mL-200 pg/mL	0.78 pg/mL
Pig glial fibrillary acidic protein,GFAP ELISA Kit	CSB-E14950p	GFAP	0.1 ng/mL-20 ng/mL	0.05 ng/mL
Mouse Glucagon-like peptide 1 receptor(GLP1R) ELISA kit	CSB-EL009514MO	GLP1R	0.45 ng/ml - 30 ng/ml	0.11 ng/ml
Human Granulins(GRN) ELISA kit	CSB-EL009939HU	GRN	7.8 pg/mL-500 pg/mL	1.95 pg/mL
Mouse β-hexosaminidase A,β-Hex A ELISA Kit	CSB-E04818m	HEXA	3.12 µU/mL-200 µU/mL	0.78 µU/mL
Human Beta-hexosaminidase subunit alpha(HEXA) ELISA kit	CSB-EL010315HU	HEXA	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human Huntingtin(HTT) ELISA kit	CSB-EL010905HU	HTT	12.5 pg/mL-800 pg/mL	3.125 pg/mL
Human iduronate sulfatase,IDS ELISA Kit	CSB-E09471h	IDS	31.25 pg/mL-2000 pg/mL	7.81 pg/mL
Rat Interferon β,IFN-β/IFNB ELISA Kit	CSB-E04845r	IFNB1	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Mouse Interferon β,IFN-β/IFNB ELISA Kit	CSB-E04945m	IFNB1	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human Interferon β,IFN-β/IFNB ELISA Kit	CSB-E09889h	IFNB1	31.25 pg/mL-2000 pg/mL	7.8 pg/mL
Bovine Interferon β,IFN-β/IFNB ELISA Kit	CSB-E13699B	IFNB1	Request Information	Request Information
Mouse Interferon γ ,IFN-γ ELISA Kit	CSB-E04578m(1)	IFNG	0.156 ng/mL-10 ng/mL	0.039 ng/mL
Human Interferon γ (IFN-γ) ELISA Kit	CSB-E04577h	IFNG	0.312 ng/ml-20 ng/ml	0.078 ng/ml
Mouse Interferon γ ,IFN-γ ELISA Kit	CSB-E04578m	IFNG	15.6 pg/mL-1000 pg/mL	3.9 pg/mL

Product Name	Code	Target	Detection Range	Sensitivity
Rat Interferon γ ,IFN-γ ELISA Kit	CSB-E04579r	IFNG	0.625 pg/mL-40 pg/mL	0.156 pg/mL
Pig Interferon γ ,IFN-γ ELISA Kit	CSB-E06794p	IFNG	15.6 pg/mL-1000 pg/mL	7.8 pg/mL
Rabbit Interferon γ ,IFN-γ ELISA Kit	CSB-E06925Rb	IFNG	62.5 pg/mL-4000 pg/mL	15.6 pg/mL
Chicken Interferon γ ,IFN-γ ELISA Kit	CSB-E08550Ch	IFNG	3.125 pg/mL-800 pg/mL	3.125 pg/mL
Bovine Interferon γ ,IFN-γ ELISA Kit	CSB-E09812b	IFNG	31.25 pg/mL-2000 pg/mL	7.8 pg/mL
Monkey Interferon γ ,IFN-γ ELISA Kit	CSB-E10049Mo	IFNG	31.25 pg/mL-2000 pg/mL	7.8 pg/mL
Dog Interferon γ ,IFN-γ ELISA Kit	CSB-E11262c	IFNG	3.12 pg/mL-200 pg/mL	0.78 pg/mL
Goat Interferon γ ,IFN-γ ELISA Kit	CSB-E12966G	IFNG	7.8 pg/mL-500 pg/mL	1.95 pg/mL
Sheep Interferon γ ,IFN-γ ELISA Kit	CSB-E14018Sh	IFNG	5 pg/mL-200 pg/mL	2.5 pg/mL
Hamster interferon, gamma (IFNG) ELISA kit	CSB-EL011050HA	IFNG	1.87 pg/mL-120 pg/mL	0.46 pg/mL
Mouse Insulin-like growth factor 1,IGF-1 ELISA Kit	CSB-E04581m	IGF1	0.156 ng/mL-10 ng/mL	0.151 ng/mL
Rat Insulin-like growth factor 1,IGF-1 ELISA Kit	CSB-E04582r	IGF1	0.156 ng/mL-10 ng/mL	0.151 ng/mL
Pig insulin-like growth factors-1,IGF-1 ELISA Kit	CSB-E06829p	IGF1	6.25 ng/mL-400 ng/mL	1.56 ng/mL
Rabbit Insulin-like growth factor 1,IGF-1 ELISA Kit	CSB-E06994Rb	IGF1	0.156 ng/mL-10 ng/mL	0.039 ng/mL
Bovine Insulin-like growth factor 1,IGF-1 ELISA Kit	CSB-E08893b	IGF1	0.75 ng/mL-60 ng/mL	0.94 ng/mL
Chicken Insulin-like growth factors 1,IGF-1 ELISA Kit	CSB-E09867Ch	IGF1	125 pg/mL-2000 pg/mL	125 pg/mL
Sheep insulin-like growth factor 1,IGF-1 ELISA Kit	CSB-E13753Sh	IGF1	15.6 ng/mL-1000 ng/mL	3.9 ng/mL
Goat insulin-like growth factor 1(IGF-1)ELISA Kit	CSB-E15909G	IGF1	80 ng/mL-1000 ng/mL	20 ng/mL
Dog insulin-like growth factor 1 (IGF1) ELISA kit	CSB-EL011086DO	IGF1	31.25 ng/mL-2000 ng/mL	7.8 ng/mL
Horse insulin-like growth factor 1 (IGF1) ELISA kit	CSB-EL011086HO	IGF1	12.5 ng/mL-800 ng/mL	3.12 ng/mL
Human insulin-like growth factors 1 receptor,IGF-1R ELISA Kit	CSB-E13766h	IGF1R	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Rat Insulin-like Growth Factors 1 Receptor,IGF-1R ELISA Kit	CSB-E13873r	IGF1R	3.12 pg/mL-200 pg/mL	0.8 pg/mL
Human insulin-like growth factors 2 receptor,IGF-2R ELISA Kit	CSB-E11267h	IGF2R	312.5 pg/mL-20000 pg/mL	78.1 pg/mL
Human Interleukin 10,IL-10 ELISA Kit	CSB-E04593h-IS	IL10	3.9 pg/ml-250 pg/ml	0.975 pg/ml
Human Interleukin 10,IL-10 ELISA KIT	CSB-E04593h	IL10	31.25 pg/ml - 2000 pg/ml	7.8 pg/ml
Mouse interleukin 10,IL-10 ELISA KIT	CSB-E04594m	IL10	62.5 pg/ml - 4000 pg/ml	15.6 pg/ml
Rat Interleukin 10,IL-10 ELISA KIT	CSB-E04595r	IL10	3.12 pg/ml-200 pg/ml.	0.78 pg/ml.
Pig interleukin 10,IL-10 ELISA Kit	CSB-E06779p	IL10	6.25 pg/mL-400 pg/mL	1.56 pg/mL

Product Name	Code	Target	Detection Range	Sensitivity
Rabbit Interleukin 10,IL-10 ELISA Kit	CSB-E06897Rb	IL10	31.25 pg/mL-2000 pg/mL	7.81 pg/mL
Dog Interleukin 10,IL-10 ELISA Kit	CSB-E11261c	IL10	1.56 pg/mL-100 pg/mL	0.39 pg/mL
Sheep Interleukin 10,IL-10 ELISA Kit	CSB-E12817Sh	IL10	5 pg/mL-80 pg/mL	5 pg/mL
Chicken Interleukin 10,IL-10 ELISA Kit	CSB-E12835C	IL10	1 pg/mL-200 pg/mL	0.5 pg/mL
Guinea pig interleukin 10,IL-10 ELISA Kit	CSB-E12915Gp	IL10	31.25 pg/mL-2000 pg/mL	7.81 pg/mL
Bovine Interleukin 10,IL-10 ELISA Kit	CSB-E12917B	IL10	5 pg/mL-1000 pg/mL	2.5 pg/mL
Hamster Interleukin-10(IL10) ELISA kit	CSB-EL011580HA	IL10	6.25 pg/mL-400 pg/mL	1.56 pg/mL
Mouse Interleukin-10 receptor subunit alpha(IL10RA) ELISA kit	CSB-EL011581MO	IL10RA	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Bovine Interleukin 12,IL-12 ELISA Kit	CSB-E14010B	IL12A	6.25 pg/mL-400 pg/mL	1.56 pg/mL
Human Interleukin 12,IL-12/P40 ELISA KIT	CSB-E04598h	IL12B	0.156 ng/ml-10 ng/ml	0.039 ng/mL
Mouse interleukin 12,IL-12/P40 ELISA KIT	CSB-E07360m	IL12B	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Pig Interleukin 12,IL-12/P40 ELISA KIT	CSB-E11341p	IL12B	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Goat interleukin 12,IL-12/P40 ELISA KIT	CSB-E12103G	IL12B	9.38 pg/mL-600 pg/mL	2.34 pg/mL
Mouse Interleukin 17,IL-17 ELISA Kit	CSB-E04608m-IS	IL17A	10.938 pg/ml-700 pg/ml	2.735 pg/ml
Mouse Interleukin 17,IL-17 ELISA Kit	CSB-E04608m	IL17A	78 pg/ml - 5000 pg/ml	19.5 pg/ml
Rat Interleukin 17,IL-17 ELISA KIT	CSB-E07451r	IL17A	93.75 pg/mL-6000 pg/mL	23.43 pg/mL
Pig interleukin 17,IL-17 ELISA Kit	CSB-E08057p	IL17A	20 pg/mL- 800 pg/mL	10 pg/mL
Human Interleukin 17A (IL-17A/IL-17) ELISA Kit	CSB-E12819h	IL17A	62.5 pg/ml - 4000 pg/ml	15.6 pg/ml
Mouse Interleukin 1 α ,IL-1 α ELISA kit	CSB-E04621m	IL1A	2.35 pg/ml - 150 pg/ml	0.58 pg/ml
Human Interleukin 1 α ,IL-1 α ELISA Kit	CSB-E04620h	IL1A	31.25 pg/mL-2000 pg/mL	7.8 pg/mL
Rat Interleukin 1 α ,IL-1 α ELISA kit	CSB-E04622r	IL1A	78 pg/ml-5000 pg/ml.	19.5 pg/ml.
Pig Interleukin 1 α ,IL-1 α ELISA kit	CSB-E09725p	IL1A	25 pg/mL-1600 pg/mL	6.25 pg/mL
Monkey Interleukin-1 alpha (IL1A) ELISA kit	CSB-E10041Mo	IL1A	7.8 pg/mL-500 pg/mL	1.95 pg/mL
Human Interleukin 1 β ,IL-1 β ELISA Kit	CSB-E08053h-IS	IL1B	2.344 pg/ml-150 pg/ml	0.586 pg/ml
Mouse Interleukin 1 β ,IL-1 β ELISA Kit	CSB-E08054m-IS	IL1B	12.5 pg/ml-800 pg/ml	3.125 pg/ml
Pig Interleukin 1 β ,IL-1 β ELISA Kit	CSB-E06782p	IL1B	3.9 pg/mL-250 pg/mL	0.975 pg/mL
Rabbit Interleukin 1 β ,IL-1 β ELISA Kit	CSB-E06900Rb	IL1B	7.8 pg/mL-500 pg/mL	1.95 pg/mL
Human Interleukin 1 β ,IL-1 β ELISA Kit	CSB-E08053h	IL1B	7.8 pg/ml - 500 pg/ml	1.95 pg/ml

Product Name	Code	Target	Detection Range	Sensitivity
Mouse Interleukin 1 β ,IL-1 β ELISA Kit	CSB-E08054m	IL1B	Request Information	Request Information
Rat Interleukin 1 β ,IL-1 β ELISA Kit	CSB-E08055r	IL1B	0.156 ng/ml-10 ng/ml.	0.039 ng/ml.
Sheep Interleukin 1 β ,IL-1 β ELISA Kit	CSB-E10115Sh	IL1B	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Bovine Interleukin 1 β (IL- 1 β) ELISA kit	CSB-E12986B	IL1B	62.5 pg/mL-4000 pg/mL	15.6 pg/mL
Dog Interleukin 1 β (IL- 1 β) ELISA Kit	CSB-E13836c	IL1B	7.8 pg/mL-500 pg/mL	1.95 pg/mL
Guinea pig Interleukin-1 beta(IL1B) ELISA kit	CSB-EL011614GU	IL1B	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human soluble interleukin-1 receptor I ,IL-1sR I ELISA kit	CSB-E04611h	IL1R1	0.156 ng/ml-10 ng/ml	0.039 ng/ml
Mouse Interleukin 2,IL-2 ELISA kit	CSB-E04627m-IS	IL2	3.9 pg/ml-250 pg/ml	0.975 pg/ml
Human Interleukin 2,IL-2 ELISA kit	CSB-E04626h	IL2	31.25 pg/ml-2000 pg/ml.	7.81 pg/ml
Mouse Interleukin 2,IL-2 ELISA kit	CSB-E04627m	IL2	31.25 pg/ml - 2000 pg/ml	7.8 pg/ml
Rat Interleukin 2,IL-2 ELISA kit	CSB-E04628r	IL2	6.25 pg/mL-400 pg/mL	1.56 pg/mL
Pig interleukin 2,IL-2 ELISA Kit	CSB-E06783p	IL2	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Rabbit Interleukin 2,IL-2 ELISA Kit	CSB-E08745Rb	IL2	31.25 pg/mL-2000 pg/mL	7.8 pg/mL
Bovine Interleukin 2,IL-2 ELISA Kit	CSB-E09810b	IL2	31.25 pg/mL-2000 pg/mL	7.8 pg/mL
Sheep Interleukin 2,IL-2 ELISA Kit	CSB-E11217Sh	IL2	31.25 pg/mL-2000 pg/mL	7.81 pg/mL
Dog Interleukin 2,IL-2 ELISA Kit	CSB-E11258c	IL2	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Goat Interleukin 2(IL-2) ELISA Kit	CSB-E12965G	IL2	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human soluble interleukin-2 receptor,IL-2sR α ELISA kit	CSB-E04629h	IL2RA	78 pg/mL-5000 pg/mL	19.5 pg/mL
Mouse soluble interleukin-2 receptor,sIL-2R α ELISA kit	CSB-E07440m	IL2RA	15.6 pg/ml - 1000 pg/ml	3.9 pg/ml
Human Interleukin 6,IL-6 ELISA KIT	CSB-E04638h	IL6	7.8 pg/mL-500 pg/mL	2.453 pg/mL
Mouse Interleukin 6,IL-6 ELISA Kit	CSB-E04639m-IS	IL6	3.9 pg/ml-250 pg/ml	0.975 pg/ml
Mouse Interleukin 6,IL-6 ELISA KIT	CSB-E04639m	IL6	Request Information	Request Information
Rat Interleukin 6,IL-6 ELISA KIT	CSB-E04640r	IL6	Request Information	Request Information
Pig Interleukin 6,IL-6 ELISA Kit	CSB-E06786p	IL6	1.25 pg/mL-80 pg/mL	0.31 pg/mL
Rabbit Interleukin 6,IL-6 ELISA Kit	CSB-E06903Rb	IL6	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Chicken Interleukin 6,IL-6 ELISA Kit	CSB-E08549Ch	IL6	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Monkey Interleukin 6,IL-6 ELISA Kit	CSB-E10050Mo	IL6	31.25 pg/mL-2000 pg/mL	7.8 pg/mL
Sheep Interleukin 6,IL-6 ELISA Kit	CSB-E10116Sh	IL6	2 pg/mL-80 pg/mL	2 pg/mL

Product Name	Code	Target	Detection Range	Sensitivity
Dog Interleukin 6,IL-6 ELISA Kit	CSB-E11260c	IL6	1.56 pg/mL-100 pg/mL	0.39 pg/mL
Bovine Interleukin 6,IL-6 ELISA Kit	CSB-E12899B	IL6	5 pg/mL-1000 pg/mL	2.5 pg/mL
Goat Interleukin 6,IL-6 ELISA Kit	CSB-E14360G	IL6	1.56 pg/mL-100 pg/mL	0.39 pg/mL
Human Interleukin 6 receptor(IL-6R)ELISA Kit	CSB-E16524h	IL6R	0.156 ng/mL-10 ng/mL	0.039 ng/mL
Rat Interleukin 6 receptor,IL-6R ELISA Kit	CSB-E11876r	IL6R	0.625 ng/mL-40 ng/mL	0.156 ng/mL
Mouse glucoprotein 130,gp130 ELISA KIT	CSB-E04572m	IL6ST	Request Information	Request Information
Human Soluble Interleukin-7 Receptor,sIL-7R ELISA Kit	CSB-E14020h	IL7R	0.78 ng/mL-50 ng/mL	0.195 ng/mL
Mouse Interleukin-7 receptor subunit alpha(IL7R) ELISA kit	CSB-EL011670MO	IL7R	78 pg/ml - 5000 pg/ml	19.5 pg/ml
Human Insulin,INS ELISA Kit	CSB-E05069h	INS	Request Information	Request Information
Rabbit Insulin,INS ELISA Kit	CSB-E06992Rb	INS	10 µIU/ml-160 µIU/ml.	2 µIU/ml
Monkey Insulin,INS ELISA Kit	CSB-E10047Mo	INS	Request Information	Request Information
Sheep Insulin(INS) ELISA Kit	CSB-E17044Sh	INS	Request Information	Request Information
Goat Insulin(INS)ELISA Kit	CSB-E17653G	INS	Request Information	Request Information
Rat Insulin,INS ELISA Kit	CSB-E05070r	INS	15.6 nIU/mL-1000 nIU/mL	3.9 nIU/mL
Mouse Insulin,INS ELISA Kit	CSB-E05071m	INS	15.6 nIU/mL-1000 nIU/mL	3.9 nIU/mL
Pig Insulin,INS ELISA Kit	CSB-E06828p	INS	2 µIU/mL-400 µIU/mL	1 µIU/mL
Fish Insulin,INS ELISA Kit	CSB-E12123Fh	INS	100 pg/mL-1600 pg/mL	100 pg/mL
Chicken Insulin ELISA kit	CSB-E13293C	INS	3.12 µIU/mL-50 µIU/mL	0.59 µIU/mL
Guinea pig insulin (INS) ELISA kit	CSB-EL011742GU	INS	12 µIU/mL-210 µIU/mL	7.5 µIU/mL
Horse insulin (INS) ELISA kit	CSB-EL011742HO	INS	1.56 nIU/mL-100 nIU/mL	0.39 nIU/mL
Human cytokeratin 20,CK20 ELISA Kit	CSB-E11710h	KRT20	0.625 ng/mL-40 ng/mL	0.156 ng/mL
Human low density lipoprotein receptor,LDLR ELISA Kit	CSB-E08950h	LDLR	31.25 pg/mL-2000 pg/mL	7.8 pg/mL
Rat low density lipoprotein receptor,LDLR ELISA Kit	CSB-E13545r	LDLR	7.8 pg/mL-500 pg/mL	1.95 pg/mL
Human Tau proteins ELISA kit	CSB-E12011h	MAPT	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Rat Tau Proteins ELISA kit	CSB-E13729r	MAPT	62.5 pg/mL-4000 pg/mL	15.6 pg/mL
Human myelin basic protein,MBP ELISA Kit	CSB-E08283h	MBP	7.8 ng/ml- 500ng/ml	1.95 ng/ml
Rat myelin basic protein,MBP ELISA Kit	CSB-E08284r	MBP	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Mouse myelin basic protein,MBP ELISA Kit	CSB-E08285m	MBP	18.75 pg/mL-1200 pg/mL	4.69 pg/mL

Product Name	Code	Target	Detection Range	Sensitivity
Human Myelin-oligodendrocyte glycoprotein(MOG) ELISA kit	CSB-EL014709HU	MOG	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human B-lymphocyte antigen (CD20) ELISA Kit	CSB-E14007h	MS4A1	3.12 ng/mL-200 ng/mL	0.78 ng/mL
Human Myostatin,MSTN ELISA Kit	CSB-E11300h	MSTN	0.625 ng/mL-20 ng/mL	0.312 ng/mL
Bovine Growth/differentiation factor 8(MSTN) ELISA kit	CSB-EL015057BO	MSTN	6.25 ng/mL-400 ng/mL	1.56 ng/mL
Rat Growth/differentiation factor 8(MSTN) ELISA kit	CSB-EL015057RA	MSTN	0.78 ng/mL-50 ng/mL	0.195 ng/mL
Human mammalian target of rapamycin,mTOR ELISA Kit	CSB-E09038h	mTOR	1.56 pg/mL-100 pg/mL	0.39 pg/mL
Human Nuclear factor-kappa B,NF-κB ELISA Kit	CSB-E12107h	NFKB1	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Mouse Nuclear factor-kappa B,NF-κB ELISA Kit	CSB-E12108m	NFKB1	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Rat Nuclear factor-kappa B(NF-κB) ELISA Kit	CSB-E13148r	NFKB1	1.56 pg/mL-100 pg/mL	0.39 pg/mL
Human Nerve growth factor,NGF ELISA kit	CSB-E04683h	NGF	6.86 pg/ml-5000 pg/ml.	6.86 pg/ml.
Mouse Nerve growth factor,NGF ELISA kit	CSB-E04684m	NGF	62.5 pg/ml-4000 pg/ml	15.6 pg/ml
Rat Nerve growth factor,NGF ELISA kit	CSB-E04685r	NGF	0.78 pg/mL-50 pg/mL	0.195 pg/mL
Human Ornithine carbamoyltransferase, mitochondrial (OTC) ELISA kit	CSB-EL017270HU	OTC	1.25 mU/mL-80 mU/mL	0.312 mU/mL
Mouse Ornithine carbamoyltransferase, mitochondrial (OTC) ELISA kit	CSB-EL017270MO	OTC	3.12 mU/mL-200 mU/mL	0.78 mU/mL
Mouse Phenylalanine-4-hydroxylase(PAH) ELISA kit	CSB-EL017396MO	PAH	23.5 pg/mL-1500 pg/mL	5.8 pg/mL
Human Protein DJ-1 (PARK7) ELISA kit	CSB-E12024h	PARK7	0.235 ng/mL-15 ng/mL	0.058 ng/mL
Mouse Protein DJ-1(PARK7) ELISA kit	CSB-EL017454MO	PARK7	62.5 pg/mL-4000 pg/mL	15.6 pg/mL
Human Proprotein convertase subtilisin/kexin type 9 (PCSK9) ELISA kit	CSB-EL017647HU	PCSK9	0.45 ng/mL-30 ng/mL	0.11 ng/mL
Rat Proprotein convertase subtilisin/kexin type 9 (PCSK9) ELISA kit	CSB-EL017647RA	PCSK9	93.75 pg/mL-6000 pg/mL	46.88 pg/mL
Monkey Proprotein convertase subtilisin/kexin type 9 (PCSK9) ELISA kit	CSB-EL017647RH	PCSK9	1.1 pg/mL-800 pg/mL	1.1 pg/mL
Mouse Programmed Death 1(PD-1)ELISA Kit	CSB-E13586m	PDCD1	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Human Programmed Death 1(PD-1)ELISA KIT	CSB-E13643h	PDCD1	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human progranulin(PGRN) ELISA Kit	CSB-E15049h	PGRN	4.7 ng/mL-300 ng/mL	1.17 ng/mL
Human Plasminogen,Plg ELISA Kit	CSB-E09177h	PLG	62.5 ng/mL-4000 ng/mL	15.6 ng/mL
Rabbit Plasminogen,Plg ELISA Kit	CSB-E11871Rb	PLG	62.5 ng/mL-4000 ng/mL	15.6 ng/mL
Rat Plasminogen (Plg) ELISA kit	CSB-E17250r	Plg	62.5 ng/mL-4000 ng/mL	15.6 ng/mL
Bovine Plasminogen(PLG) ELISA kit	CSB-EL018188BO	PLG	62.5 ng/mL-4000 ng/mL	15.6 ng/mL
Mouse Plasminogen(PLG) ELISA kit	CSB-EL018188MO	PLG	62.5 ng/mL-4000 ng/mL	15.63 ng/mL

Product Name	Code	Target	Detection Range	Sensitivity
Human Apoptosis-associated speck-like protein containing a CARD(PYCARD) ELISA kit	CSB-EL019114HU	PYCARD	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Mouse Apoptosis-associated speck-like protein containing a CARD(PYCARD) ELISA kit	CSB-EL019114MO	PYCARD	7.81 pg/mL-500 pg/mL	1.95 pg/mL
Human Reticulon-4(RTN4) ELISA kit	CSB-EL020572HU	RTN4	46.88 pg/mL-3000 pg/mL	11.75 pg/mL
Mouse P-Selectin ELISA kit	CSB-E04709m	SELP	31.25 pg/mL-2000 pg/mL	7.8 pg/mL
Rat P-Selectin ELISA kit	CSB-E07399r	SELP	9.38 ng/mL-600 ng/mL	2.34 ng/mL
Dog P-selectin(SELP) ELISA kit	CSB-EL020978DO	SELP	12.5 ng/mL-800 ng/mL	3.12 ng/mL
Human soluble Cluster of differentiation 100(sCD100)ELISA kit	CSB-E15792h	SEMA4D	125 pg/mL-8000 pg/mL	31.25 pg/mL
Human NAD-dependent deacetylase sirtuin-3, mitochondrial (SIRT3) ELISA kit	CSB-EL021341HU	SIRT3	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human Alpha-synuclein (SNCA) ELISA Kit	CSB-E18033h	SNCA	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Rat Alpha-synuclein(SNCA) ELISA kit	CSB-EL021912RA	SNCA	0.625 ng/mL-40 ng/mL	0.156 ng/mL
Mouse Super Oxidase Dimutase,SOD ELISA Kit	CSB-E08556m	SOD1	5.4 pg/mL-4000 pg/mL	5.4 pg/mL
Human Superoxide dismutase [Cu-Zn] (SOD1) ELISA kit	CSB-E16845h	SOD1	62.5 pg/ml - 4000 pg/ml	15.6 pg/ml
Bovine Cu/Zn-Superoxide Dismutase,Cu/Zn-SOD ELISA Kit	CSB-E14090B	SOD1	15.6 ng/mL-1000 ng/mL	3.9 ng/mL
Rat superoxide dismutase 1, soluble (SOD1) ELISA kit	CSB-EL022397RA	SOD1	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human Sortilin(SORT1) ELISA kit	CSB-EL022412HU	SORT1	46.88 pg/mL-3000 pg/mL	11.72 pg/mL
Mouse Sortilin(SORT1) ELISA kit	CSB-EL022412MO	SORT1	6.25 pg/mL-400 pg/mL	1.56 pg/mL
Human TAR DNA-binding protein 43 (TARDBP/TDP43) ELISA kit	CSB-E17007h	TARDBP	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Mouse transferrin receptor,TFR ELISA Kit	CSB-E08389m	TFRC	0.625 ng/mL-40 ng/mL	0.156 ng/mL
Human Transforming Growth factor β 1,TGF- β 1 ELISA kit	CSB-E04725h	TGFB1	0.78 ng/mL-50 ng/mL	0.747 ng/mL
Dog Transforming growth factor β 1,TGF- β 1 ELISA Kit	CSB-E04796c	TGFB1	0.781 ng/ml - 50 ng/ml	0.747 ng/ml
Mouse Transforming Growth factor β 1,TGF- β 1 ELISA kit	CSB-E04726m	TGFB1	0.78 ng/mL-50 ng/mL	0.2 ng/mL
Rat Transforming Growth Factor β 1(TGF- β 1) ELISA kit	CSB-E04727r	TGFB1	6.25 pg/mL-400 pg/mL	1.56 pg/mL
Pig transforming growth factor β 1(TGF- β 1)ELISA Kit	CSB-E06843p	TGFB1	0.78 ng/mL-50 ng/mL	0.195 ng/mL
Chicken Transforming Growth factor β 1,TGF- β 1 ELISA kit	CSB-E09875Ch	TGFB1	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Bovine transforming growth factor β 1,TGF- β 1 ELISA Kit	CSB-E14208B	TGFB1	6.25 ng/mL-400 ng/mL	1.56 ng/mL
Fish transforming growth factor, beta 1 (TGFB1) ELISA kit	CSB-EL023446FI	TGFB1	6.25 ng/mL-400 ng/mL	3.12 ng/mL
Horse Transforming growth factor beta-1(TGFB1) ELISA kit	CSB-EL023446HO	TGFB1	3.12 pg/mL-200 pg/mL	1.56 pg/mL
Sheep Transforming growth factor beta-1(TGFB1) ELISA kit	CSB-EL023446SH	TGFB1	15.6 ng/mL-1000 ng/mL	3.9 ng/mL

Product Name	Code	Target	Detection Range	Sensitivity
Human tyrosine hydroxylase,TH ELISA Kit	CSB-E09661h	TH	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Rat Tyrosine Hydroxylase(TH) ELISA Kit	CSB-E13102r	TH	78 pg/mL-5000 pg/mL	19.5 pg/mL
Mouse tyrosine hydroxylase(TH)ELISA Kit	CSB-E17645m	TH	62.5 pg/mL-4000 pg/mL	15.6 pg/mL
Human toll-like receptor 2,TLR2 ELISA Kit	CSB-E11781h	TLR2	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Human toll-like receptor 4,TLR4 ELISA Kit	CSB-E12954h	TLR4	0.156 ng/mL-10 ng/mL	0.039 ng/mL
Mouse Soluble Toll-like Receptor 4,sTLR4 ELISA Kit	CSB-E14275m	TLR4	0.156 ng/mL-10 ng/mL	0.039 ng/mL
Rat toll-like receptor 4 (TLR4) ELISA Kit	CSB-E15822r	TLR4	0.625 ng/mL-40 ng/mL	0.156 ng/mL
Pig Tumor necrosis factor (TNF/TNFA/TNFSF2) ELISA kit	CSB-E16980p	TNF	0.047 ng/mL-3 ng/mL	0.040 ng/mL
Mouse Tumor necrosis factor α ,TNF- α ELISA KIT	CSB-E04741m(1)	TNF	62.5 pg/mL-4000 pg/mL	15.6 pg/mL
Mouse Tumor necrosis factor α ,TNF- α ELISA Kit	CSB-E04741m-IS	TNF	7.8 pg/ml-500 pg/ml	1.95 pg/ml
Human Tumor necrosis factor α ,TNF- α ELISA KIT	CSB-E04740h	TNF	7.8 pg/mL-500 pg/mL	1.95 pg/mL
Mouse Tumor necrosis factor α ,TNF- α ELISA KIT	CSB-E04741m	TNF	15.6 pg/ml-1000 pg/ml	3.9 pg/ml
Guinea pig Tumor necrosis factor α ,TNF- α ELISA Kit	CSB-E06772Gu	TNF	3.12 pg/mL-200 pg/mL	0.78 pg/mL
Rabbit Tumor necrosis factor α ,TNF α ELISA Kit	CSB-E06998Rb	TNF	78 pg/mL-5000 pg/mL	19.5 pg/mL
Goat Tumor necrosis factor α ,TNF- α ELISA KIT	CSB-E09811g	TNF	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Monkey Tumor necrosis factor α ,TNF- α ELISA KIT	CSB-E10039Mo	TNF	1.56 pg/mL-100 pg/mL	0.39 pg/mL
Dog Tumor necrosis factor α ,TNF- α ELISA Kit	CSB-E11737c	TNF	0.156 ng/mL-10 ng/mL	0.039 ng/mL
Rat TNF- α ELISA kit	CSB-E11987r	TNF	6.25 pg/mL-400 pg/mL	1.56 pg/mL
Bovine TNF- α ELISA kit	CSB-E12020B	TNF	0.1 ng/mL-20 ng/mL	0.05 ng/mL
Sheep Tumor Necrosis Factor α ,TNF- α ELISA Kit	CSB-E13853Sh	TNF	3.125 pg/mL-50 pg/mL	3.125 pg/mL
Horse Tumor necrosis factor α ,TNF- α ELISA Kit	CSB-E14031Hs	TNF	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human a proliferation inducing ligand(APRIL)ELISA Kit	CSB-E15012h	TNFSF13	0.78 ng/mL-50 ng/mL	0.195 ng/mL
Mouse Tumor necrosis factor ligand superfamily member 13(TNFSF13) ELISA kit	CSB-EL023989MO	TNFSF13	93.75 pg/mL-6000 pg/mL	23.43 pg/mL
Human Tumor necrosis factor ligand superfamily member 9(TNFSF9) ELISA kit	CSB-EL023997HU	TNFSF9	18.75 pg/mL-1200 pg/mL	4.7 pg/mL
Mouse Tumor necrosis factor ligand superfamily member 9(TNFSF9) ELISA kit	CSB-EL023997MO	TNFSF9	23.5 pg/mL-1500 pg/mL	5.8 pg/mL
Human p53/tumor protein,p53/TP53 ELISA Kit	CSB-E08334h	TP53	9.38 pg/mL-600 pg/mL	2.34 pg/mL
Rat p53/tumor protein,p53/TP53 ELISA Kit	CSB-E08336r	TP53	12.5 pg/mL-800 pg/mL	3.12 pg/mL
Human Triggering receptor expressed on myeloid cells 2(TREM2) ELISA kit	CSB-EL024405HU	TREM2	0.312 ng/mL-20 ng/mL	0.078 ng/mL

Product Name	Code	Target	Detection Range	Sensitivity
Mouse Triggering receptor expressed on myeloid cells 2 (TREM2) ELISA kit	CSB-EL024405MO	TREM2	0.156 ng/mL-10 ng/mL	0.039 ng/mL
Rat triggering receptor expressed on myeloid cells 2 (TREM2) ELISA kit	CSB-EL024405RA	TREM2	15.6 pg/mL-1000 pg/mL	3.9 pg/mL
Human transthyretin(TTR) ELISA kit	CSB-E11169h	TTR	1.56 µg/mL-25 µg/mL	0.78 µg/mL
Bovine Transthyretin(TTR) ELISA kit	CSB-EL025270BO	TTR	312.5 ng/mL-5000 ng/mL	312.5 ng/mL
Mouse Transthyretin(TTR) ELISA kit	CSB-EL025270MO	TTR	3.12 ng/mL-200 ng/mL	0.78 ng/mL
Human Transitional endoplasmic reticulum ATPase(VCP) ELISA kit	CSB-EL025813HU	VCP	23.44 pg/mL-1500 pg/mL	5.86 pg/mL
Mouse Vascular Endothelial cell Growth Factor,VEGF ELISA KIT	CSB-E04756m	VEGFA	3.906 pg/mL-250 pg/mL	0.857 pg/mL
Human Vascular Endothelial cell Growth Factor(VEGF) ELISA KIT	CSB-E11718h	VEGFA	31.25 pg/mL-2000 pg/mL	25.297 pg/mL
Rat Vascular Endothelial cell Growth Factor,VEGF ELISA Kit	CSB-E04757r	VEGFA	3.9 pg/mL-250 pg/mL	0.97 pg/mL
Dog Vascular Endothelial cell Growth Factor,VEGF ELISA KIT	CSB-E11744c	VEGFA	3.12 pg/mL-200 pg/mL	0.78 pg/mL
Pig Vascular Endothelial cell Growth Factor,VEGF ELISA KIT	CSB-E12053p	VEGFA	4.69 pg/mL-300 pg/mL	1.17 pg/mL
Bovine Vascular endothelial growth factor A(VEGFA) ELISA kit	CSB-EL025833BO	VEGFA	47 pg/mL-3000 pg/mL	11.75 pg/mL
Human von Willebrand Factor,vWF ELISA Kit	CSB-E08437h	VWF	6.25 ng/ml-400 ng/ml.	1.56 ng/ml.
Rat von Willebrand Factor,vWF ELISA Kit	CSB-E08438r	VWF	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Mouse von Willebrand Factor,vWF ELISA Kit	CSB-E08439m	VWF	0.78 ng/mL-50 ng/mL	0.195 ng/mL
Rabbit von Willebrand Factor,vWF ELISA Kit	CSB-E08525Rb	VWF	0.312 ng/mL-20 ng/mL	0.078 ng/mL
Pig von Willebrand Factor,vWF ELISA Kit	CSB-E09395p	VWF	0.625 ng/mL-40 ng/mL	0.156 ng/mL
Dog von Willebrand Factor,vWF ELISA Kit	CSB-E12759c	VWF	0.625 ng/mL-40 ng/mL	0.32 ng/mL

CUSABIO Featured Products Related to Neuroscience

● Cusabio Cytokines for Neural Cell Culture

Cytokine	Role	CUSABIO Product			
		Product Code	Organism	Source	Tag Info
BDNF	Brain-derived neurotrophic factor, supports neuronal survival and promotes maturation, essential for learning and memory processes.	CSB-AP003781HU	Human	E.coli	Tag-Free
CNTF	Ciliary neurotrophic factor, which promotes nerve cell survival and has potential therapeutic effects in neurodegenerative diseases.	CSB-AP002841HU	Human	E.Coli	Tag-Free
EGF	Epidermal growth factor, which promotes the growth and proliferation of many cell types, especially epithelial cells.	CSB-AP003711HU	Human	E.coli	Tag-Free
		CSB-AP004091MO	Mouse	E.coli	C-terminal 6xHis-tagged
aFGF (FGF1)	Acidic Fibroblast Growth Factor , promotes cell proliferation and differentiation, especially in the repair of damaged tissues.	CSB-AP003761HU	Human	E.coli	Tag-Free
bFGF (FGF2)	Basic Fibroblast Growth Factor ,promotes the proliferation of neural stem cells and precursor cells while contributing to wound healing and angiogenesis	CSB-AP003751HU	Human	E.coli	Tag-Free
		CSB-AP003831HU	Human	E.coli	Tag-Free
		CSB-AP003841HU	Human	E.coli	Tag-Free
		CSB-AP004061MO	Mouse	E.coli	Tag-Free
		CSB-AP004171RA	Rat	E.coli	Tag-Free
GDNF	Glial cell-derived neurotrophic factor, important for the survival and regeneration of dopaminergic neurons.	CSB-AP002881HU	Human	E.Coli	Tag-Free
IGF1	Insulin-like growth factor, which promotes cell proliferation and differentiation and has a positive effect on the development and regeneration of the nervous system.	CSB-AP003701HU	Human	E.coli	Tag-Free
		CSB-AP003721HU	Human	E.coli	Tag-Free
		CSB-AP003731HU	Human	E.coli	Tag-Free
NGF	Nerve growth factor, promotes nerve cell growth and survival, especially in the development of sensory and sympathetic neurons.	CSB-AP003771HU	Human	E.coli	Tag-Free
		CSB-AP003951HU	Human	Mammalian cell	Tag-Free
		CSB-AP004081MO	Mouse	E.coli	Tag-Free